

INFORMATION

Customer Name: Date:
 Address: PO#:
 City: Phone:
 Account:
 Authorizer Name:

STANDARD NO CHARGE SPECIFICATIONS

SEAT WIDTH: 14" – 18"
 SEAT DEPTH [ADJUSTABLE]
 SEAT HEIGHT [ADJUSTABLE]: 12.25" – 18.25"
 BACK ANGLE [ADJUSTABLE]: 0°, 5°, 10°, 15° [DEGREES]
 WEIGHT OF CHAIR: 60 LBS [APPROXIMATELY]
 SINGLE RELEASE, DUAL CYLINDERS
 ADJUSTABLE STROLLER HANDLE: 0° – 180° [DEGREES]
 TILT ADJUSTMENTS SET AT 35° [DEGREES]
 STANDARD LOCKOUT COLLARS
 WEIGHT CAPACITY: 250 LBS
 TWO-WAY LOCKING GAS STRUT
 STANDARD T STYLE ARMS
 REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

FRAME TYPE 01

- STPB Tilt Frame – Black \$2,699
 Tilt Adjustment: 0–35° 0–45° N/C
- STPINTS Independant Power Tilt System (A1) \$2,000

CLIENT WEIGHT 02

- STPWA Up to 250Lbs N/C
- STPWAMF Exceeds 250Lbs \$730
- STPWAMG Exceeds 350Lbs \$900

SEAT WIDTH 03

- STPSW14 Seat Width 14" N/C
- STPSW15 Seat Width 15" N/C
- STPSW16 Seat Width 16" N/C
- STPSW17 Seat Width 17" N/C
- STPSW18 Seat Width 18" N/C
- STPSW19 Seat Width 19" (A5) – W111 \$200
- STPSW20 Seat Width 20" (A5) – W111 \$200
- STPSW21–24 Specify: _____ (A5) – W111 \$450
- STPSW25–32 Specify: _____ \$450

FOOT REST & FOOT PLATES 04

- STPOFST Offset Foot Rest N/C
- STPF6A 60 Degree Hanger Pivot Pin N/C
- STPF7A 70 Degree Hanger Pivot Pin N/C
- STPF8A 80 Degree Hanger Pivot Pin N/C
- STPEA Elevating Legrest Pivot Pin (A2) – W106 \$404
- STPAL Articulating Legrest Pivot Pin (A2) – W106 \$450
- STPNR Non-removable Pin Style Legrest N/C
- STPEB Elevating Legrest System \$450
- STPPPL Plastic Footplates N/C
- STPPXL Extra Large Footplates \$150
- STPPAF Angle Adjustable Footplates – W107 \$150
- STPIFB One Piece Foot Board Padded \$80
- STPFRE Foot Rest Extensions Short Medium Long N/C

SEAT DEPTH 05

- STPSD16 Seat Depth 16" N/C
- STPSD17 Seat Depth 17" N/C
- STPSD18 Seat Depth 18" N/C
- STPSD19 Seat Depth 19" \$200
- STPSD20 Seat Depth 20" \$200
- STPSD21–24 Specify: _____ \$450

BACK TYPE 06

- STPAB1 Angle Adjustable Back – Preset N/C
 0° 5° 10° 15°
- STPENDO Endomorph Canes W () N/C
- STPRBK Recline Kit \$560

SEAT TO FLOOR SPECIFICATIONS (± 0.5")

(LOW SEAT SETTING)

STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 13.00"	20"	5"	XL SHORT
<input type="checkbox"/> 13.75"	20"	5"	XL SHORT
<input type="checkbox"/> 14.50"	20"	6"	SHORT
<input type="checkbox"/> 14.50"	22"	5"	SHORT
<input type="checkbox"/> 14.50"	22"	6"	SHORT
<input type="checkbox"/> 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/> 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/> 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/> 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/> 16.00"	22"	6"	LONG

(HIGH SEAT SETTING)

STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 16.75"	22"	6"	LONG
<input type="checkbox"/> 17.50"	22"	6"	LONG
<input type="checkbox"/> 17.50"	22"	6"	LONG
<input type="checkbox"/> 17.50"	24"	6"	LONG
<input type="checkbox"/> 18.25"	24"	8"	LONG

ARM REST Select one below **07**

- STPTA T-Arm Ht. Adj./ Removable (A3) N/C
- STPFA Flip Back Ht. Adj./ Removable N/C
- STPTAS T-Arm Ht. Adj./ Removable (A3) N/C
- STPFAS Flip Back Ht. Adj./ Removable N/C
- STPGWAPF Gel Waterfall Arm Pads (A1) \$150

REAR WHEELS Select one below **08**

- STPRWTO 20" Solid Treaded Tire N/C
- STPRWT2 22" Solid Treaded Tire N/C
- STPRWT4 24" Solid Treaded Tire N/C
- STPRWT12I 12" Solid Tire Inside Mount 120
- STPRWT16 16" Solid Tire N/C
- STPAQR Axles – Quick Release \$120
- STPHRT Narrow Handrim Config. (1/2" less) N/C
- STPPCHR Plastic Coated Hand Rims N/C

CASTER Select one below **09**

- STPFC1D 3" X 1" Caster N/C
- STPFC1A 5" X 1.25" Caster N/C
- STPFC1B 6" X 1.25" Caster N/C
- STPFC1C 8" X 1.25" Caster N/C
- STPFC2A 5" X 1.75" Caster (Soft Roll) \$80
- STPFC2B 6" X 1.75" Caster \$80
- STPFC2C 8" X 1.75" Caster \$80
- STPFC4B 6" X 2 Caster (Steel Fork) \$80
- STPFC4C 8" X 2 Caster (Steel Fork) \$80

WHEEL LOCK Select one below **10**

- STPALPH Aluminium Push To Lock N/C
- STPALPL Aluminium Pull To Lock N/C
- STPSTPH Steel Push To Lock N/C
- STPSTPL Steel Pull To Lock N/C
- STPBRXC 4" Brake Extension \$40
- STPBRXA 6" Brake Extension \$40
- STPBRXB 9" Brake Extension \$130

SEATING BENIFITS AND ACCESSORIES

Kindly refer to new basic seating benifit order form

PLUS PACKAGES Select one below **15**

- PLUS PACKAGE 1:** CLIENT 250 – 349 LBS - \$730
AADL CODE: W101 INCLUDES THE FOLLOWING:
 STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS,
 HEAVY DUTY STRUTS AND SEAT PAN
 SEAT WIDTH INCREASE FROM 16" TO 22"
 SEAT DEPTH INCREASE FROM 16" TO 24"
 5,6 OR 8 X 1 1/4" HD CASTERS
- PLUS PACKAGE 2:** CLIENT 350 – 449 LBS - \$900
AADL CODE: W102 INCLUDES THE FOLLOWING:
 STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS,
 HEAVY DUTY STRUTS AND SEAT PAN
 SEAT WIDTH INCREASE FROM 18" TO 26"
 SEAT DEPTH INCREASE FROM 16" TO 24"
 5,6 OR 8 X 1 1/4" HD CASTERS
 HD STEEL FORKS
- PLUS PACKAGE 3:** CLIENT 450 – 600 LBS - \$1,110
AADL CODE: W103 INCLUDES THE FOLLOWING:
 STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS,
 HEAVY DUTY STRUTS AND SEAT PAN
 SEAT WIDTH INCREASE FROM 20" TO 30"
 SEAT DEPTH INCREASE FROM 16" TO 24"
 5,6 OR 8 X 1 1/4" HD CASTERS
 HD STEEL FORKS

70 Degree Pivot legrest are recommended for all PLUS packages. If 60 Degree or 80 Degree are required, Pin Style legrest are recommended.

PLEASE NOTE:

- A1 =** Client's Cost
 - A2 =** Clinical Justification
 - A3 =** T-Arm recommended for seat depth & angle back adjustment.
 - A4 =** To be mounted on inside of frame with inside wheels option.
 - A5 =** Included when Heavy duty packages are ordered
- For more Seat To Floor combinations, visit our website
<https://powerplusbility.com/stf>
 Please confirm with customer service (csr@powerplusbility.com)

NOTES:

INFORMATION

Customer Name: Date:
 Address: PO#:
 City: Phone:
 Account:
 Authorizer Name:

BACKRESTS		01
<input type="checkbox"/> N590	Relax Form Fit (No Lateral)	\$690
<input type="checkbox"/> N590	Relax Gentle Fit (<3" Lateral)	\$690
<input type="checkbox"/> N590	Relax Gentle Form Fit (<3" Lateral)	\$690
<input type="checkbox"/> N590	Relax Gentle Spiral Fit (<3" Lateral)	\$690
<input type="checkbox"/> N593	Relax Snug Form Fit (>3" Lateral)	\$790
<input type="checkbox"/> N593	Relax Snug Spiral Fit (>3" Lateral)	\$790
<input type="checkbox"/> N593	Relax Full Form Fit (>3" Lateral)	\$790
<input type="checkbox"/> N593	Relax Full Spiral Fit (>3" Lateral)	\$790
<input type="checkbox"/> N593	Relax Snug Fit (>3" Lateral)	\$790
<input type="checkbox"/> N593	Relax Full Fit (>3" Lateral)	\$790
Dimensions (in inches):		
Width: _____ x Height: _____		
ADDITIONAL OPTIONS		02
<input type="checkbox"/> N500	Heavy-Duty Backrest (>20" Wide)	\$300
<input type="checkbox"/> V571	Replacement Back Cover	\$183
PRESSURE REDUCTION OPTIONS		03
Limited to 3 pieces per backrest.		
<input type="checkbox"/> V853	Lumbar Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V854	Lateral Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V855	Pelvic Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V856	Wedge Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V901	Thoracic Gel Pressure Reduction QTY: ____ Location: ____	\$69
<input type="checkbox"/> V902	Full Spinal Gel Pressure Reduction QTY: ____ Location: ____	\$103.50
<input type="checkbox"/> V903	Full Back Gel Pressure Reduction QTY: ____ Location: ____	\$137.08

HEADREST PAD		07
<input type="checkbox"/> V081	Standard Pad	\$190
<input type="checkbox"/> N582	Headrest with Lateral Support	\$255
HEADREST HARDWARE		08
<input type="checkbox"/> V338	Standard Straight Hardware	\$165
<input type="checkbox"/> N533	Multi-Axis Hardware (L-Bend)	\$200
<input type="checkbox"/> N535	Gooseneck Multi-Axis (S-Bend)	\$280
<input type="checkbox"/> V085	Mounting Plate <input type="checkbox"/> 1/2" or <input type="checkbox"/> 3/4"	\$189.70
<input type="checkbox"/> V570	Replacement Headrest Cover	\$120

ACCESSORIES		06
<input type="checkbox"/> V577	Arm Trough Tray Assembly	\$432
<input type="checkbox"/> V476	Wheelchair Tray (Up to 18" Wide)	\$257.60
<input type="checkbox"/> V542	Wheelchair Tray (Over 18" Wide)	\$310.20
<input type="checkbox"/> V762	Elbow Stops	\$90
<input type="checkbox"/> V766	Lateral Support	\$90
<input type="checkbox"/> V769	Hand Block	\$90
<input type="checkbox"/> N814	Lateral Thigh Support	\$150
<input type="checkbox"/> V619	Modified Foot Support (Padded Footbox)	\$359.10
<input type="checkbox"/> V581	Amputee Board (Up to 20" Wide)	\$350
<input type="checkbox"/> V583	Heavy Duty Amputee Board	\$412.50
<input type="checkbox"/> V579	Pad for Amputee Board	\$62.50
<input type="checkbox"/> V572	Modular Arm Support (Left Side)	\$62.50
<input type="checkbox"/> V572	Modular Arm Support (Right Side)	\$62.50

NOTES:

CUSHIONS		04
<input type="checkbox"/> V980	Relax Comfy (Up to 20"x20")	\$690
<input type="checkbox"/> V980HD	Relax Firm (Over 20"x20")	\$690
<input type="checkbox"/> V980	Relax AIRFLOW (Up to 20"x20")	\$790
<input type="checkbox"/> V980HD	Relax AIRFLOW (Over 20"x20")	\$790
Dimensions (in inches):		
Width: _____ x Depth: _____		
ADDITIONAL OPTIONS		05
<input type="checkbox"/> V101	Replacement Cushion Cover	\$165
<input type="checkbox"/> V006	Incontinence Cover	\$135

