

VGP Voyager Wheelchair - Folding Wheelchair



INFORMATION

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR	Select one below	1	ADP
<input type="checkbox"/> VGPS	Frame - Silver	\$1,690	WM0001063
<input type="checkbox"/> VGPB	Frame - Black	\$1,690	WM0001063
<input type="checkbox"/> VGPC	Frame - Colour : _____	\$2,090	WM0001063

CLIENT WEIGHT	Select one below	2	ADP
<input type="checkbox"/> VGPWA	Up to 250Lbs	N/C	
<input type="checkbox"/> VGPWAMF	Exceeds 250Lbs	\$730	WAMF
<input type="checkbox"/> VGPWAMG	Exceeds 350Lbs	\$992	WAMG
<input type="checkbox"/> VGPWAMH	Exceeds 400Lbs (Max 600Lbs)	\$1,110	WAMH

SEAT WIDTH	Select one below	3	ADP
<input type="checkbox"/> VGPSW14	Seat Width 14"	N/C	
<input type="checkbox"/> VGPSW15	Seat Width 15"	N/C	
<input type="checkbox"/> VGPSW16	Seat Width 16"	N/C	
<input type="checkbox"/> VGPSW17	Seat Width 17"	N/C	
<input type="checkbox"/> VGPSW18	Seat Width 18"	N/C	
<input type="checkbox"/> VGPSW19	Seat Width 19"	\$399	WAMA
<input type="checkbox"/> VGPSW20	Seat Width 20"	\$399	WAMA
<input type="checkbox"/> VGPSWS	Custom Seat Width: _____	\$399	WAMA

SEAT DEPTH	Select one below	4	ADP
<input type="checkbox"/> VGPSD14	Seat Depth 14"	N/C	
<input type="checkbox"/> VGPSD15	Seat Depth 15"	N/C	
<input type="checkbox"/> VGPSD16	Seat Depth 16"	N/C	
<input type="checkbox"/> VGPSD17	Seat Depth 17"	N/C	
<input type="checkbox"/> VGPSD18	Seat Depth 18"	N/C	
<input type="checkbox"/> VGPSD19	Seat Depth 19"	\$265	WAMB
<input type="checkbox"/> VGPSD20	Seat Depth 20"	\$265	WAMB
<input type="checkbox"/> VGPSDS	Custom Seat Depth: _____	\$265	WAMB

BACK HEIGHT	Select one below	5	ADP
<input type="checkbox"/> VGPBH15	Back Height 15"	N/C	
<input type="checkbox"/> VGPBH16	Back Height 16"	N/C	
<input type="checkbox"/> VGPBH17	Back Height 17"	N/C	
<input type="checkbox"/> VGPBH18	Back Height 18"	N/C	
<input type="checkbox"/> VGPBH19	Back Height 19"	N/C	
<input type="checkbox"/> VGPBH20	Back Height 20"	N/C	
<input type="checkbox"/> VGPBH21	Back Height 21"	N/C	
<input type="checkbox"/> VGPBHC	Back Height Custom : _____	\$120	

STANDARD NO CHARGE SPECIFICATIONS
HEAVY DUTY STEEL FRAME AND CROSS BRACES
STANDARD 18" BACK HEIGHT
20", 22", OR 24" NYLON MOULDED REAR WHEELS
EXTRUDED ALUMINUM FORKS
CHROME HAND RIMS
5", 6", OR 8" CASTERS
BLACK NYLON UPHOLSTERY
POWDER COATED FRAME
FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
PRECISION BEARINGS
MOULDED PLASTIC FOOTPLATES
SEAT WIDTH 14" TO 18"
VERTICAL/HORIZONTAL REAR AXLE PLATE ADJUSTABILITY
ANTI-TIP ROLLERS
6" BRAKE EXTENSIONS
DOUBLE CROSS BRACE
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	12.75"	20"	5"	2H-N
<input type="checkbox"/>	13"	20"	5"	3H,2H,1H
<input type="checkbox"/>	13.75"	20"	5"	3H,2H
<input type="checkbox"/>	14.25"	20"	6"	5H,3H,2H
<input type="checkbox"/>	14.5"	22"	5"	5H,3H
<input type="checkbox"/>	14.25"	22"	6"	5H,3H,2H
<input type="checkbox"/>	15.25"	20"	5"	5H
<input type="checkbox"/>	15"	20"	6"	5H,3H
<input type="checkbox"/>	15"	22"	5"	5H
<input type="checkbox"/>	15"	22"	6"	5H,3H
<input type="checkbox"/>	16.25"	22"	6"	5H
<input type="checkbox"/>	16.5"	22"	8"	5H
<input type="checkbox"/>	17.25"	22"	8"	5H

(HEIGHT EXTENDED CASTER)				
	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	18.25"	24"	8"	5H,3H

For more Seat To Floor combinations, visit our website

Notes:

FOOT REST AND FOOT PLATES		Select one below	5	ADP	SEATING - BACKS		Custom: select back & fill custom	12	ADP
<input type="checkbox"/>	VGPF6A	60 Degree Hanger Pivot	N/C		<input type="checkbox"/>	VGPBK01	RELAX Gentle Fit Back	\$544	SEBWM0020
<input type="checkbox"/>	VGPF7A	70 Degree Hanger Pivot	N/C		<input type="checkbox"/>	VGPBK02	RELAX Full Fit Back	\$608	SEBWM0010
<input type="checkbox"/>	VGPF8A	80 Degree Hanger Pivot	N/C		<input type="checkbox"/>	VGPBK03	RELAX Snug Fit Back	\$700	SEBWM0027
<input type="checkbox"/>	VGPEA	Elevating Legrest Pivot	\$404	WAMK	<input type="checkbox"/>	VGPBK04	RELAX Form Fit Back No Laterals	\$820	SEBWM0026
<input type="checkbox"/>	VGPAL	Articulating Leg Rests Pivot	\$450	WAMK	<input type="checkbox"/>	VGPBK06	RELAX Form Fit Back Gentle Laterals	\$820	SEBWM0026
<input type="checkbox"/>	VGPPPL	Plastic	N/C		<input type="checkbox"/>	VGPBK07	RELAX Form Fit Back Full Laterals	\$820	SEBWM0026
<input type="checkbox"/>	VGPPXL	Extra Large Foot Rest	\$150	SE0000322	<input type="checkbox"/>	VGPBK05	RELAX Form Fit Back Snug Laterals	\$820	SEBWM0026
<input type="checkbox"/>	VGPPAF	Angle Adjustable Footplates	\$240	WAMJ	<input type="checkbox"/>	VGPBK10	RELAX Spiral Fit Back Gentle Laterals	\$832	SEBWM0030
<input type="checkbox"/>	VGPIFB	One Piece Foot Board	\$80	SEFND0015	<input type="checkbox"/>	VGPBK09	RELAX Spiral Fit Back Snug Laterals	\$832	SEBWM0030
<input type="checkbox"/>	VGPIFBP	One Piece Foot Board Padding	\$44	SEFND0035	<input type="checkbox"/>	VGPBK11	RELAX Spiral Fit Back Full Laterals	\$832	SEBWM0030
<input type="checkbox"/>	VGPNLS	Nylon Leg Strap Padded	\$60	SEFND0010	<input type="checkbox"/>	VGPTAB	RELAX Tension Back Upto 18"	\$140	WAMD
<input type="checkbox"/>	VGPNLSC	Nylon Leg Strap Custom Size ()	\$120	SEFCF107L	<input type="checkbox"/>	VGPTAB+	RELAX Tension Back Over 18"	\$240	WAME
<input type="checkbox"/>	VGPCSL	Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style	\$240	SEFND0020	<input type="checkbox"/>	VGPBK12	* Custom Width : _____	\$176	SEBWM0035
BACK TYPE		Select one below	6	ADP	<input type="checkbox"/>	VGPBK13	* Custom Height : _____	\$176	SEBWM0040
<input type="checkbox"/>	VGAB1	Angle Adjustable Back	N/C		*Note: Please refer to the Seating Booklet for custom sizes				
ARM REST		Select one below	8	ADP	<input type="checkbox"/>	VGPBGL1	Gel Full Overlay	\$348	SEBND0020 X 3
(Ht. 10"-14")	<input type="checkbox"/>	VGPF6A	Flip Back Ht. Adj./ Removable	\$340	<input type="checkbox"/>	VGPBGL2	Gel Center Vertebrae	\$116	SEBND0020
(Ht. 10"-14")	<input type="checkbox"/>	VGPTA	T-Arm Ht. Adj./ Removable	\$240	<input type="checkbox"/>	VGPBGL3	Gel Left Lateral	\$116	SEBND0020
(Ht. 7"-11")	<input type="checkbox"/>	VGPFAS	Flip Back Ht. Adj./ Removable Short	N/C	<input type="checkbox"/>	VGPBGL4	Gel Right Lateral	\$116	SEBND0020
(Ht. 7"-11")	<input type="checkbox"/>	VGPTAS	T-Arm Ht. Adj./ Removable Short	\$340	<div style="border: 1px solid red; padding: 2px; display: inline-block;"> Maximum 1 piece additional pressure relief gel </div>				
<input type="checkbox"/>	VGPGWAPS	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030					
* Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required									
REAR WHEEL		Select one below	9	ADP	SEATING - CUSHIONS		Select one below	13	ADP
<input type="checkbox"/>	VGPRWT0	20" Solid Treaded Tire	N/C		<input type="checkbox"/>	VGPC01	RELAX Soft Cushion	\$480	SESWM0050
<input type="checkbox"/>	VGPRWT2	22" Solid Treaded Tire	N/C		<input type="checkbox"/>	VGPC02	RELAX Comfy Cushion	\$624	SESWM0030
<input type="checkbox"/>	VGPRWT4	24" Solid Treaded Tire	N/C		<input type="checkbox"/>	VGPC03	RELAX Firm Cushion	\$784	SESWM0040
<input type="checkbox"/>	VGPRWA2	22" Air Tire	\$120		<input type="checkbox"/>	VGPC04	RELAX AIR FLOW Cushion	\$587	SE0001669
<input type="checkbox"/>	VGPRWA4	24" Air Tire	\$120		<input type="checkbox"/>	VGPCGL2	Full Overlay Gel (2 Pieces)	\$258	SESND1045
<input type="checkbox"/>	VGPAQR	Axles – Quick Release	\$120	WAMU	<input type="checkbox"/>	VGPCCL	Contoured Leg Channel (Bolster) Std.	\$180	SEMCF203L
CASTER		Select one below	10	ADP	<input type="checkbox"/>	VGPCCLA	Contoured Leg Channel (Bolster) with ABS	\$260	SESDND1035 SESDND1015
<input type="checkbox"/>	VGPF1A	5" X 1.25" Caster	N/C		<input type="checkbox"/>	VGPCCLC	* Contoured Leg Channel (Bolster) Custom	\$320	SEMCF209L
<input type="checkbox"/>	VGPF1B	6" X 1.25" Caster	N/C		*Note: Please refer to the Seating Booklet for custom sizes				
<input type="checkbox"/>	VGPF1C	8" X 1.25" Caster	N/C		<input type="checkbox"/>	VGPC00	Extra Cushion Cover	\$99	SESDND1055
<input type="checkbox"/>	VGPF2A	5" X 1.75" Caster (Soft Roll)	\$80		ACCESSORIES			14	ADP
<input type="checkbox"/>	VGPF2B	6" X 1.75" Caster	\$80		<input type="checkbox"/>	VGPSBA	Seat Belt (Auto Style)	\$26	
<input type="checkbox"/>	VGPF2C	8" X 1.75" Caster	\$80		<input type="checkbox"/>	VGPSBAP	Seat Belt (Padded)	\$56	SERND0030
<input type="checkbox"/>	VGPF4B	6" X 2 Caster	\$80		<input type="checkbox"/>	VGPSBA4P	Seat Belt (4 Point Padded)	\$130	SERND0020
<input type="checkbox"/>	VGPF4C	8" X 2 Caster	\$80		<input type="checkbox"/>	VGPTRA	Full Padded Tray Locking	\$280	SETND0005 SETND0035 SETND0005 SETND0035
WHEEL LOCK		Select one below	11	ADP	<input type="checkbox"/>	VGPTRB	Full Padded Tray Slide On	\$280	SETND0035
<input type="checkbox"/>	VGPALPH	Aluminium Push To Lock	N/C		<input type="checkbox"/>	VGPTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	VGPALPL	Aluminium Pull To Lock	N/C		<input type="checkbox"/>	VGPTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	VGPOSTPH	Steel Push To Lock	N/C		<input type="checkbox"/>	VGPTRD	Half Padded Left Tray	\$240	SETND0015 SETND0035
<input type="checkbox"/>	VGPOSTPL	Steel Pull To Lock	N/C		<input type="checkbox"/>	VGPTRE	Half Padded Right Tray	\$240	SETND0015 SETND0035
<input type="checkbox"/>	VGPRXC	4" Brake Extension	N/C		<input type="checkbox"/>	VGPTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>	VGPRXA	6" Brake Extension	N/C		<input type="checkbox"/>	VGPTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>	VGPRXB	9" Brake Extension	\$130		<input type="checkbox"/>	VGPTDFB	Transport Tie Downs Front	\$85	
					<input type="checkbox"/>	VGPTDRB	Transport Tie Downs Rear	\$85	
					<input type="checkbox"/>	VGPOXY	Oxygen Tank Holder	\$195	WEPX
					<input type="checkbox"/>	VGPIVA	Adjustable I V Pole	\$130	
					<input type="checkbox"/>	VGPCHT-L	Left Arm Trough Tray	\$400	SETND0015 SETND0035
					<input type="checkbox"/>	VGPCHT-R	Right Arm Trough Tray	\$400	SETND0030
					<input type="checkbox"/>	VGPLATB	Lateral Bite out	\$180	SELND0015
					<input type="checkbox"/>	VGPAADD8	8" Adductor	\$164	SELND0005
					<input type="checkbox"/>	VGPAADD10	10" Adductor	\$164	SELND0005
					<input type="checkbox"/>	VGPAADD12	12" Adductor	\$164	SELND0005
					<input type="checkbox"/>	VGPFRC	Padded Footplates	\$88	SEMCF803L