

# MAGIC PLUS


FOLDING WHEELCHAIR



**ADP CODE:** WM0001062

The MAGIC PLUS is a folding chair designed with an aluminum light weight composition for easy maneuverability. This chair also offers one of the lowest seat to floor height settings available.





Vertical Rear Axle Adjustability  
Ranging From 12.5" To 18.25"



Angle Adjustable Back Canes 0°, 5°,  
10° Or 15°




Dual cross braces





Two Point Flip Back Arm



Height Adjustable Arm



Swing In And Swing Out Front Rigging

## CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# MAGIC PLUS



## STANDARD FEATURES

- Light Weight Aluminum Frame & Cross Braces
- Multi Axle Position - Vertical
- 20", 22" or 24" Urethane Molded Rear Wheels
- New Light Weight 20", 22" or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- T-Arm Height Adjustable/ Removable
- New Flip Back Arms (Allows Angle Back Adjustment)
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 20"
- Black Padded Arm Pads
- Anti-Tip Rollers
- 6" Brake Extensions
- Overall Wheelchair Width: 27" (18" Seat Width Unfolded)  
18" (18" Seat Width folded)

SPECIFICATIONS	MAGIC PLUS
SEAT WIDTH:	14" TO 20"
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.5" TO 18.25"
REAR WHEEL:	20", 22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	5", 6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWING AWAY, DUAL SWING AWAY OR ELEVATING LEG RESTS
BACK ANGLE (ADJUSTABLE):	0°, 5°, 10° OR 15°
BACK HEIGHT RANGE:	15" TO 24"
WEIGHT CAPACITY:	250 lbs. (Up to 600 lbs. With HD package)
WEIGHT OF CHAIR:	31 lbs.

MAGIC PLUS  
FOLDING WHEELCHAIR



Learn more about this product by scanning the QR Code

