

RELAX SERIES

RELAX

HEADREST


<p>MAHS-F Dia.1</p>  <p>Headrest with Contoured Pad Large pad with angled edges.</p>	<p>MAHS-KP Dia.2</p>  <p>Headrest with Medium sized pad with similar features as the contoured pad.</p>	<p>MAHS-FT Dia.3</p>  <p>Headrest with Flat Pad Small pad with a slight curve at the two ends.</p>
<p>MAHS-FXL Dia.4</p>  <p>Headrest with Flex Pad Large pad that can be flexed to form a custom position.</p>	<p>MAHS-FXS Dia.5</p>  <p>Headrest with Flex Pad Medium pad that can be flexed to form a custom position.</p>	<p>MAHS-PL Dia.6</p>  <p>Headrest with Poly Urethane Pad Large pad that is injection molded and water resistant.</p>
<p>MAHS-C Dia.7</p>  <p>Headrest with Square Pad Large pad that has thicker edges and narrow in the center.</p>	<p>UHRA Dia.8</p>  <p>Universal Headrest Plate - 1/2"</p>	<p>UHRA3-4 Dia.9</p>  <p>Universal Headrest Plate - 3/4"</p>

Quick Reference Guide

Accessories - Headrests and Hardware

Part Number	Description	Dia	Product ADP Code
MAHS-F	Muti-Axis Headrest S Bend with Large Contoured Pad	1	SE0001542
MAHS-KP	Muti-Axis Headrest S Bend with Medium Contoured Pad	2	SE0001542
MAHS-FT	Muti-Axis Headrest S Bend with Small Flat Pad	3	SE0001542
MAHS-FXL	Muti-Axis Headrest S Bend with Large Flex Pad	4	SE0001542
MAHS-FXS	Muti-Axis Headrest S Bend with Medium Flex Pad	5	SE0001542
MAHS-PL	Muti-Axis Headrest S Bend with Large Poly Urethane Pad	6	SE0001542
MAHS-C	Muti-Axis Headrest S Bend with Large Square Pad	7	SE0001542
UHRA	Universal Headrest Plate - 1/2"	8	SEHND0020
UHRA3-4	Universal Headrest Plate - 3/4"	9	SEHND0020

Headrest Neck Chart



Headrest Pad Chart

PAD CODE	DESCRIPTION	INSIDE MEASUREMENT (W X H X T)	OUTSIDE MEASUREMENT (W X H X T)
C	Large Square Pad	12" X 8" X 2.25" (3.5" edges)	12" X 8" X 2.25" (3.5" edges)
F	Large Contoured Pad	14" X 6" X 1.5"	17" X 6" X 1.5"
KP	Medium Contoured Pad	11" X 6" X 1.5"	13" X 6" X 1.5"
FT	Small Flat Pad	10" X 6" X 1.5"	11" X 6" X 1.5"
FXL	Large Flex Pad	16" X 6" X 1.5"	16" X 6" X 1.5"
FXS	Small Flex Pad	11" X 6" X 1.5"	13" X 6" X 1.5"
PL	Large Poly Urethane Pad	12.5" X 5" X 1.5"	16" X 5" X 1.5"

ORDERING A HEADREST SET

Step 1: Choose a Neck from the Headrest Neck Chart Above

Step 2: Choose a Pad from the Headrest Pad Chart Above

Step 3: Add the Headrest Neck Code and the Pad Code together when submitting an order.