

VGP VOYAGER Wheelchair (Manual Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR *Select one below* **1**

- VOYAGER PLUS FRAME - Black 1450

SEAT WIDTH *Select one below* **2**

- SW14"-20" Specify _____ N/C
(Sw14" - 20" Light Weight Aluminium Frame)
- SW21"-28" Specify _____
(Sw21" -28" Must Include Plus Package 1, 2 Or 3)

SEAT DEPTH *Select one below* **3**

- SDF18" (By Frame) N/C
- SD16" -18" Specify _____ (By Upholstery) N/C
- SDF20" (By Frame) 350
- SD18" -20" Specify _____ (By Upholstery) N/C
- SDF22" (By Frame) 350
- SD20" -22" Specify _____ (By Upholstery) N/C

BACK *Select one below* **4**

- BH15" - 20" Specify _____ N/C
- BH21" + Specify _____ 120
- AB Angle Adj. Back N/C
- TAB Tension Adj Back Uph. (250Lbs) 180
- ENDO* Endomorph Canes **W ()** N/C

Specify Width

Endomorph Canes require TA or TAS Arm Rest

STROLLER HANDLE *Select one below* **5**

- SHS Removable Stroller Handle N/C
- SHR Removable Raised Stroller Handle N/C

FOOT REST *Select one below* **6**

- F6 60 Degree Hanger N/C
- F7 70 Degree Hanger N/C
- F8 80 Degree Hanger N/C
- EA Elevating Legrest 404
- AL Articulating Legrest 450

FOOT REST AND FOOT PLATES *Select one below* **7**

- PPL Plastic N/C
- PAF Angle Adjustable - Extra Large 150
- PXL Extra Large – Alum (F6A,F7A,F8A) 150
- 1FB One Piece Foot Board 80

SPECIFICATIONS

LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES
STANDARD 18" BACK HEIGHT
20", 22", OR 24" NYLON MOULDED REAR WHEELS
EXTRUDED ALUMINUM FORKS
CHROME HAND RIMS
5", 6", OR 8" CASTERS
BLACK NYLON UPHOLSTERY
POWDER COATED FRAME
FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
PRECISION BEARINGS
MOULDED PLASTIC FOOTPLATES
SEAT WIDTH 14" TO 18"
BLACK PADDED ARM PADS
ANTI-TIP ROLLERS
6" BRAKE EXTENSIONS
DOUBLE CROSS BRACE

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	

(STANDARD SEAT SETTING)

<input type="checkbox"/> MH 16.75	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	8"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	24"	8"	LONG
<input type="checkbox"/> TH18.25"	24"	8"	LONG
<input type="checkbox"/> TT 19.00"*	24"	8"	LONG

**TT 19.00' requires taller caster journal*

ARM REST <i>Select one below</i>		8	SEATING - BACKS <i>Select one below</i>		13
<input type="checkbox"/> FA	Flip Back Ht. Adj./ Removable	240		W () x H ()	
<input type="checkbox"/> TA*	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/> ADLBK01	RELAX Gentle Fit Flat Back	350
<input type="checkbox"/> FAS	Flip Back (Short) Ht. Adj./ Remove	240			
<input type="checkbox"/> TAS*	T-Arm (Short) Ht. Adj./ Removable	240	SEATING - CUSHION <i>Select one below</i>		14
<input type="checkbox"/> GWAPS	Gel Waterfall Arm Pads	150		W () x D ()	
<i>*14" - 20" Wide Light Weight Aluminum Frame Only</i>			<input type="checkbox"/> ADLC01	RELAX Soft Cushion	350
REAR WHEELS <i>Select one below</i>		9	PLUS PACKAGES <i>Select one below</i>		15
<input type="checkbox"/> RWT0	20" Solid Treaded Tire	N/C	<input type="checkbox"/> PLUS PACKAGE 1:		
<input type="checkbox"/> RWT2	22" Solid Treaded Tire	N/C	CLIENT 250 TO 349 LBS - \$730.00		
<input type="checkbox"/> RWT4	24" Solid Treaded Tire	N/C	INCLUDES THE FOLLOWING:		
<input type="checkbox"/> RWS0*	20" Spoked Treaded Tire	N/C	STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY		
<input type="checkbox"/> RWS2*	22" Spoked Treaded Tire	N/C	SEAT WIDTH INCREASE FROM 16" TO 24"		
<input type="checkbox"/> RWS4*	24" Spokedtreaded Tire	N/C	SEAT DEPTH INCREASE FROM 16" TO 22"		
<input type="checkbox"/> RWP2	22" Pyramid Tire	80	5,6 OR 8 X 1 3/4" HD CASTERS		
<input type="checkbox"/> RWP4	24" Pyramid Tire	80	WELDED CASTER JOURNAL		
<input type="checkbox"/> RWA2	22" Air Tire	120	<input type="checkbox"/> PLUS PACKAGE 2:		
<input type="checkbox"/> RWA4	24" Air Tire	120	CLIENT 350 TO 449 LBS - \$992.00		
<input type="checkbox"/> AQR	Axles – Quick Release (349 Lbs)	N/C	INCLUDES THE FOLLOWING:		
<input type="checkbox"/> PCHR	Plastic Coated Hand Rims	N/C	STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY		
<i>*14" - 20" Wide Light Weight Aluminum Frame Only</i>			SEAT WIDTH INCREASE FROM 18" TO 26"		
CASTER <i>Select one below</i>		10	SEAT DEPTH INCREASE FROM 16" TO 22"		
<input type="checkbox"/> FC1A	5" X 1.25" Caster	N/C	5,6 OR 8 X 1 3/4" HD CASTERS		
<input type="checkbox"/> FC1B	6" X 1.25" Caster	N/C	HD STEEL FORKS		
<input type="checkbox"/> FC1C	8" X 1.25" Caster	N/C	WELDED CASTER JOURNAL		
<input type="checkbox"/> FC2A	5" X 1.75" Caster	80	<input type="checkbox"/> PLUS PACKAGE 3:		
<input type="checkbox"/> FC2B	6" X 1.75" Caster	80	CLIENT 450 TO 600 LBS - \$1110.00		
<input type="checkbox"/> FC2C	8" X 1.75" Caster	80	INCLUDES THE FOLLOWING:		
<input type="checkbox"/> FC3A	5" X 1.25" Air Caster	80	STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY		
<input type="checkbox"/> FC3B	6" X 1.25" Air Caster	80	SEAT WIDTH INCREASE FROM 20" TO 28"		
<input type="checkbox"/> FC3C	8" X 1.25" Air Caster	80	SEAT DEPTH INCREASE FROM 16" TO 22"		
WHEEL LOCK <i>Select one below</i>		11	5,6 OR 8 X 1 3/4" HD CASTERS		
<input type="checkbox"/> ALPL	Aluminium Pull To Lock	N/C	HD STEEL FORKS		
<input type="checkbox"/> ALPH	Aluminium Push To Lock	N/C	WELDED CASTER JOURNAL		
<input type="checkbox"/> STPL	Steel Pull To Lock	N/C	<input type="checkbox"/> FIXED REINFORCED ONE PIECE BACK POST \$260.00		
<input type="checkbox"/> STPH	Steel Push To Lock	N/C	(SPECIFY ANGLE _____)		
<input type="checkbox"/> BRXA	6" Brake Extension	N/C	<input type="checkbox"/> HEAVY DUTY TENSION ADJ. BACK UPH. \$220.00		
<input type="checkbox"/> BRXB	9" Brake Extension	60			
ACCESSORIES		12	Notes:		
<input type="checkbox"/> SBA	Seat Belt (Auto Style)	26			
<input type="checkbox"/> SBB	Seat Belt (Velcro Style)	26			
<input type="checkbox"/> SBC	Seat Belt (Plastic Style)	26			
<input type="checkbox"/> CNC	Crutch And Cane Holder (A1)	60			
<input type="checkbox"/> IVA	Adjustable I V Pole (A1)	130			
<input type="checkbox"/> NLSP	Nylon Leg Strap Padded (A1)	50			
<i>Please Note : A1 = Clients Responsibility</i>					