

**NVT NVT Wheelchair (Tilt Wheelchair)**



**POWER PLUS MOBILITY**

*We will truly move you*

Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



**FRAME COLOUR** *Select one below* **1**

- NVTB Tilt Frame - Black 3299
- NVTS Tilt Frame - Silver 3299
- NVTC Tilt Frame - Custom : \_\_\_\_\_ 3699

**CLIENT WEIGHT** *Select one below* **2**

- NVTWA Up to 250Lbs N/C
- NVTWAMF Exceeds 250Lbs 730
- NVTWAMG Exceeds 350Lbs 992
- NVTWAMH Exceeds 400Lbs (Max 600Lbs) 1110

**SEAT WIDTH** *Select one below* **3**

- NVTSW16 Seat Width 16" N/C
- NVTSW17 Seat Width 17" N/C
- NVTSW18 Seat Width 18" N/C
- NVTSW19 Seat Width 19" 399
- NVTSW20 Seat Width 20" 399
- NVTSWC Custom Seat Width: \_\_\_\_\_ 399

**SEAT DEPTH** *Select one below* **4**

- NVTSD16 Seat Depth 16" N/C
- NVTSD17 Seat Depth 17" N/C
- NVTSD18 Seat Depth 18" N/C
- NVTSD19 Seat Depth 19" 265
- NVTSD20 Seat Depth 20" 265
- NVTSDC Custom Seat Depth: \_\_\_\_\_ 265

**FOOT REST AND FOOT PLATES** *Select one below* **5**

- NVTOFST Offset Foot Rest N/C
- NVTF6A 60 Degree Hanger  Pivot  Pin N/C
- NVTF7A 70 Degree Hanger  Pivot  Pin N/C
- NVTF8A 80 Degree Hanger  Pivot  Pin N/C
- NVTEA Elevating Legrest  Pivot  Pin 404
- NVTAL Articulating L/R  Pivot  Pin 450
- NVTEB Elevating Legrest System N/C
- NVTPPL Plastic N/C
- NVTPXL Extra Large 150
- NVTPAF Angle Adjustable Footplates 240
- NVTFRCH Padded Footplates 88
- NVT1FB One Piece Foot Board 80
- NVTNLS Nylon Leg Strap Padded 60
- NVTNLS Custom Size ( ) 120
- NVTCSL Contracture Sling 240

**STANDARD NO CHARGE SPECIFICATIONS**

SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS

**Tilt Options** **6**

- NVTINTS Independent Power Tilt System 2000

**SEAT TO FLOOR SPECIFICATIONS**

**(LOW SEAT SETTING)**

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	

**(STANDARD SEAT SETTING)**

<input type="checkbox"/> MH 16.75"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	8"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	24"	8"	LONG
<input type="checkbox"/> A 18.25"	24"	8"	LONG

**Notes:**

BACK TYPE		7	SEATING - BACKS		For custom, choose one back then fill custom	13
<input type="checkbox"/>	NVTAB1	Angle Adjustable Back	N/C	<input type="checkbox"/>	NVTBK01	<b>RELAX</b> Gentle Fit Back 544
<i>Specify Width</i> <input type="checkbox"/>	NVTENDO	Endomorph Canes W ( )	N/C	<input type="checkbox"/>	NVTBK02	<b>RELAX</b> Full Fit Back 608
<input type="checkbox"/>	NVTRBK	Recline Kit	N/C	<input type="checkbox"/>	NVTBK03	<b>RELAX</b> Snug Fit Back 700
HEADREST		Select one below	8	<input type="checkbox"/>	NVTBK04	<b>RELAX</b> Form Fit Back No Laterals 820
<input type="checkbox"/>	NVMAHF	Multi-Axis Hardware Complete with Large Pad	512	<input type="checkbox"/>	NVTBK06	<b>RELAX</b> Form Fit Back Gentle Laterals 820
<input type="checkbox"/>	NVMAHFXL	Multi-Axis Hardware Complete with Flex Large	512	<input type="checkbox"/>	NVTBK07	<b>RELAX</b> Form Fit Back Full Laterals 820
<input type="checkbox"/>	NVMAHKP	Multi-Axis Hardware Complete with Medium Pad	512	<input type="checkbox"/>	NVTBK05	<b>RELAX</b> Form Fit Back Snug Laterals 820
<input type="checkbox"/>	NVMAHFXS	Multi-Axis Hardware Complete with Flex Medium	512	<input type="checkbox"/>	NVTBK10	<b>RELAX</b> Spiral Fit Back Gentle Laterals 832
<input type="checkbox"/>	NVMAHFT	Multi-Axis Hardware Complete with Small Flat	512	<input type="checkbox"/>	NVTBK09	<b>RELAX</b> Spiral Fit Back Snug Laterals 832
<i>(1/2" Tubing)</i> <input type="checkbox"/>	NVUPHT219	Contour Headrest Complete with Swivel Hardware	433	<input type="checkbox"/>	NVTBK11	<b>RELAX</b> Spiral Fit Back Full Laterals 832
<i>(3/4" Tubing)</i> <input type="checkbox"/>	NVUPHT218	Contour Headrest with Std. Assembly (1/2")	292	<input type="checkbox"/>	NVTTAB	<b>RELAX</b> Tension Back Upto 18" 140
<input type="checkbox"/>	NVUPHT217	Contour Headrest with Std. Assembly (3/4")	292	<input type="checkbox"/>	NVTTAB+	<b>RELAX</b> Tension Back Over 18" 240
<input type="checkbox"/>	NVUHA	Universal Hardware Adaptor (Adaptor only)	57	<input type="checkbox"/>	NVTBGL1	Half Overlay Gel (3 Piece) 348
<input type="checkbox"/>	NVMAHGA-P	Multi Axis with Headrest Plate (No Pad)	255	<input type="checkbox"/>	NVTBGL2	Full Overlay Gel (6 Pieces) 696
ARM REST		Select one below	9	<input type="checkbox"/>	NVTNPVT	Pin Extension (1" Extension) 40
<i>(Std 10"-14")</i> <input type="checkbox"/>	NVTTA	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	NVTBK12	*Custom Width : _____ +176
<i>(Std 10"-14")</i> <input type="checkbox"/>	NVTFA	Flip Back Ht. Adj./ Removable	340	<input type="checkbox"/>	NVTBK13	*Custom Height : _____ +176
<i>(Short 7"-11")</i> <input type="checkbox"/>	NVTTAS	T-Arm Ht. Adj./ Removable	240	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>		
<i>(Short 7"-11")</i> <input type="checkbox"/>	NVTFAS	Flip Back Ht. Adj./ Removable	340	SEATING - CUSHION		Select one below
<input type="checkbox"/>	NVTGWAPF	Gel Waterfall Arm Pads	180			14
REAR WHEELS		Select one below	10			
<input type="checkbox"/>	NVTRWTO	20" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTC01	<b>RELAX</b> Soft Cushion 480
<input type="checkbox"/>	NVTRWT2	22" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTC02	<b>RELAX</b> Comfy Cushion 624
<input type="checkbox"/>	NVTRWT4	24" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTC03	<b>RELAX</b> Firm Cushion 784
<i>(STF 13-19")</i> <input type="checkbox"/>	NVTRWT12	12" Solid Tire	N/C	<input type="checkbox"/>	NVTCGL1	Half Overlay Gel (1 Piece) 116
<i>(STF 17-19")</i> <input type="checkbox"/>	NVTRWT12I	12" Solid Tire Inside Mount	120	<input type="checkbox"/>	NVTCGL2	Full Overlay Gel (2 Pieces) 232
<input type="checkbox"/>	NVTRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	NVTCLL	Contoured Leg Channel (Bolster) Std. 180
<input type="checkbox"/>	NVTRWA2	22" Air Tire	120	<input type="checkbox"/>	NVCCLA	Contoured Leg Channel (Bolster) with ABS 260
<input type="checkbox"/>	NVTRWA4	24" Air Tire	120	<input type="checkbox"/>	NVTCLC	Contoured Leg Channel (Bolster) Custom 320
<input type="checkbox"/>	NVTAQR	Axles – Quick Release	120	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>		
<input type="checkbox"/>	NVTHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	NVTC00	Extra Cushion Cover 99
CASTER		Select one below	11	ACCESSORIES		15
<input type="checkbox"/>	NVTFC1A	5" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTSBA	Seat Belt (Auto Style) 26
<input type="checkbox"/>	NVTFC1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTSBAP	Seat Belt (Padded) 56
<input type="checkbox"/>	NVTFC1C	8" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTSBA4P	Seat Belt (4 Point Padded) 130
<input type="checkbox"/>	NVTFC2A	5" X 1.75" Caster (Soft Roll)	80	<input type="checkbox"/>	NVTTRA	Full Padded Tray Locking 280
<input type="checkbox"/>	NVTFC2B	6" X 1.75" Caster	80	<input type="checkbox"/>	NVTTRB	Full Padded Tray Slide On 280
<input type="checkbox"/>	NVTFC2C	8" X 1.75" Caster	80	<input type="checkbox"/>	NVTTRC	Full Plastic Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 280
<input type="checkbox"/>	NVTFC4B	6" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	NVTTRD	Half Padded Left Tray 240
<input type="checkbox"/>	NVTFC4C	8" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	NVTTRE	Half Padded Right Tray 240
WHEEL LOCK		Select one below	12	<input type="checkbox"/>	NVTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	NVTALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	NVTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	NVTALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	NVTDFB	Transport Tie Downs Front 85
<input type="checkbox"/>	NVTNVTH	Steel Push To Lock	N/C	<input type="checkbox"/>	NVTDRB	Transport Tie Downs Rear 85
<input type="checkbox"/>	NVTNVTL	Steel Pull To Lock	N/C	<input type="checkbox"/>	NVTOXY	Oxygen Tank Holder 195
<input type="checkbox"/>	NVTBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	NVTIVA	Adjustable I V Pole 130
<input type="checkbox"/>	NVTBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	NVTCHPT-L	Left Arm Trough Tray 400
<input type="checkbox"/>	NVTBRXB	9" Brake Extension	130	<input type="checkbox"/>	NVTCHPT-R	Right Arm Trough Tray 400
<input type="checkbox"/>	NVTATTBRX	Attendant Brakes	240	<input type="checkbox"/>	NVTLATB	Lateral Bite out 180
				<input type="checkbox"/>	NVTADD8	8" Adductor 164
				<input type="checkbox"/>	NVTADD10	10" Adductor 164
				<input type="checkbox"/>	NVTADD12	12" Adductor 164