


MULTI-AXIS HEADREST SYSTEM

Toolless multi axis adjustment



MAHS-F	MAHS-KP	MAHS-FT	MAHS-FXL	MAHS-FXS	MAHS-PL	MAHS-C
Headrest with Contoured Pad Large pad with angled edges.	Headrest with KP Pad Medium sized pad with similar features as the contoured pad.	Headrest with Flat Pad Small pad with a slight curve at the two ends.	Headrest with Flex Pad Large pad that can be flexed to form a custom position.	Headrest with Flex Pad Medium pad that can be flexed to form a custom position.	Headrest with Poly Urethane Pad Large pad that is injection molded and water resistant.	Headrest with Square Pad Large pad that has thicker edges and narrow in the center.

FLEX HEADREST PAD



Malleable metal

Easily bend the headrest pad inwards from both ends to get a custom firm fit or simply bend the pad outwards for a wider pad.

RELAX Multi Axis Headrests provide the best support for your head and neck when using a wheelchair that tilts. This allows for additional support to hold your head up while using your wheelchair. The **RELAX** Multi Axis Headrest comes with different options for headrest pad as well as different types of materials.


We customize anything!

Power Plus Mobility own advance machinery specializes in manufacturing custom parts for customer's specification from either a drawing or a sample part. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities.

Multi-Axis Headrest ADP: **SE0001542**

HEADREST REFERENCE GUIDE

Types of headrest neck:



Types of headrest pad:

HEADREST PAD SIZING CHART				
PAD CODE (SIZE)		DESCRIPTION	INSIDE MEASUREMENT (WIDTH X HEIGHT X THICKNESS)	OUTSIDE MEASUREMENT (WIDTH X HEIGHT X THICKNESS)
C	LARGE	SQUARE PAD	12"X8"X2.25" (3.5" EDGES)	12"X8"X2.25" (3.5" EDGES)
F	LARGE	CONTOURED PAD	14" X 6" X 1.5"	17" X 6" X 1.5"
KP	MEDIUM	CONTOURED PAD	11" X 6" X 1.5"	13" X 6" X 1.5"
FT	SMALL	FLAT PAD	10" X 6" X 1.5"	11" X 6" X 1.5"
FXL	LARGE	FLEX PAD	16" X 6" X 1.5"	16" X 6" X 1.5"
FXS	SMALL	FLEX PAD	11" X 6" X 1.5"	13" X 6" X 1.5"
PL	LARGE	POLY URETHANE PAD	12.5"X5"X1.5"	16"X5"X1.5"

Types of headrest mount:



How to order a headrest system?

To order a headrest with a specific pad, please refer to the neck chart and the pad chart above. Choose the part number of the neck required and part number for pad as required.
e.g. **MAH-FXL (NECK CODE) - (PAD CODE)**

