

DU RELAX Full Fit Back - Seating



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



STANDARD SIZES	PRICING	COVER OPTIONS	PRICING
<input type="checkbox"/> DU14X18 Full Fit Back W 14" 18" H	608	<input type="checkbox"/> DUCOV Standard	N/C
<input type="checkbox"/> DU15X18 Full Fit Back W 15" 18" H	608	<input type="checkbox"/> DUCT4 TEK4 Infection Control Cover	N/C
<input type="checkbox"/> DU16X18 Full Fit Back W 16" 18" H	608		
<input type="checkbox"/> DU17X18 Full Fit Back W 17" 18" H	608		
<input type="checkbox"/> DU18X18 Full Fit Back W 18" 18" H	608		
<input type="checkbox"/> DU20X18 Full Fit Back W 20" 18" H	608		

CUSTOM SIZES	Choose one below, then enter size	PRICING
<input type="checkbox"/> DUC	Full Fit Back Custom	608
<input type="checkbox"/> DUCW	Width: _____	+176
<input type="checkbox"/> DUCH	Height: _____	+176

EXTRA COVER	PRICING
<input type="checkbox"/> DUCC Standard Cover	99
<input type="checkbox"/> DUOT4 TEK4 Infection Control Cover	99
<input type="checkbox"/> DUCD Custom Cover <input type="checkbox"/> Standard <input type="checkbox"/> TEK4	134
<input type="checkbox"/> DUCDW Width: _____	
<input type="checkbox"/> DUCDH Height: _____	

*Note: Please enter the size above for standard or custom cover


ADDITIONAL OVERLAYS	PRICING
<input type="checkbox"/> DUOLF Foam Pad	116
<input type="checkbox"/> DUOLV1 Visco Overlay Half (3 Pieces)	348
<input type="checkbox"/> DUOLV2 Visco Overlay Full (6 pieces)	696
<input type="checkbox"/> DUOLG1 Gel Overlay Half (3 Pieces)	348
<input type="checkbox"/> DUOLG2 Gel Overlay Full (6 Pieces)	696

PIN OPTIONS	PRICING
<input type="checkbox"/> DUP1 1" Long Pins	40
(Standard) <input type="checkbox"/> DUP34 3/4" Long Pins	N/C

Key Features

- Provides aggressive laterals for complex posture support
- Comfortable High Resilience Foam
- 3lbs Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Adjustable/Removable Lumbar Support
- Easy to Install and Remove
- Height Adjustable

BACK SHELL SPECIFICATIONS



FULL LATERAL
 X = Standard shell height 18"
 Y = Lateral depth 5" finished
 Z = Lateral height 10"

DIAGRAM