

## MDU RELAX Snug Fit Back - Seating



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



STANDARD SIZES	PRICING	COVER OPTIONS	PRICING
<input type="checkbox"/> MDU16X18 Snug Fit Back W 16" 18" H	700	<input type="checkbox"/> MDUCOV Standard	N/C
<input type="checkbox"/> MDU18X18 Snug Fit Back W 18" 18" H	700	<input type="checkbox"/> MDUCT4 TEK4 Infection Control Cover	N/C
<input type="checkbox"/> MDU20X18 Snug Fit Back W 20" 18" H	700		

CUSTOM SIZES	Choose one below, then enter size	PRICING	EXTRA COVER	PRICING
<input type="checkbox"/> MDUC	Snug Fit Back Custom	700	<input type="checkbox"/> MDUCC Standard Cover	99
MDUCW	Width: _____	+176	<input type="checkbox"/> MDUOT4 TEK4 Infection Control Cover	99
MDUCH	Height: _____	+ 176	<input type="checkbox"/> MDUCD Custom Cover <input type="checkbox"/> Standard <input type="checkbox"/> TEK4	134
			<input type="checkbox"/> MDUCDW Width: _____	
			<input type="checkbox"/> MDUCDH Height: _____	

\*Note: Please enter the size above for standard or custom cover

ADDITIONAL OVERLAYS	PRICING
<input type="checkbox"/> MDUOLF Foam Pad	116
<input type="checkbox"/> MDUOLV1 Visco Overlay Half (3 Pieces)	348
<input type="checkbox"/> MDUOLV2 Visco Overlay Full (6 pieces)	696
<input type="checkbox"/> MDUOLG1 Gel Overlay Half (3 Pieces)	348
<input type="checkbox"/> MDUOLG2 Gel Overlay Full (6 Pieces)	696


PIN OPTIONS	PRICING
<input type="checkbox"/> MDUP1 1" Long Pins	40
(Standard) <input type="checkbox"/> MDUP34 3/4" Long Pins	N/C

### Key Features

- Provides aggressive laterals for complex posture support
- Comfortable High Resilience Foam
- 3lbs Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Adjustable/Removable Lumbar Support
- Easy to Install and Remove
- Height Adjustable

### BACK SHELL SPECIFICATIONS

### DIAGRAM



#### SNUG LATERAL

- W = Standard shell height 18"
- X = Lateral depth 5"
- Y = Space below lateral 4" finished
- Z = Lateral height 6"