

RELAX AIR FLOW Cushion - Seating

By  **SYSTEM**



INFORMATION

Customer Name:

Address:

City:

Account #:

Phone #:

Authorizer Name:

PO #:

Date:



P302C - Polyair Cushion Dual Valve

☐	ACD1616	AIR FLOW Cushion - Dual Valve W 16" X L 16"	\$704
☐	ACD1618	AIR FLOW Cushion - Dual Valve W 16" X L 18"	\$704
☐	ACD1816	AIR FLOW Cushion - Dual Valve W 18" X L 16"	\$704
☐	ACD1818	AIR FLOW Cushion - Dual Valve W 18" X L 18"	\$704
☐	ACD2018	AIR FLOW Cushion - Dual Valve W 20" X L 18"	\$704
☐	ACD2020	AIR FLOW Cushion - Dual Valve W 20" X L 20"	\$704
☐	ACD2218	AIR FLOW Cushion - Dual Valve W 22" X L 18"	\$704
☐	ACD2420	AIR FLOW Cushion - Dual Valve W 24" X L 20"	\$704

KEY FEATURES

- Four way stretchable POLYMAILLE® cover allows the air cells to conform to the clients body
- Removable POLYMAILLE® Cover
- Durable cover with a soft fabric
- Laxprene® air pockets for maximum comfort and softness
- Bottom of the cushion cover is based on a non-skid material to ensure secure placement when user is seated

EXTRA PRESSURE PUMP

☐	ACD624A	AIRFLOW Cushion Extra Pump	\$120
--------------------------	---------	-----------------------------------	-------


EXTRA COVER

☐	ACDCOV3	AIRFLOW Cushion Extra Removable POLYMAILLE® Cover	\$120
	ACDDW	Width:	
	ACDDH	Length:	


*Note: Please enter the size above for standard or custom cover

PRESSURE GRAPH

Pelvic Circumference Vs Weight



Example on how to use chart



1. Measure user's waist (Pelvic Circumference).
2. Measure user's weight.
3. Reference the chart and create a intersecting point between the weight and the pelvic circumference.
4. User's setup pressure level is where both points connect.
e.g: Pressure point is 5 for Pelvic Circumference 32 Inches and Weight 198 lbs.

NOTE: The chart is only a reference guide for the MINIMUM cushion pressure before you setup the user. Please ensure the user has 1 to 1.5 inch gap from the bottox to the seat pan for optimal comfort.