

VGP VOYAGER PLUS Wheelchair (Manual Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR		Select one below	1
<input type="checkbox"/>	VOYAGER PLUS FRAME - Black		1450

SEAT WIDTH		Select one below	2
<input type="checkbox"/>	SW14"-20"	Specify _____	N/C

SW14"-20" Light Weight Aluminum Frame

<input type="checkbox"/>	SW21"-28"	Specify _____	N/C
--------------------------	-----------	---------------	-----

SW21"-28" Must Include Plus Package 1, 2 Or 3

SEAT DEPTH		Select one below	3
------------	--	------------------	---

<input type="checkbox"/>	SDF18"	(By Frame)	N/C
<input type="checkbox"/>	SD16"-18"	Specify _____ (By Upholstery)	N/C
<input type="checkbox"/>	SDF20"	(By Frame)	350
<input type="checkbox"/>	SD18"-20"	Specify _____ (By Upholstery)	N/C
<input type="checkbox"/>	SDF22"	(By Frame)	350
<input type="checkbox"/>	SD20"-22"	Specify _____ (By Upholstery)	N/C

BACK		Select one below	4
------	--	------------------	---

<input type="checkbox"/>	BH15"-20"	Specify _____	N/C
<input type="checkbox"/>	BH21"+	Specify _____	120
<input type="checkbox"/>	AB	Angle Adj. Back	N/C
<input type="checkbox"/>	ENDO*	Endomorph Canes W ()	N/C

Specify Width

Endomorph Canes require TA or TAS Arm Rest

STROLLER HANDLE		Select one below	5
-----------------	--	------------------	---

<input type="checkbox"/>	SHS	Removable Stroller Handle	N/C
<input type="checkbox"/>	SHR	Removable Raised Stroller Handle	N/C

FOOT REST		Select one below	6
-----------	--	------------------	---

<input type="checkbox"/>	F6	60 Degree Hanger	N/C
<input type="checkbox"/>	F7	70 Degree Hanger	N/C
<input type="checkbox"/>	F8	80 Degree Hanger	N/C
<input type="checkbox"/>	EA	Elevating Legrest (A2)	404
<input type="checkbox"/>	AL	Articulating Legrest (A2)	450

FOOT REST AND FOOT PLATES		Select one below	7
---------------------------	--	------------------	---

<input type="checkbox"/>	PPL	Plastic	N/C
<input type="checkbox"/>	PAF	Angle Adjustable - Extra Large	150
<input type="checkbox"/>	PXL	Extra Large – Alum (F6A,F7A,F8A)	150
<input type="checkbox"/>	1FB	One Piece Foot Board	80

ARM REST		Select one below	8
----------	--	------------------	---

<input type="checkbox"/>	FA	Flip Back Ht. Adj./ Removable	240
<input type="checkbox"/>	TA*	T-Arm Ht. Adj./ Removable	240
<input type="checkbox"/>	FAS	Flip Back (Short) Ht. Adj./ Remove	240
<input type="checkbox"/>	TAS*	T-Arm (Short) Ht. Adj./ Removable	240
<input type="checkbox"/>	GWAPS	Gel Waterfall Arm Pads (A1)	150

14" - 20" Wide Light Weight Aluminum Frame Only

Please Note : **A1** = Clients Responsibility

Please Note : **A2** = Clinical Justification

SPECIFICATIONS

- LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- BLACK PADDED ARM PADS
- DOUBLE CROSS BRACE

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	SL 12.25	20"	5"	XL SHORT
<input type="checkbox"/>	SL 13.00"	20"	5"	XL SHORT
<input type="checkbox"/>	LL 13.75"	20"	5"	XL SHORT
<input type="checkbox"/>	L 14.50"	20"	6"	SHORT
<input type="checkbox"/>	L 14.50"	22"	5"	SHORT
<input type="checkbox"/>	L14.50"	22"	6"	SHORT
<input type="checkbox"/>	LH 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/>	SH 16.00"	22"	6"	LONG

(STANDARD SEAT SETTING)

<input type="checkbox"/>	MH 16.75	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	8"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	24"	8"	LONG
<input type="checkbox"/>	TH18.25"	24"	8"	LONG
<input type="checkbox"/>	TT 19.00"*	24"	8"	LONG

*TT 19.00' requires taller caster journal

REAR WHEELS <i>Select one below</i>			9
<input type="checkbox"/>	RWT0	20" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWT2	22" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWT4	24" MAG Solid Treaded Tire	N/C
<input type="checkbox"/>	RWP2	22" Pyramid Tire	80
<input type="checkbox"/>	RWP4	24" Pyramid Tire	80
<input type="checkbox"/>	RWA2	22" Air Tire	120
<input type="checkbox"/>	RWA4	24" Air Tire	120
<input type="checkbox"/>	AQR	Axles – Quick Release (349 Lbs)	N/C
<input type="checkbox"/>	PCHR	Plastic Coated Hand Rims	N/C

14" - 20" Wide Light Weight Aluminium Frame Only

CASTER <i>Select one below</i>			10
<input type="checkbox"/>	FC1A	5" X 1.25" Caster	N/C
<input type="checkbox"/>	FC1B	6" X 1.25" Caster	N/C
<input type="checkbox"/>	FC1C	8" X 1.25" Caster	N/C
<input type="checkbox"/>	FC2A	5" X 1.75" Caster	80
<input type="checkbox"/>	FC2B	6" X 1.75" Caster	80
<input type="checkbox"/>	FC2C	8" X 1.75" Caster	80
<input type="checkbox"/>	FC3A	5" X 1.25" Air Caster	80
<input type="checkbox"/>	FC3B	6" X 1.25" Air Caster	80
<input type="checkbox"/>	FC3C	8" X 1.25" Air Caster	80

WHEEL LOCK <i>Select one below</i>			11
<input type="checkbox"/>	ALPL	Aluminium Pull To Lock	N/C
<input type="checkbox"/>	ALPH	Aluminium Push To Lock	N/C
<input type="checkbox"/>	STPL	Steel Pull To Lock	N/C
<input type="checkbox"/>	STPH	Steel Push To Lock	N/C
<input type="checkbox"/>	BRXA	6" Brake Extension	40
<input type="checkbox"/>	BRXB	9" Brake Extension	60

SEATING - BACKS <i>Select one below</i>			12
W () x H ()			
<input type="checkbox"/>	TAB	Tension Adjustable Back	180

ACCESSORIES			14
<input type="checkbox"/>	SBA	Seat Belt (Auto Style)	22
<input type="checkbox"/>	SBB	Seat Belt (Velcro Style)	15.75
<input type="checkbox"/>	CNC	Crutch And Cane Holder (A1)	60
<input type="checkbox"/>	IVA	Adjustable I V Pole (A1)	130
<input type="checkbox"/>	NLSP	Nylon Leg Strap Padded (A1)	50
<input type="checkbox"/>	VGPAMPL	Amputee Support Left	45
<input type="checkbox"/>	VGPAMPR	Amputee Support Right	45
<input type="checkbox"/>	VGPOXY	Oxygen Tank Holder	130
<input type="checkbox"/>	VGPATR	Anti Tippers	88

Please Note : A1 = Clients Responsibility

PLUS PACKAGES <i>Select one below</i>		15
<input type="checkbox"/>	PLUS PACKAGE 1: AADL CODE: W101 CLIENT 250 TO 349 LBS - \$880.00 INCLUDES THE FOLLOWING: STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 16" TO 24" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/>	PLUS PACKAGE 2: AADL CODE: W102 CLIENT 350 TO 449 LBS - \$1140.00 INCLUDES THE FOLLOWING: STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 18" TO 26" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS HD STEEL FORKS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/>	PLUS PACKAGE 3: AADL CODE: W103 CLIENT 450 TO 600 LBS - \$1300.00 INCLUDES THE FOLLOWING: STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 20" TO 28" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS HD STEEL FORKS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/>	FIXED REINFORCED ONE PIECE BACK POST	\$260.00
	(SPECIFY ANGLE _____)	
<input type="checkbox"/>	HEAVY DUTY TENSION ADJ. BACK UPH.	\$220.00

Please Note : Reference Page 12 of the AADL Wheelchair APL for options and catalogue numbers for A wheelchairs.

Notes:

AADL CODE	PRODUCT	MSRP	DETAILS
BACKS			
N590	Relax Personal Fit - LESS THAN 3" LATERAL	□ \$690.00	
N590	Relax Gentle Fit - LESS THAN 3" LATERAL	□ \$690.00	
N593	Relax Snug Fit - GREATER THAN 3" LATERAL	□ \$790.00	
N593	Relax Full Fit - GREATER THAN 3" LATERAL	□ \$790.00	
N590	Relax Form Fit GENTLE - LESS THAN 3" LATERAL	□ \$690.00	
N593	Relax Form Fit SNUG - GREATER THAN 3" LATERAL	□ \$790.00	
N593	Relax Form Fit FULL - GREATER THAN 3" LATERAL	□ \$790.00	
N590	Relax Spiral Fit GENTLE - LESS THAN 3" LATERAL	□ \$690.00	
N593	Relax Spiral Fit SNUG - GREATER THAN 3" LATERAL	□ \$790.00	
N593	Relax Spiral Fit FULL - GREATER THAN 3" LATERAL	□ \$790.00	
N500	CUSTOM WIDTH CHARGE WHERE APPLICABLE - GREATER THAN 20" WIDE	□ \$300.00	
V853	Lumbar/Extra Pressure Reduction - Maximum 3 pieces per back	□ \$48.45	
V854	Lateral Pressure Reduction - Maximum 3 Pieces Per Back	□ \$48.45	
V855	Pelvic Pressure Reduction - Maximum 3 Pieces Per Back	□ \$48.45	
V856	Wedge Pressure Reduction - Maximum 3 Pieces Per Back	□ \$48.45	
V901	Thoracic Gel Pressure Reduction - Maximum 3 Pieces Per Back	□ \$69.00	
V902	Full Spinal Area Gel Pressure Reduction - Maximum 3 Pieces Per Back	□ \$103.50	
V903	Full Back Area Gel Pressure Reduction - Maximum 3 Pieces Per Back	□ \$137.08	
V571	REPLACEMENT BACK COVER FOR ALL SIZES	□ \$183.00	
CUSHIONS			
V980	Relax Comfy CUSHION - Full Time User Sizes Up-to 20"x20"	□ \$589.00	
V980 + V529	Relax Firm CUSHION - Full Time User Sizes Over 20"x20" (V529 = \$190)	□ \$779.00	
N570	Relax AIRFLOW CUSHION - Adjustable Positioning Up-to 20"x20"	□ \$719.00	
N570 + V529	Relax AIRFLOW CUSHION - Adjustable Positioning Over 20"x20" (V529 = \$190)	□ \$909.00	
V101	REPLACEMENT CUSHION COVER FOR ALL SIZES	□ \$165.00	
HEADRESTS			
V081	STANDARD PAD - F1 Contoured Pad Large	□ \$190.00	
N582	HEADREST WITH LATERAL SUPPORT - FXL1 Flex Pad Large	□ \$255.00	
V338	STANDARD STRAIGHT HARDWARE - LAH12 Vertical Horizontal Headrest Fixed 1/2"	□ \$165.00	
N533	MULTI-AXIS HARDWARE - MAHL Multi-Axis Headrest (L-Bend) 1/2"	□ \$200.00	
N535	GOOSENECK STYLE MULTI-AXIS - MAH Multi-Axis Headrest (S-Bend) 1/2"	□ \$280.00	
V085	MOUNTING PLATE - UHRA Universal Headrest Plate - 1/2" OR 3/4"	□ \$189.70	
V570	REPLACEMENT HEADREST COVER FOR ALL SIZES	□ \$120.00	
ACCESSORIES			
V577	ARM TROUGH TRAY ASSEMBLY - CHPTL-P = Left Tray, CHPTR-P = Right Tray	□ \$432.00	
V476	WHEELCHAIR TRAY - Standard or Padded, Full or Half Up-to 18" Wide	□ \$257.60	
V542	WHEELCHAIR TRAY - Standard or Padded, Full or Half Over 18" Wide	□ \$310.20	
V762	ELBOW STOPS	□ \$90.00	
V766	LATERAL SUPPORT	□ \$90.00	
V769	HAND BLOCK	□ \$90.00	
N814	LATERAL THIGH SUPPORT	□ \$150.00	
V619	MODIFIED FOOT SUPPORT - Footbox, Padded Footbox, Contracture Sling (V619 \$179.55 X 2)	□ \$359.10	
V581	AMPUTEE BOARD - Up-to 20" Wide	□ \$350.00	
V583	HEAVY DUTY AMPUTEE BOARD - Up-to 28" x 22"	□ \$412.50	
V579	PAD FOR AMPUTEE BOARD	□ \$62.50	

NOTES: