

NVT NV Dynamic Tilt Wheelchair



INFORMATION

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



DYNAMIC FRAME		1	ADP
<input checked="" type="checkbox"/>	NVTB Dynamic Tilt Frame - Black	2699	WM0001076
<input checked="" type="checkbox"/>	NVTWANA Standard Manual Wheelchair Frame with Dynamic	945	WANA
<input checked="" type="checkbox"/>	NVTWAMV Clothing Guard	99	WAMV

*Note: All three ADP codes apply when you order the NV Tilt Wheelchair

CLIENT WEIGHT		Select one below	2	ADP
<input type="checkbox"/>	NVTWA Up to 250Lbs		N/C	
<input type="checkbox"/>	NVTWAMF Exceeds 250Lbs		\$730	WAMF
<input type="checkbox"/>	NVTWAMG Exceeds 350Lbs		\$992	WAMG
<input type="checkbox"/>	NVTWAMH Exceeds 400Lbs (Max 600Lbs)		\$1,110	WAMH

SEAT WIDTH		Select one below	3	ADP
<input type="checkbox"/>	NVTSW18 Seat Width 18"		N/C	
<input type="checkbox"/>	NVTSW19 Seat Width 19"		\$399	WAMA
<input type="checkbox"/>	NVTSW20 Seat Width 20"		\$399	WAMA
<input type="checkbox"/>	NVTSWC Custom Seat Width: _____		\$399	WAMA

SEAT DEPTH		Select one below	4	ADP
<input type="checkbox"/>	NVTS16 Seat Depth 16"		N/C	
<input type="checkbox"/>	NVTS17 Seat Depth 17"		N/C	
<input type="checkbox"/>	NVTS18 Seat Depth 18"		N/C	
<input type="checkbox"/>	NVTS19 Seat Depth 19"		\$265	WAMB
<input type="checkbox"/>	NVTS20 Seat Depth 20"		\$265	WAMB
<input type="checkbox"/>	NVTSDC Custom Seat Depth: _____		\$265	WAMB

STANDARD NO CHARGE SPECIFICATIONS	
SEAT WIDTH: 14" – 18"	
SEAT DEPTH (ADJUSTABLE)	
SEAT HEIGHT (ADJUSTABLE): 16.75", 17.75" and 18.75"	
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]	
WEIGHT CAPACITY: 250 LBS	
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)	
TILT ADJUSTMENTS 0 TO 30 (DEGREES)	
TWO-WAY LOCKING GAS STRUT	
SINGLE RELEASE, DUAL CYLINDERS	
STANDARD LOCKOUT COLLARS	
STANDARD T STYLE ARMS	
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS	

OTHER FRAME COLOUR		1	ADP
<input type="checkbox"/>	NVTS Tilt Frame - Silver	N/C	WM0001040
<input type="checkbox"/>	NVTBL Tilt Frame - Blue	N/C	WM0001040
<input type="checkbox"/>	NVTBU Tilt Frame - Burgundy	N/C	WM0001040
<input type="checkbox"/>	NVTC Tilt Frame - Custom : _____	\$400	WM0001040

FOOT REST AND FOOT PLATES		Select one below	5	ADP
<input type="checkbox"/>	NVTOFST Offset Foot Rest		N/C	
<input type="checkbox"/>	NVTF6A 60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTF7A 70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTF8A 80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTEA Elevating Legrest <input type="checkbox"/> Pivot		\$404	WAMK
<input type="checkbox"/>	NVTAL Articulating L/R <input type="checkbox"/> Pivot		\$450	WAMK
<input type="checkbox"/>	NVTEB Elevating Legrest System		\$450	WAMK
<input type="checkbox"/>	NVTPPL Plastic Footplates		N/C	
<input type="checkbox"/>	NVTPXL Extra Large Footplates		\$150	SE0000322
<input type="checkbox"/>	NVTPAF Angle Adjustable Footplates		\$240	WAMI
<input type="checkbox"/>	NVT1FB One Piece Foot Board		\$80	SEFND0015
<input type="checkbox"/>	NVT1FBP One Piece Foot Board Padding		\$44	SEFND0035
<input type="checkbox"/>	NVTNLS Nylon Leg Strap Padded		\$60	SEFND0010
<input type="checkbox"/>	NVTNLS Custom Size ()		\$120	SEFCF107L
<input type="checkbox"/>	NVTCSL Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style		\$240	SEFND0020

TILT OPTION		Select one below	6	ADP
<input type="checkbox"/>	NVTINTS Independant Power Tilt System		\$2,000	

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/>	16.75"	20"	6-8"	LONG
<input type="checkbox"/>	17.75"	22"	6-8"	LONG
<input type="checkbox"/>	18.75"	24"	6-8"	LONG

BACK TYPE		Select one below	7	ADP
<input type="checkbox"/>	NVTAB1 Angle Adjustable Back		N/C	
<input type="checkbox"/>	NVTENDO Endomorph Canes W () <input type="checkbox"/> Standard		N/C	
<input type="checkbox"/>	NVTRBK Recline Kit		\$560	WAMI

HEADREST		Select one below	8	ADP
<input type="checkbox"/>	NVTMAHF Multi-Axis Hardware Complete with Large Pad		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFXL Multi-Axis Hardware Complete with Flex Large		\$512	SE0001542
<input type="checkbox"/>	NVTMAHKP Multi-Axis Hardware Complete with Medium Pad		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFXS Multi-Axis Hardware Complete with Flex Medium		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFT Multi-Axis Hardware Complete with Small Flat		\$512	SE0001542
<input type="checkbox"/>	NVTUPHT219 Contour Headrest Complete with Swivel Hardware		\$433	SEHND0005 SEHCF050L
<input type="checkbox"/>	NVTUPHT218 Contour Headrest with Std. Assembly (1/2" Tube)		\$292	SEHND0020
<input type="checkbox"/>	NVTUPHT217 Contour Headrest with Std. Assembly (3/4" Tube)		\$221	SEHWM0005
<input type="checkbox"/>	NVTUHRA Universal Hardware Adaptor (Adaptor only)		\$86	SEHND0020

ARM REST		Select one below	9	ADP	SEATING - BACKS		Custom: select back & fill custom	13	ADP
(Ht. 10"-14")	<input type="checkbox"/> NVTTA	T-Arm Ht. Adj./ Removable	\$240		<input type="checkbox"/> NVTBK01	RELAX Gentle Fit Back	W () x H ()	\$544	SEBWM0020
(Ht. 10"-14")	<input type="checkbox"/> NVTTAS	T-Arm Ht. Adj./ Removable Short	\$240		<input type="checkbox"/> NVTBK02	RELAX Full Fit Back		\$608	SEBWM0010
	<input type="checkbox"/> NVTGWAPF	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030	<input type="checkbox"/> NVTBK03	RELAX Snug Fit Back		\$700	SEBWM0027
*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required									
REAR WHEEL		Select one below	10	ADP	<input type="checkbox"/> NVTBK04	RELAX Form Fit Back No Laterals		\$820	SEBWM0026
<input type="checkbox"/> NVTRWT0		20" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK06	RELAX Form Fit Back Gentle Laterals		\$820	SEBWM0026
<input type="checkbox"/> NVTRWT2		22" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK07	RELAX Form Fit Back Full Laterals		\$820	SEBWM0026
<input type="checkbox"/> NVTRWT4		24" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK05	RELAX Form Fit Back Snug Laterals		\$820	SEBWM0026
<input type="checkbox"/> NVTRWA2		22" Air Tire	\$120		<input type="checkbox"/> NVTBK10	RELAX Spiral Fit Back Gentle Laterals		\$832	SEBWM0030
<input type="checkbox"/> NVTRWA4		24" Air Tire	\$120		<input type="checkbox"/> NVTBK09	RELAX Spiral Fit Back Snug Laterals		\$832	SEBWM0030
<input type="checkbox"/> NVTAQR		Axles – Quick Release	\$120	WAMU	<input type="checkbox"/> NVTBK11	RELAX Spiral Fit Back Full Laterals		\$832	SEBWM0030
(Standard)	<input type="checkbox"/> NVTHRT	Narrow Handrim Config. (1/2" less)	N/C		<input type="checkbox"/> NVTTAB	RELAX Tension Back Upto 18"		\$140	WAMD
*Note: Quick Release is not available for wheels smaller than 20"									
CASTER		Select one below	11	ADP	<input type="checkbox"/> NVTTAB+	RELAX Tension Back Over 18"		\$240	WAME
<input type="checkbox"/> NVTFC1B		6" X 1.25" Caster	N/C		<input type="checkbox"/> NVTBK12	* Custom Width : _____		\$176	SEBWM0035
<input type="checkbox"/> NVTFC1C		8" X 1.25" Caster	N/C		<input type="checkbox"/> NVTBK13	* Custom Height : _____		\$176	SEBWM0040
<input type="checkbox"/> NVTFC2B		6" X 1.75" Caster	\$80		*Note: Please refer to the Seating Booklet for custom sizes				
<input type="checkbox"/> NVTFC2C		8" X 1.75" Caster	\$80		<input type="checkbox"/> NVTBGL1	Gel Full Overlay	Maximum 1 piece additional pressure relief gel	\$348	SEBND0020 X 3
<input type="checkbox"/> NVTFC4B		6" X 2 Caster	\$80		<input type="checkbox"/> NVTBGL2	Gel Center Vertebrae		\$116	SEBND0020
<input type="checkbox"/> NVTFC4C		8" X 2 Caster	\$80		<input type="checkbox"/> NVTBGL3	Gel Left Lateral		\$116	SEBND0020
					<input type="checkbox"/> NVTBGL4	Gel Right Lateral		\$116	SEBND0020
WHEEL LOCK		Select one below	12	ADP	SEATING - CUSHIONS		Select one below	14	ADP
<input type="checkbox"/> NVTALPH		Aluminium Push To Lock	N/C			W () x D ()			
<input type="checkbox"/> NVTALPL		Aluminium Pull To Lock	N/C		<input type="checkbox"/> NVTC01	RELAX Soft Cushion		\$480	SESWM0050
<input type="checkbox"/> NVTBRXC		4" Brake Extension	N/C		<input type="checkbox"/> NVTC02	RELAX Comfy Cushion		\$624	SESWM0030
<input type="checkbox"/> NVTBRXA		6" Brake Extension	N/C		<input type="checkbox"/> NVTC03	RELAX Firm Cushion		\$784	SESWM0040
<input type="checkbox"/> NVTBRXB		9" Brake Extension	\$130		<input type="checkbox"/> NVTC04	RELAX AIR FLOW Cushion		\$587	SE0001669
<input type="checkbox"/> NVTATTBRX		Attendant Brakes	\$240		<input type="checkbox"/> NVTCGL2	Full Overlay Gel		\$258	SESND1045
<input type="checkbox"/> NVTULB		Unilateral Hand Brake <input type="checkbox"/> Left <input type="checkbox"/> Right	\$120		<input type="checkbox"/> NVTCCL	Contoured Leg Channel (Bolster) Std.		\$180	SEMCFF203L
					<input type="checkbox"/> NVTCCLA	Contoured Leg Channel (Bolster) with ABS		\$260	SESND1035 SESND1015
					<input type="checkbox"/> NVTCCLC	* Contoured Leg Channel (Bolster) Custom		\$320	SEMCFF209L
*Note: Please refer to the Seating Booklet for custom sizes									
					<input type="checkbox"/> NVTC00	Extra Cushion Cover		\$99	SESND1055
ACCESSORIES			15	ADP	<input type="checkbox"/> NVTSBA	Seat Belt (Auto Style)		\$26	
<input type="checkbox"/> NVTSBA		Seat Belt (Auto Style)	\$26		<input type="checkbox"/> NVTSBAP	Seat Belt (Padded)		\$56	SERND0030
<input type="checkbox"/> NVTSBA4P		Seat Belt (4 Point Padded)	\$130		<input type="checkbox"/> NVTSBA4P	Seat Belt (4 Point Padded)		\$130	SERND0020
<input type="checkbox"/> NVTTRA		Full Padded Tray Locking	\$280	SETND0005 SETND0035 SETND0005	<input type="checkbox"/> NVTTRB	Full Padded Tray Slide On		\$280	SETND0035
<input type="checkbox"/> NVTTRB		Full Padded Tray Slide On	\$280	SETND0035	<input type="checkbox"/> NVTTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$280	SETND0010
<input type="checkbox"/> NVTTRCL		Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010	<input type="checkbox"/> NVTTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$280	SETND0010
<input type="checkbox"/> NVTTRC		Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010	<input type="checkbox"/> NVTTRD	Half Padded Left Tray		\$240	SETND0015 SETND0035 SETND0015
<input type="checkbox"/> NVTTRD		Half Padded Left Tray	\$240	SETND0015 SETND0035 SETND0015	<input type="checkbox"/> NVTTRE	Half Padded Right Tray		\$240	SETND0035
<input type="checkbox"/> NVTTRE		Half Padded Right Tray	\$240	SETND0035	<input type="checkbox"/> NVTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$240	SETND0015
<input type="checkbox"/> NVTTRF		Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015	<input type="checkbox"/> NVTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$240	SETND0015
<input type="checkbox"/> NVTTRG		Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015	<input type="checkbox"/> NVTTDFB	Transport Tie Downs Front		\$85	
<input type="checkbox"/> NVTTDFB		Transport Tie Downs Front	\$85		<input type="checkbox"/> NVTTDRB	Transport Tie Downs Rear		\$85	
<input type="checkbox"/> NVTTDRB		Transport Tie Downs Rear	\$85		<input type="checkbox"/> NVTOXY	Oxygen Tank Holder		\$195	WEPX
<input type="checkbox"/> NVTOXY		Oxygen Tank Holder	\$195	WEPX	<input type="checkbox"/> NVTIVA	Adjustable I V Pole		\$130	
<input type="checkbox"/> NVTIVA		Adjustable I V Pole	\$130		<input type="checkbox"/> NVTCHPT-L	Left Arm Trough Tray		\$400	SETND0015 SETND0035 SETND0030
<input type="checkbox"/> NVTCHPT-L		Left Arm Trough Tray	\$400	SETND0015 SETND0035 SETND0030	<input type="checkbox"/> NVTCHPT-R	Right Arm Trough Tray		\$400	SETND0030
<input type="checkbox"/> NVTCHPT-R		Right Arm Trough Tray	\$400	SETND0030	<input type="checkbox"/> NVTLATB	Lateral Bite out		\$180	SELND0015
<input type="checkbox"/> NVTLATB		Lateral Bite out	\$180	SELND0015	<input type="checkbox"/> NVTADD8	8" Adductor		\$164	SELND0005
<input type="checkbox"/> NVTADD8		8" Adductor	\$164	SELND0005	<input type="checkbox"/> NVTADD10	10" Adductor		\$164	SELND0005
<input type="checkbox"/> NVTADD10		10" Adductor	\$164	SELND0005	<input type="checkbox"/> NVTADD12	12" Adductor		\$164	SELND0005
<input type="checkbox"/> NVTADD12		12" Adductor	\$164	SELND0005	<input type="checkbox"/> NVTFRCH	Padded Footplates		\$88	SEMCFF803L
<input type="checkbox"/> NVTFRCH		Padded Footplates	\$88	SEMCFF803L					