

NV TILT

NARROW VERSION TILT

DYNAMIC TILT WHEELCHAIR



ADP CODE: WM0001076

NV Tilt reimagines mobility with its innovative design as a narrow-width ***Dynamic Tilt*** wheelchair, expertly crafted for easy maneuverability and increased mobility through doorways and tight spaces. Engineered for both practicality and comfort, this wheelchair's slender frame with a wide seating surface that doesn't compromise on stability or user experience.




Adjustable Stroller Handle



Upper Frame Back Cane, and Front Rigging Adjustable to Balance Center of Gravity



Narrow Wheel Width



Seat Pan Rigid Skirt Guards and Endomorph Canes



Rear Axle Adjustability Maximizing Hand to Wheel Contact



Removable T-Style Arms

NV Tilt reimagines mobility with its innovative design as a narrow-width **Dynamic Tilt** wheelchair, expertly crafted for easy maneuverability and increased mobility through doorways and tight spaces. Engineered for both practicality and comfort, this wheelchair's slender frame with a wide seating surface that doesn't compromise on stability or user experience. The NV Tilt prioritizes accessibility, making it an ideal choice for individuals seeking enhanced freedom and independence in their daily lives. Its sleek and ergonomic construction embodies an era of wheelchair mobility, where seamless movement through various environments is a reality.

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to users' specifications from either a drawing, pictures or a sample. The custom parts we make can be simple or complex and can be produced in small or large quantities based on what is requested. Contact your local sales representative with your specifications.

NV TILT

NARROW VERSION TILT

DYNAMIC TILT WHEELCHAIR



STANDARD FEATURES

- Narrow Overall Width with Wide Seating Surface
- Depth Adjustable Frame for Center of Gravity Adjustment
- ABS Seat Pan with Rigid Skirt Guard
- Dual Gas Struts For Easy Tilt (Up To 30°)
- Tilt Collar lock, Locked Out at Your Position
- Recessed Lower Frame For increased maneuverability
- Dual Gas Cylinders for Smooth & Easy Tilting
- Soft Set Gas Cylinder Release Mechanism
- Black Frame Colour - Standard
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals

SPECIFICATIONS	STP
SEAT WIDTH:	18" TO 36"
SEAT DEPTH:	16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 16.75", 17.75" and 18.75"
TILT ANGLE:	0° TO 30°
BACK ANGLE:	(ADJUSTABLE): 0° , 5° , 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS

WELCOME HOME
TILT
TILT WHEELCHAIR



Learn more about this product by scanning the QR Code



NVT NV Dynamic Tilt Wheelchair



INFORMATION

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



DYNAMIC FRAME		1	ADP
<input checked="" type="checkbox"/>	NVTB Dynamic Tilt Frame - Black	2699	WM0001076
<input checked="" type="checkbox"/>	NVTWANA Standard Manual Wheelchair Frame with Dynamic	945	WANA
<input checked="" type="checkbox"/>	NVTWAMV Clothing Guard	99	WAMV

*Note: All three ADP codes apply when you order the NV Tilt Wheelchair

CLIENT WEIGHT		Select one below	2	ADP
<input type="checkbox"/>	NVTWA Up to 250Lbs		N/C	
<input type="checkbox"/>	NVTWAMF Exceeds 250Lbs		\$730	WAMF
<input type="checkbox"/>	NVTWAMG Exceeds 350Lbs		\$992	WAMG
<input type="checkbox"/>	NVTWAMH Exceeds 400Lbs (Max 600Lbs)		\$1,110	WAMH

SEAT WIDTH		Select one below	3	ADP
<input type="checkbox"/>	NVTSW18 Seat Width 18"		N/C	
<input type="checkbox"/>	NVTSW19 Seat Width 19"		\$399	WAMA
<input type="checkbox"/>	NVTSW20 Seat Width 20"		\$399	WAMA
<input type="checkbox"/>	NVTSWC Custom Seat Width: _____		\$399	WAMA

SEAT DEPTH		Select one below	4	ADP
<input type="checkbox"/>	NVTS16 Seat Depth 16"		N/C	
<input type="checkbox"/>	NVTS17 Seat Depth 17"		N/C	
<input type="checkbox"/>	NVTS18 Seat Depth 18"		N/C	
<input type="checkbox"/>	NVTS19 Seat Depth 19"		\$265	WAMB
<input type="checkbox"/>	NVTS20 Seat Depth 20"		\$265	WAMB
<input type="checkbox"/>	NVTSDC Custom Seat Depth: _____		\$265	WAMB

OTHER FRAME COLOUR		1	ADP
<input type="checkbox"/>	NVTS Tilt Frame - Silver	N/C	WM0001040
<input type="checkbox"/>	NVTBL Tilt Frame - Blue	N/C	WM0001040
<input type="checkbox"/>	NVTBU Tilt Frame - Burgundy	N/C	WM0001040
<input type="checkbox"/>	NVTC Tilt Frame - Custom : _____	\$400	WM0001040

FOOT REST AND FOOT PLATES		Select one below	5	ADP
<input type="checkbox"/>	NVTOFST Offset Foot Rest		N/C	
<input type="checkbox"/>	NVTF6A 60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTF7A 70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTF8A 80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/>	NVTEA Elevating Legrest <input type="checkbox"/> Pivot		\$404	WAMK
<input type="checkbox"/>	NVTAL Articulating L/R <input type="checkbox"/> Pivot		\$450	WAMK
<input type="checkbox"/>	NVTEB Elevating Legrest System		\$450	WAMK
<input type="checkbox"/>	NVTPPL Plastic Footplates		N/C	
<input type="checkbox"/>	NVTPXL Extra Large Footplates		\$150	SE0000322
<input type="checkbox"/>	NVTPAF Angle Adjustable Footplates		\$240	WAMI
<input type="checkbox"/>	NVT1FB One Piece Foot Board		\$80	SEFND0015
<input type="checkbox"/>	NVT1FBP One Piece Foot Board Padding		\$44	SEFND0035
<input type="checkbox"/>	NVTNLS Nylon Leg Strap Padded		\$60	SEFND0010
<input type="checkbox"/>	NVTNLS Custom Size ()		\$120	SEFCF107L
<input type="checkbox"/>	NVTCSL Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style		\$240	SEFND0020

TILT OPTION		Select one below	6	ADP
<input type="checkbox"/>	NVTINTS Independant Power Tilt System		\$2,000	

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/>	16.75"	20"	6-8"	LONG
<input type="checkbox"/>	17.75"	22"	6-8"	LONG
<input type="checkbox"/>	18.75"	24"	6-8"	LONG

BACK TYPE		Select one below	7	ADP
<input type="checkbox"/>	NVTAB1 Angle Adjustable Back		N/C	
<input type="checkbox"/>	NVTENDO Endomorph Canes W () <input type="checkbox"/> Standard		N/C	
<input type="checkbox"/>	NVTRBK Recline Kit		\$560	WAMI

HEADREST		Select one below	8	ADP
<input type="checkbox"/>	NVTMAHF Multi-Axis Hardware Complete with Large Pad		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFXL Multi-Axis Hardware Complete with Flex Large		\$512	SE0001542
<input type="checkbox"/>	NVTMAHKP Multi-Axis Hardware Complete with Medium Pad		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFXS Multi-Axis Hardware Complete with Flex Medium		\$512	SE0001542
<input type="checkbox"/>	NVTMAHFT Multi-Axis Hardware Complete with Small Flat		\$512	SE0001542
<input type="checkbox"/>	NVTUPHT219 Contour Headrest Complete with Swivel Hardware		\$433	SEHND0005 SEHCF050L
<input type="checkbox"/>	NVTUPHT218 Contour Headrest with Std. Assembly (1/2" Tube)		\$292	SEHND0020
<input type="checkbox"/>	NVTUPHT217 Contour Headrest with Std. Assembly (3/4" Tube)		\$221	SEHWM0005
<input type="checkbox"/>	NVTUHRA Universal Hardware Adaptor (Adaptor only)		\$86	SEHND0020

ARM REST		Select one below	9	ADP	SEATING - BACKS		Custom: select back & fill custom	13	ADP
(Ht. 10"-14")	<input type="checkbox"/> NVTTA	T-Arm Ht. Adj./ Removable	\$240		<input type="checkbox"/> NVTBK01	RELAX Gentle Fit Back	W () x H ()	\$544	SEBWM0020
(Ht. 10"-14")	<input type="checkbox"/> NVTTAS	T-Arm Ht. Adj./ Removable Short	\$240		<input type="checkbox"/> NVTBK02	RELAX Full Fit Back		\$608	SEBWM0010
	<input type="checkbox"/> NVTGWAPF	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030	<input type="checkbox"/> NVTBK03	RELAX Snug Fit Back		\$700	SEBWM0027
*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required									
REAR WHEEL		Select one below	10	ADP	<input type="checkbox"/> NVTBK04	RELAX Form Fit Back No Laterals		\$820	SEBWM0026
<input type="checkbox"/>	NVTRWTO	20" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK06	RELAX Form Fit Back Gentle Laterals		\$820	SEBWM0026
<input type="checkbox"/>	NVTRWT2	22" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK07	RELAX Form Fit Back Full Laterals		\$820	SEBWM0026
<input type="checkbox"/>	NVTRWT4	24" Solid Treaded Tire	N/C		<input type="checkbox"/> NVTBK05	RELAX Form Fit Back Snug Laterals		\$820	SEBWM0026
<input type="checkbox"/>	NVTRWA2	22" Air Tire	\$120		<input type="checkbox"/> NVTBK10	RELAX Spiral Fit Back Gentle Laterals		\$832	SEBWM0030
<input type="checkbox"/>	NVTRWA4	24" Air Tire	\$120		<input type="checkbox"/> NVTBK09	RELAX Spiral Fit Back Snug Laterals		\$832	SEBWM0030
<input type="checkbox"/>	NVTAQR	Axles – Quick Release	\$120	WAMU	<input type="checkbox"/> NVTBK11	RELAX Spiral Fit Back Full Laterals		\$832	SEBWM0030
(Standard)	<input type="checkbox"/> NVTHRT	Narrow Handrim Config. (1/2" less)	N/C		<input type="checkbox"/> NVTTAB	RELAX Tension Back Upto 18"		\$140	WAMD
*Note: Quick Release is not available for wheels smaller than 20"									
CASTER		Select one below	11	ADP	<input type="checkbox"/> NVTTAB+	RELAX Tension Back Over 18"		\$240	WAME
<input type="checkbox"/>	NVTFC1B	6" X 1.25" Caster	N/C		<input type="checkbox"/> NVTBK12	* Custom Width : _____		\$176	SEBWM0035
<input type="checkbox"/>	NVTFC1C	8" X 1.25" Caster	N/C		<input type="checkbox"/> NVTBK13	* Custom Height : _____		\$176	SEBWM0040
<input type="checkbox"/>	NVTFC2B	6" X 1.75" Caster	\$80		*Note: Please refer to the Seating Booklet for custom sizes				
<input type="checkbox"/>	NVTFC2C	8" X 1.75" Caster	\$80		<input type="checkbox"/> NVTBGL1	Gel Full Overlay	Maximum 1 piece additional pressure relief gel	\$348	SEBND0020 X 3
<input type="checkbox"/>	NVTFC4B	6" X 2 Caster	\$80		<input type="checkbox"/> NVTBGL2	Gel Center Vertebrae		\$116	SEBND0020
<input type="checkbox"/>	NVTFC4C	8" X 2 Caster	\$80		<input type="checkbox"/> NVTBGL3	Gel Left Lateral		\$116	SEBND0020
<input type="checkbox"/>					<input type="checkbox"/> NVTBGL4	Gel Right Lateral		\$116	SEBND0020
WHEEL LOCK		Select one below	12	ADP	SEATING - CUSHIONS		Select one below	14	ADP
<input type="checkbox"/>	NVTALPH	Aluminium Push To Lock	N/C			W () x D ()			
<input type="checkbox"/>	NVTALPL	Aluminium Pull To Lock	N/C		<input type="checkbox"/> NVTC01	RELAX Soft Cushion		\$480	SESWM0050
<input type="checkbox"/>	NVTBRXC	4" Brake Extension	N/C		<input type="checkbox"/> NVTC02	RELAX Comfy Cushion		\$624	SESWM0030
<input type="checkbox"/>	NVTBRXA	6" Brake Extension	N/C		<input type="checkbox"/> NVTC03	RELAX Firm Cushion		\$784	SESWM0040
<input type="checkbox"/>	NVTBRXB	9" Brake Extension	\$130		<input type="checkbox"/> NVTC04	RELAX AIR FLOW Cushion		\$587	SE0001669
<input type="checkbox"/>	NVTATTBRX	Attendant Brakes	\$240		<input type="checkbox"/> NVTCGL2	Full Overlay Gel		\$258	SESND1045
<input type="checkbox"/>	NVTULB	Unilateral Hand Brake <input type="checkbox"/> Left <input type="checkbox"/> Right	\$120		<input type="checkbox"/> NVTCCL	Contoured Leg Channel (Bolster) Std.		\$180	SEMCFF203L
*Note: Please refer to the Seating Booklet for custom sizes									
<input type="checkbox"/>					<input type="checkbox"/> NVTCCLA	Contoured Leg Channel (Bolster) with ABS		\$260	SESND1035 SESND1015
<input type="checkbox"/>					<input type="checkbox"/> NVTCCLC	* Contoured Leg Channel (Bolster) Custom		\$320	SEMCFF209L
<input type="checkbox"/>					<input type="checkbox"/> NVTC00	Extra Cushion Cover		\$99	SESND1055
ACCESSORIES			15	ADP	<input type="checkbox"/> NVTSBA	Seat Belt (Auto Style)		\$26	
<input type="checkbox"/>	NVTSBA	Seat Belt (Auto Style)	\$26		<input type="checkbox"/> NVTSBAP	Seat Belt (Padded)		\$56	SERND0030
<input type="checkbox"/>	NVTSBA4P	Seat Belt (4 Point Padded)	\$130		<input type="checkbox"/> NVTTRA	Full Padded Tray Locking		\$280	SERND0020
<input type="checkbox"/>	NVTTRB	Full Padded Tray Slide On	\$280	SETND0005 SETND0035 SETND0005 SETND0035	<input type="checkbox"/> NVTTRC	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$280	SETND0010
<input type="checkbox"/>	NVTTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280		<input type="checkbox"/> NVTTRD	Half Padded Left Tray		\$240	SETND0015 SETND0035 SETND0015 SETND0035
<input type="checkbox"/>	NVTTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280		<input type="checkbox"/> NVTTRE	Half Padded Right Tray		\$240	SETND0015
<input type="checkbox"/>	NVTTRD	Half Padded Left Tray	\$240		<input type="checkbox"/> NVTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$240	SETND0015
<input type="checkbox"/>	NVTTRE	Half Padded Right Tray	\$240		<input type="checkbox"/> NVTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear		\$240	SETND0015
<input type="checkbox"/>	NVTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240		<input type="checkbox"/> NVTTDFB	Transport Tie Downs Front		\$85	
<input type="checkbox"/>	NVTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240		<input type="checkbox"/> NVTTDRB	Transport Tie Downs Rear		\$85	
<input type="checkbox"/>	NVTTDFB	Transport Tie Downs Front	\$85		<input type="checkbox"/> NVTOXY	Oxygen Tank Holder		\$195	WEPX
<input type="checkbox"/>	NVTTDRB	Transport Tie Downs Rear	\$85		<input type="checkbox"/> NVTIVA	Adjustable I V Pole		\$130	
<input type="checkbox"/>	NVTOXY	Oxygen Tank Holder	\$195		<input type="checkbox"/> NVTCHPT-L	Left Arm Trough Tray		\$400	SETND0015 SETND0035 SETND0030
<input type="checkbox"/>	NVTIVA	Adjustable I V Pole	\$130		<input type="checkbox"/> NVTCHPT-R	Right Arm Trough Tray		\$400	SETND0015 SETND0035 SETND0030
<input type="checkbox"/>	NVTCHPT-L	Left Arm Trough Tray	\$400		<input type="checkbox"/> NVTLATB	Lateral Bite out		\$180	SELND0015
<input type="checkbox"/>	NVTCHPT-R	Right Arm Trough Tray	\$400		<input type="checkbox"/> NVTADD8	8" Adductor		\$164	SELND0005
<input type="checkbox"/>	NVTLATB	Lateral Bite out	\$180		<input type="checkbox"/> NVTADD10	10" Adductor		\$164	SELND0005
<input type="checkbox"/>	NVTADD8	8" Adductor	\$164		<input type="checkbox"/> NVTADD12	12" Adductor		\$164	SELND0005
<input type="checkbox"/>	NVTADD10	10" Adductor	\$164		<input type="checkbox"/> NVTFRCH	Padded Footplates		\$88	SEMCFF803L
<input type="checkbox"/>	NVTADD12	12" Adductor	\$164						
<input type="checkbox"/>	NVTFRCH	Padded Footplates	\$88						