

**STP STP Wheelchair (Tilt Wheelchair)**



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_

FRAME	Select one below	1	STANDARD NO CHARGE SPECIFICATIONS			
<input type="checkbox"/> STPB	Tilt Frame - Black	2499	SEAT WIDTH: 14" – 18"			
	Tilt Adjustment: <input type="checkbox"/> 0-35° <input type="checkbox"/> 0-45°	N/C	SEAT DEPTH (ADJUSTABLE)			
CLIENT WEIGHT	Select one below	2	SEAT HEIGHT (ADJUSTABLE): 10.25" – 18.25"			
<input type="checkbox"/> STPWA	Up to 250Lbs	N/C	BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]			
<input type="checkbox"/> STPWAMF	Exceeds 250Lbs	730	ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]			
<input type="checkbox"/> STPWAMG	Exceeds 350Lbs	900	WEIGHT CAPACITY: 250 LBS			
<input type="checkbox"/> STPWAMF	Exceeds 450Lbs	1110	WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)			
SEAT WIDTH	Select one below	3	TILT ADJUSTMENTS 0 – 35 [DEGREES]			
<input type="checkbox"/> STPSW14	Seat Width 14"	N/C	TWO-WAY LOCKING GAS STRUT			
<input type="checkbox"/> STPSW15	Seat Width 15"	N/C	SINGLE RELEASE, DUAL CYLINDERS			
<input type="checkbox"/> STPSW16	Seat Width 16"	N/C	STANDARD LOCKOUT COLLARS			
<input type="checkbox"/> STPSW17	Seat Width 17"	N/C	REAR HEIGHT ADJUSTABLE ANTI TIPPERS			
<input type="checkbox"/> STPSW18	Seat Width 18"	N/C				
<input type="checkbox"/> STPSW19	Seat Width 19"	200				
<input type="checkbox"/> STPSW20	Seat Width 20"	200				
<small>21-30" Width</small> <input type="checkbox"/> STPSW21-30	Specify _____	450				
SEAT DEPTH	Select one below	4	SEAT TO FLOOR SPECIFICATIONS			
<input type="checkbox"/> STPSD16	Seat Depth 16"	N/C	LOW PIVOT POINT SETTING			
<input type="checkbox"/> STPSD17	Seat Depth 17"	N/C				
<input type="checkbox"/> STPSD18	Seat Depth 18"	N/C	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> STPSD19	Seat Depth 19"	200	<input type="checkbox"/> 10.50"	16"	3"	2 HOLE
<input type="checkbox"/> STPSD20	Seat Depth 20"	200	<input type="checkbox"/> 13.00"	20"	5"	2 HOLE
<small>21-24" Depth</small> <input type="checkbox"/> STPSD21-24	Specify _____	450	<input type="checkbox"/> 13.75"	20"	5"	2 HOLE
			<input type="checkbox"/> 14.50"	20"	6"	3 HOLE
			<input type="checkbox"/> 14.50"	22"	5"	3 HOLE
			<input type="checkbox"/> 14.50"	22"	6"	3 HOLE
			<input type="checkbox"/> 15.25"	20"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
			<input type="checkbox"/> 15.25"	20"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
			<input type="checkbox"/> 15.25"	22"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
			<input type="checkbox"/> 15.25"	22"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
			<input type="checkbox"/> 16.00"	22"	6"	5 HOLE
FOOT REST AND FOOT PLATES	Select one below	5	HIGH PIVOT POINT SETTING			
<input type="checkbox"/> STPOFST	Offset Foot Rest	N/C	<input type="checkbox"/> 16.75"	22"	6"	5 HOLE
<input type="checkbox"/> STPF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	<input type="checkbox"/> 17.50"	22"	6"	5 HOLE
<input type="checkbox"/> STPF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	<input type="checkbox"/> 17.50"	24"	6"	5 HOLE
<input type="checkbox"/> STPF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	<input type="checkbox"/> 18.25"	24"	8"	5 HOLE
<input type="checkbox"/> STPEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin (A2)	404	For more Seat To Floor combinations, visit our website <a href="https://powerplumobility.com/stf">https://powerplumobility.com/stf</a> Please confirm with customer service (csr@powerplumobility.com)			
<input type="checkbox"/> STPAL	Articulating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin (A2)	450				
<input type="checkbox"/> STPNR	Non-removable Pin Style Legrest	N/C				
<input type="checkbox"/> STPEB	Elevating Legrest System	450				
<input type="checkbox"/> STPPPL	Plastic Footplates	N/C				
<input type="checkbox"/> STPPXL	Extra Large Footplates	150				
<input type="checkbox"/> STPPAF	Angle Adjustable Footplates	150				
<input type="checkbox"/> STP1FB	One Piece Foot Board Padded	80				
BACK TYPE	Select one below	6	<b>Please Note :</b> A1 = Client's Responsibility A2 = Clinical Justification A3 = T-Arm recommended for seat depth & angle back adjustment.			
<input type="checkbox"/> STPAB1	Angle Adjustable Back - Preset	N/C				
	<input type="checkbox"/> 0° <input type="checkbox"/> 5° <input type="checkbox"/> 10° <input type="checkbox"/> 15°					
<small>Specify Width</small> <input type="checkbox"/> STPENDO	Endomorph Canes W ( )	N/C				
<input type="checkbox"/> STPRBK	Recline Kit	560				
ARM REST	Select one below	7				
<small>Std 10"-14"</small> <input type="checkbox"/> STPTA	T-Arm Ht. Adj./ Removable (A3)	N/C				
<small>Std 10"-14"</small> <input type="checkbox"/> STPFA	Flip Back Ht. Adj./ Removable	N/C				
<small>Short 7"-11"</small> <input type="checkbox"/> STPTAS	T-Arm Ht. Adj./ Removable (A3)	N/C				
<small>Short 7"-11"</small> <input type="checkbox"/> STPFAS	Flip Back Ht. Adj./ Removable	N/C				
<input type="checkbox"/> STPGWAPF	Gel Waterfall Arm Pads (A1)	150				

HEADREST		Select one below	8
<input type="checkbox"/>	STPMAHKP	Multi Axis Complete With Small Pad	N/C
<input type="checkbox"/>	STPMAHF	Multi Axis Complete With Large Pad	N/C
<input type="checkbox"/>	STPMAHFLAT	Multi Axis Complete Flat Pad (malleable)	100

REAR WHEELS		Select one below	9	
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C	
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C	
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C	
STF 17-19"	<input type="checkbox"/>	STPRWT12I	12" Solid Tire Inside Mount	120
<input type="checkbox"/>	STPRWT16	16" Solid Tire	N/C	
<input type="checkbox"/>	STPAQR	Axles – Quick Release	120	
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C	
<input type="checkbox"/>	STPPCHR	Plastic Coated Hand Rims	N/C	

CASTER		Select one below	10
<input type="checkbox"/>	STPFC1D	3" X 1" Caster	N/C
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster (Soft Roll)	80
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	80
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	80
<input type="checkbox"/>	STPFC4B	6" X 2 Caster (Steel Fork)	80
<input type="checkbox"/>	STPFC4C	8" X 2 Caster (Steel Fork)	80

WHEEL LOCK		Select one below	11
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C
<input type="checkbox"/>	STPBRXC	4" Brake Extension	40
<input type="checkbox"/>	STPBRXA	6" Brake Extension	40
<input type="checkbox"/>	STPBRXB	9" Brake Extension	130

SEATING - BACKS		Select one below	12
<input type="checkbox"/>	STPBK01	<b>RELAX</b> Gentle Fit Back	N/C

SEATING - CUSHION		Select one below	13	
<input type="checkbox"/>	STPC01	<b>RELAX</b> Soft Cushion (Up to 250 lbs)	375	
14-20" wide	<input type="checkbox"/>	STPC02	<b>RELAX</b> Comfy Cushion (up to 250 lbs)	575
21-28" wide	<input type="checkbox"/>	STPC03	<b>RELAX</b> Firm Cushion (up to 600 lbs)	575

ACCESSORIES AND POWER OPTION		Select one below	14
<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style)	23
<input type="checkbox"/>	STPCNC	Crutch And Cane Holder (A1)	60
<input type="checkbox"/>	STPIVA	Adjustable I V Pole (A1)	130
<input type="checkbox"/>	STPNLSP	Nylon Leg Strap Padded	50
<input type="checkbox"/>	STPAMPL	Amputee Support Left	45
<input type="checkbox"/>	STPAMPR	Amputee Support Right	45
<input type="checkbox"/>	STPOXY	Oxygen Tank Holder	135
<input type="checkbox"/>	STPINTS	Independent Power Tilt System (A1)	2000
<input type="checkbox"/>	STPTDFB	Transit Tie-down Front (Non-occupied) (A4)	N/C
<input type="checkbox"/>	STPTDRB	Transit Tie-down Rear (Non-occupied) (A4)	N/C

PLUS PACKAGES	Select one below	15
<input type="checkbox"/>	<b>PLUS PACKAGE 1:</b> <b>CLIENT 250 TO 349 LBS - \$730.00</b> INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS, HEAVY DUTY STRUTS AND SEAT PAN SEAT WIDTH INCREASE FROM 16" TO 22" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 1/4" HD CASTERS	
<input type="checkbox"/>	<b>PLUS PACKAGE 2:</b> <b>CLIENT 350 TO 449 LBS - \$900.00</b> INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS, HEAVY DUTY STRUTS AND SEAT PAN SEAT WIDTH INCREASE FROM 18" TO 26" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS	
<input type="checkbox"/>	<b>PLUS PACKAGE 3:</b> <b>CLIENT 450 TO 600 LBS - \$1110.00</b> INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK, WELDED CASTER JOURNALS, HEAVY DUTY STRUTS AND SEAT PAN SEAT WIDTH INCREASE FROM 20" TO 30" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS	

70 Degree Pivot legrest are recommended for all PLUS packages. If 60 Degree or 80 Degree are required, Pin Style legrest are recommended.

**Notes:**

**Please Note :**  
**A1** = Client's Responsibility  
**A2** = Clinical Justification  
**A3** = T-Arm recommended for seat depth & angle back adjustment.  
**A4** = To be mounted on inside of frame with inside wheels option.

AADL CODE	PRODUCT	AADL PRICE	MSRP
<b>SEATING OFF CHAIR</b>			
<b>BACKS</b>			
N590	Relax Personal Fit - LESS THAN 3" LATERAL LEVEL B	\$690.00	\$690.00
N590	Relax Gentle Fit - LESS THAN 3" LATERAL LEVEL B	\$690.00	\$690.00
N593	Relax Snug Fit - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N593	Relax Full Fit - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N590	Relax Form Fit GENTLE - LESS THAN 3" LATERAL LEVEL B	\$690.00	\$690.00
N593	Relax Form Fit SNUG - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N593	Relax Form Fit FULL - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N590	Relax Spiral Fit GENTLE - LESS THAN 3" LATERAL LEVEL B	\$690.00	\$690.00
N593	Relax Spiral Fit SNUG - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N593	Relax Spiral Fit FULL - GREATER THAN 3" LATERAL LEVEL B	\$790.00	\$790.00
N500	Oversized / Heavy Duty Option for Backrests - GREATER THAN 20" WIDE LEVEL B	\$300.00	\$300.00
V853	Lumbar/Extra Pressure Reduction - Maximum 3 pieces per back LEVEL A	\$48.45	\$48.45
V854	Lateral Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45	\$48.45
V855	Pelvic Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45	\$48.45
V856	Wedge Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45	\$48.45
V901	Thoracic Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$69.00	\$69.00
V902	Full Spinal Area Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$103.50	\$103.50
V903	Full Back Area Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$137.08	\$137.08
V571	REPLACEMENT BACK COVER FOR ALL SIZES LEVEL A OR LEVEL B	\$183.00	\$183.00
<b>CUSHIONS</b>			
V980	Relax Comfy CUSHION - Full Time User Sizes Up-to 20"x20" LEVEL A	\$589.00	\$589.00
V980 + V529	Relax Firm CUSHION - Full Time User Sizes Over 20"x20" (V529 = \$190) LEVEL A	\$779.00	\$779.00
V980	Relax <b>AIRFLOW</b> CUSHION - Adjustable Positioning Up-to 20"x20" LEVEL A	\$589.00	\$589.00
V980 + V529	Relax <b>AIRFLOW</b> CUSHION - Adjustable Positioning Over 20"x20" (V529 = \$190) LEVEL A	\$779.00	\$779.00
V101	REPLACEMENT CUSHION COVER FOR ALL SIZES LEVEL A OR LEVEL B	\$165.00	\$165.00
<b>HEADRESTS</b>			
V081	STANDARD PAD - <b>F1</b> Contoured Pad Large LEVEL A	\$190.00	\$190.00
N582	HEADREST WITH LATERAL SUPPORT - <b>FXL1</b> Flex Pad Large LEVEL B	\$255.00	\$255.00
V338	STANDARD STRAIGHT HARDWARE - <b>LAH12</b> Vertical Horizontal Headrest Fixed 1/2" LEVEL A	\$165.00	\$165.00
N533	MULTI-AXIS HARDWARE - <b>MAHL</b> Multi-Axis Headrest (L-Bend) 1/2" LEVEL B	\$200.00	\$200.00
N535	GOOSENECK STYLE MULTI-AXIS - <b>MAH</b> Multi-Axis Headrest (S-Bend) 1/2" LEVEL B	\$280.00	\$280.00
V085	MOUNTING PLATE - <b>UHRA</b> Universal Headrest Plate - 1/2" OR 3/4" LEVEL A	\$189.70	\$189.70
V570	REPLACEMENT HEADREST COVER FOR ALL SIZES LEVEL A	\$120.00	\$120.00
<b>ACCESSORIES</b>			
V577	ARM TROUGH TRAY ASSEMBLY - <b>CHPTL-P</b> = Left Tray, <b>CHPTR-P</b> = Right Tray LEVEL A	\$432.00	\$432.00
V476	WHEELCHAIR TRAY - Standard or Padded, Full or Half Up-to 18" Wide LEVEL A	\$257.60	\$257.60
V542	WHEELCHAIR TRAY - Standard or Padded, Full or Half Over 18" Wide LEVEL A	\$310.20	\$310.20
V762	ELBOW STOPS LEVEL A	\$90.00	\$90.00
V766	LATERAL SUPPORT LEVEL A	\$90.00	\$90.00
V769	HAND BLOCK LEVEL A	\$90.00	\$90.00
N814	LATERAL THIGH SUPPORT LEVEL B 150	\$150.00	\$150.00
V619	MODIFIED FOOT SUPPORT - Footbox, Padded Footbox, Contracture Sling (V619 \$179.55 X 2) LEVEL A	\$359.10	\$359.10
V581	AMPUTEE BOARD - Up-to 20" Wide LEVEL A	\$350.00	\$350.00
V583	HEAVY DUTY AMPUTEE BOARD - Up-to 28" x 22" LEVEL A	\$412.50	\$412.50
V579	PAD FOR AMPUTEE BOARD LEVEL A	\$62.50	\$62.50