

INFORMATION

Customer Name: Date:
 Address: PO#:
 City: Phone:
 Account:
 Authorizer Name:

STANDARD NO CHARGE SPECIFICATIONS

LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES	5", 6", OR 8" CASTERS	PRECISION BEARINGS
STANDARD 18" BACK HEIGHT	BLACK NYLON UPHOLSTERY	MOULDED PLASTIC FOOTPLATES
20", 22", OR 24" NYLON MOULDED REAR WHEELS	POWDER COATED FRAME	SEAT WIDTH 14" – 18"
EXTRUDED ALUMINUM FORKS	HEIGHT ADJUSTABLE REMOVABLE ARMS	BLACK PADDED ARM PADS
CHROME HAND RIMS	DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS	DOUBLE CROSS BRACE

FRAME TYPE 01

MGP PLUS FRAME – Black \$1,690

SEAT WIDTH 02

MGPSW14"-20" Specify: N/C
 SW14"-20" Light Weight Aluminum Frame

MGPSW21"-32" Specify: N/C
 SW21"-28" Must Include Plus Package 1, 2 Or 3. 28"+ requires regid frame

SEAT DEPTH 03

MGPSDF18" Seat Depth 18" (By frame) N/C

MGPSD16"-18" Specify: (By upholstery) N/C

MGPSDF20" Seat Depth 20" (By frame) \$350

MGPSD18"-20" Specify: (By upholstery) N/C

MGPSDF22" Seat Depth 22" (By frame) \$350

MGPSD20"-22" Specify: (By upholstery) N/C

BACK 04

MGPBH15"-20" Specify: N/C

MGPBH21"+ Specify: \$120

MG PAB Angle Adjustable Back – Preset N/C

0° 5° 10° 15° N/C

MGPEUDO* Endomorph Canes **W ()** N/C
 Endomorph Canes require TA or TAS Arm Rest

STROLLER HANDLE 05

MGPSHS Removable Stroller Handle N/C

MGPSHR Removable Raised Stroller Handle N/C

FOOT REST 06

MGPF6 60 Degree Hanger Pivot Pin N/C

MGPF7 70 Degree Hanger Pivot Pin N/C

MGPF8 80 Degree Hanger Pivot Pin N/C

MGPEA Elevating Legrest Pivot Pin (A2) \$404

MG PAL Articulating Legrest Pivot Pin (A2) \$450

MGPNR Non-removable Pin Style Legrest N/C

FOOT REST & FOOT PLATE 07

MGPPPL Plastic N/C

MGPPAF Angle Adjustable – Extra Large \$150

MGPPXL Extra Large – Alum (F6A,F7A,F8A) \$150

MGPIFB One Piece Foot Board \$80

BACK TYPE 08

MG PFA Flip Back Ht. Adj./ Removable \$240

MGPTA T-Arm Ht. Adj./ Removable (A3) \$240

MGPFAS Flip Back (Short) Ht. Adj./ Removable \$240

MGPTAS T-Arm Ht. Adj./ Removable (A3) \$240

MGPGWAPS Gel Waterfall Arm Pads (A1) \$150

SEAT TO FLOOR SPECIFICATIONS (± 0.5")

(LOW SEAT SETTING)

STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 12.25"	20"	5"	<input type="checkbox"/> 2 HOLE
<input type="checkbox"/> 13.00"	20"	5"	<input type="checkbox"/> 2 HOLE
<input type="checkbox"/> 13.75"	20"	5"	<input type="checkbox"/> 2 HOLE
<input type="checkbox"/> 14.50"	20"	6"	<input type="checkbox"/> 3 HOLE
<input type="checkbox"/> 14.50"	22"	5"	<input type="checkbox"/> 3 HOLE
<input type="checkbox"/> 14.50"	22"	6"	<input type="checkbox"/> 3 HOLE
<input type="checkbox"/> 15.25"	20"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	20"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	22"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	22"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 16.00"	22"	6"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 16.75"	22"	6"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 17.50"	22"	6"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 17.50"	22"	8"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 17.50"	22"	6"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 17.50"	24"	8"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 18.25"	24"	8"	<input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 19.00"+	24"	8"	<input type="checkbox"/> 5 HOLE

NOTE: *TT 19.00' requires taller caster journal

REAR WHEELS	Select one below	09
<input type="checkbox"/> MGPRWT0	20" MAG Solid Treaded Tire	N/C
<input type="checkbox"/> MGPRWT2	22" MAG Solid Treaded Tire	N/C
<input type="checkbox"/> MGPRWT4	24" MAG Solid Treaded Tire	N/C
<input type="checkbox"/> MGPRWP2	22" Pyramid Tire	\$80
<input type="checkbox"/> MGPRWP4	24" Pyramid Tire	\$80
<input type="checkbox"/> MGPRWA2	22" Air Tire	\$120
<input type="checkbox"/> MGPRWA4	24" Air Tire	\$120
<input type="checkbox"/> MGPAQR	Axles – Quick Release (349 Lbs)	N/C
<input type="checkbox"/> MGPPCHR	Plastic Coated Hand Rims	N/C
<input type="checkbox"/> MGPOADS	One Arm Drive System <input type="checkbox"/> Left <input type="checkbox"/> Right	\$1,100

Available sizes 16" – 24" with 22" or 24" wheel only

CASTERS	Select one below	10
<input type="checkbox"/> MGPFCA1	5" X 1.25" Caster	N/C
<input type="checkbox"/> MGPFCA1B	6" X 1.25" Caster	N/C
<input type="checkbox"/> MGPFCA1C	8" X 1.25" Caster	N/C
<input type="checkbox"/> MGPFCA2	5" X 1.75" Caster	\$80
<input type="checkbox"/> MGPFCA2B	6" X 1.75" Caster	\$80
<input type="checkbox"/> MGPFCA2C	8" X 1.75" Caster	\$80
<input type="checkbox"/> MGPFCA3A	5" X 1.25" Air Caster	\$80
<input type="checkbox"/> MGPFCA3B	6" X 1.25" Air Caster	\$80
<input type="checkbox"/> MGPFCA3C	8" X 1.25" Air Caster	\$80

WHEEL LOCK	Select one below	11
<input type="checkbox"/> MGPALPL	Aluminium Pull To Lock	N/C
<input type="checkbox"/> MGPALPH	Aluminium Push To Lock	N/C
<input type="checkbox"/> MGPSTPL	Steel Pull To Lock	N/C
<input type="checkbox"/> MGPSTPH	Steel Push To Lock	N/C
<input type="checkbox"/> MGPBAXA	6" Brake Extension	\$40
<input type="checkbox"/> MGPBAXB	9" Brake Extension	\$60

SEATING BENIFITS AND ACCESSORIES

Kindly refer to new basic seating benefit order form

PLUS PACKAGES	Select one below	14
<input type="checkbox"/> MGP PLUS PACKAGE 1:	CLIENT 250 – 349 LBS - \$880.00 INCLUDES THE FOLLOWING: REMOVABLE STROLLER HANDLE, STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 16" TO 24" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/> MGP PLUS PACKAGE 2:	LIENT 350 – 449 LBS - \$1,140.00 INCLUDES THE FOLLOWING: REMOVABLE STROLLER HANDLE, STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 18" TO 26" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS HD STEEL FORKS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/> MGP PLUS PACKAGE 3:	CLIENT 450 – 600 LBS - \$1,300.00 INCLUDES THE FOLLOWING: REMOVABLE STROLLER HANDLE, STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY SEAT WIDTH INCREASE FROM 20" TO 28" SEAT DEPTH INCREASE FROM 16" TO 22" 5,6 OR 8 X 1 ¼" HD CASTERS HD STEEL FORKS *** WELDED CASTER JOURNAL ***	
<input type="checkbox"/> FIXED REINFORCED ONE PIECE BACK POST	- \$260.00 (SPECIFY ANGLE: _____)	
<input type="checkbox"/> HEAVY DUTY TENSION ADJ. BACK UPH	- \$220.00	
<input type="checkbox"/> HEAVY DUTY SOLID SEA T INSERT (V556)	- \$143.00 Please specify: <input type="checkbox"/> IN SEAT UPHOLESTERY POCKET <input type="checkbox"/> IN CUSHION COVER	

PLEASE NOTE:

- A1 =** Client's Cost
 - A2 =** Clinical Justification
 - A3 =** T-Arm recommended for seat depth & angle back adjustment.
 - A4 =** To be mounted on inside of frame with inside wheels option.
- For more Seat To Floor combinations, visit our website**
<https://powerplusmobility.com/stf>
Please confirm with customer service (csr@powerplusmobility.com)

NOTES:

INFORMATION

Customer Name: Date:

Address: PO#:

City: Phone:

Account:

Authorizer Name:

BACKRESTS		O1
<input type="checkbox"/>	N590 Relax Form Fit (No Lateral)	\$690
<input type="checkbox"/>	N590 Relax Gentle Fit (<3" Lateral)	\$690
<input type="checkbox"/>	N590 Relax Gentle Form Fit (<3" Lateral)	\$690
<input type="checkbox"/>	N590 Relax Gentle Spiral Fit (<3" Lateral)	\$690
<input type="checkbox"/>	N593 Relax Snug Form Fit (>3" Lateral)	\$790
<input type="checkbox"/>	N593 Relax Snug Spiral Fit (>3" Lateral)	\$790
<input type="checkbox"/>	N593 Relax Full Form Fit (>3" Lateral)	\$790
<input type="checkbox"/>	N593 Relax Full Spiral Fit (>3" Lateral)	\$790
<input type="checkbox"/>	N593 Relax Snug Fit (>3" Lateral)	\$790
<input type="checkbox"/>	N593 Relax Full Fit (>3" Lateral)	\$790
Dimensions (in inches):		
Width: _____ x Height: _____		

ADDITIONAL OPTIONS		O2
<input type="checkbox"/>	N500 Heavy-Duty Backrest (>20" Wide)	\$300
<input type="checkbox"/>	V571 Replacement Back Cover	\$183

PRESSURE REDUCTION OPTIONS		O3
Limited to 3 pieces per backrest.		
<input type="checkbox"/>	V853 Lumbar Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/>	V854 Lateral Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/>	V855 Pelvic Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/>	V856 Wedge Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/>	V901 Thoracic Gel Pressure Reduction QTY: ____ Location: ____	\$69
<input type="checkbox"/>	V902 Full Spinal Gel Pressure Reduction QTY: ____ Location: ____	\$103.50
<input type="checkbox"/>	V903 Full Back Gel Pressure Reduction QTY: ____ Location: ____	\$137.08

CUSHIONS		O4
<input type="checkbox"/>	V980 Relax Comfy (Up to 20"x20")	\$690
<input type="checkbox"/>	V980HD Relax Firm (Over 20"x20")	\$690
<input type="checkbox"/>	V980 Relax AIRFLOW (Up to 20"x20")	\$790
<input type="checkbox"/>	V980HD Relax AIRFLOW (Over 20"x20")	\$790
Dimensions (in inches):		
Width: _____ x Depth: _____		

ADDITIONAL OPTIONS		O5
<input type="checkbox"/>	V101 Replacement Cushion Cover	\$165
<input type="checkbox"/>	VO06 Incontinence Cover	\$135

HEADREST PAD		O7
<input type="checkbox"/>	V081 Standard Pad	\$190
<input type="checkbox"/>	N582 Headrest with Lateral Support	\$255

HEADREST HARDWARE		O8
<input type="checkbox"/>	V338 Standard Straight Hardware	\$165
<input type="checkbox"/>	N533 Multi-Axis Hardware (L-Bend)	\$200
<input type="checkbox"/>	N535 Gooseneck Multi-Axis (S-Bend)	\$280
<input type="checkbox"/>	V085 Mounting Plate <input type="checkbox"/> 1/2" or <input type="checkbox"/> 3/4"	\$189.70
<input type="checkbox"/>	V570 Replacement Headrest Cover	\$120

ACCESSORIES		O6
<input type="checkbox"/>	V577 Arm Trough Tray Assembly	\$432
<input type="checkbox"/>	V476 Wheelchair Tray (Up to 18" Wide)	\$257.60
<input type="checkbox"/>	V542 Wheelchair Tray (Over 18" Wide)	\$310.20
<input type="checkbox"/>	V762 Elbow Stops	\$90
<input type="checkbox"/>	V766 Lateral Support	\$90
<input type="checkbox"/>	V769 Hand Block	\$90
<input type="checkbox"/>	N814 Lateral Thigh Support	\$150
<input type="checkbox"/>	V619 Modified Foot Support (Padded Footbox)	\$359.10
<input type="checkbox"/>	V581 Amputee Board (Up to 20" Wide)	\$350
<input type="checkbox"/>	V583 Heavy Duty Amputee Board	\$412.50
<input type="checkbox"/>	V579 Pad for Amputee Board	\$62.50
<input type="checkbox"/>	V572 Modular Arm Support (Left Side)	\$62.50
<input type="checkbox"/>	V572 Modular Arm Support (Right Side)	\$62.50

NOTES:

