

EXT Extreme Wheelchair (Tilt Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

- EXTB Dynamic Frame - Black 3899
- EXTC Dynamic Frame - Colour : _____ 4299

CLIENT WEIGHT *Select one below* **2**

- EXTWA Up to 250Lbs N/C
- EXTWAMF Exceeds 250Lbs 730
- EXTWAMG Exceeds 350Lbs 992
- EXTWAMH Exceeds 400Lbs (Max 600Lbs) 1110

SEAT WIDTH *Select one below* **3**

- EXTSW16 Seat Width 16" N/C
- EXTSW17 Seat Width 17" N/C
- EXTSW18 Seat Width 18" N/C
- EXTSW19 Seat Width 19" 399
- EXTSW20 Seat Width 20" 399
- EXTSWC Custom Seat Width: _____ 399

SEAT DEPTH *Select one below* **4**

- EXTSD16 Seat Depth 16" N/C
- EXTSD17 Seat Depth 17" N/C
- EXTSD18 Seat Depth 18" N/C
- EXTSD19 Seat Depth 19" 265
- EXTSD20 Seat Depth 20" 265
- EXTSDC Custom Seat Depth: _____ 265

FOOT REST AND FOOT PLATES *Select one below* **5**

- EXTF6A 60 Degree Hanger Pivot Pin N/C
- EXTF7A 70 Degree Hanger Pivot Pin N/C
- EXTF8A 80 Degree Hanger Pivot Pin N/C
- EXTEA Elevating Legrest Pivot Pin 404
- EXTAL Articulating Leg Rests Pivot Pin 450
- EXTEB Elevating Legrest System 450
- EXTPPL Plastic N/C
- EXTPXL Extra Large 150
- EXTPAF Angle Adjustable Footplates 240
- EXTFRCH Padded Footplates 88
- EXT1FB One Piece Foot Board 80
- EXT1FBP One Piece Foot Board Padding 44
- EXTNLS Nylon Leg Strap Padded 60
- EXTNLSC Nylon Leg Strap Custom **Size ()** 120
- EXTCSL Contracture Sling Standard M Style 240
- EXTOFST Offset Foot Rest Receivers N/C

STANDARD NO CHARGE SPECIFICATIONS

- SEAT WIDTH: 14" – 18"
- SEAT DEPTH (ADJUSTABLE): 16" – 18"
- SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
- SEAT ANGLE: 0 – 30 [DEGREES]
- BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
- ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
- WEIGHT CAPACITY: 250 LBS
- WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
- DYNAMIC TILT ADJUSTMENTS 0 TO 25 (DEGREES)
- TWO-WAY LOCKING GAS STRUT
- SINGLE RELEASE, DUAL CYLINDERS
- STANDARD LOCKOUT COLLARS
- REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

STF	R. WHEEL	CASTER	FORK
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT
<input type="checkbox"/> L 14.50"	20"	6"	SHORT
<input type="checkbox"/> L 14.50"	22"	5"	SHORT
<input type="checkbox"/> L14.50"	22"	6"	SHORT
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/> SH 16.00"	22"	6"	LONG

(STANDARD SEAT SETTING)

<input type="checkbox"/> MH 16.75"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	8"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	24"	6"	LONG
<input type="checkbox"/> A 18.25"	24"	8"	LONG

Notes:

BACK TYPE		6	SEATING - BACKS <i>For custom, choose one back then fill custom</i>		12		
<input type="checkbox"/>	EXTAB1	Angle Adjustable Back	N/C	W () x H ()			
<i>Specify Width</i> <input type="checkbox"/>	EXTENDO	Endomorph Canes	W ()	<input type="checkbox"/>	EXTBK01	RELAX Gentle Fit Back	544
<input type="checkbox"/>	EXTRBK	Recline Kit	560	<input type="checkbox"/>	EXTBK02	RELAX Full Fit Back	608
HEADREST <i>Select one below</i>			7	<input type="checkbox"/>	EXTBK03	RELAX Snug Fit Back	700
<input type="checkbox"/>	EXTMAHF	Multi-Axis Hardware Complete with Large Pad	512	<input type="checkbox"/>	EXTBK04	RELAX Form Fit Back No Laterals	820
<input type="checkbox"/>	EXTMAHFSL	Multi-Axis Hardware Complete with Flex Large	512	<input type="checkbox"/>	EXTBK06	RELAX Form Fit Back Gentle Laterals	820
<input type="checkbox"/>	EXTMAHKP	Multi-Axis Hardware Complete with Medium Pad	512	<input type="checkbox"/>	EXTBK07	RELAX Form Fit Back Full Laterals	820
<input type="checkbox"/>	EXTMAHFXS	Multi-Axis Hardware Complete with Flex Medium	512	<input type="checkbox"/>	EXTBK05	RELAX Form Fit Back Snug Laterals	820
<input type="checkbox"/>	EXTMAHFT	Multi-Axis Hardware Complete with Small Flat	512	<input type="checkbox"/>	EXTBK10	RELAX Spiral Fit Back Gentle Laterals	832
<input type="checkbox"/>	EXTUPHT219	Contour Headrest Complete with Swivel Hardware	433	<input type="checkbox"/>	EXTBK09	RELAX Spiral Fit Back Snug Laterals	832
<i>(1/2" Tubing)</i> <input type="checkbox"/>	EXTUPHT218	Contour Headrest with Std. Assembly (1/2")	292	<input type="checkbox"/>	EXTBK11	RELAX Spiral Fit Back Full Laterals	832
<i>(3/4" Tubing)</i> <input type="checkbox"/>	EXTUPHT217	Contour Headrest with Std. Assembly (3/4")	292	<input type="checkbox"/>	EXTTAB	RELAX Tension Back Upto 18"	140
<input type="checkbox"/>	EXTUHA	Universal Hardware Adaptor (Adaptor only)	57	<input type="checkbox"/>	EXTTAB+	RELAX Tension Back Over 18"	240
<input type="checkbox"/>	EXTMAHGA-P	Multi Axis with Headrest Plate (No Pad)	255	<input type="checkbox"/>	EXTBGL1	Half Overlay Gel (3 Piece)	348
ARM REST <i>Select one below</i>			8	<input type="checkbox"/>	EXTBGL2	Full Overlay Gel (6 Pieces)	696
<i>(Std 10"-14")</i> <input type="checkbox"/>	EXTTA	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	EXTPEXT	Pin Extension (1" Extension)	40
<input type="checkbox"/>	EXTGWAPS	* Gel Waterfall Arm Pads	180	<input type="checkbox"/>	EXTBK12	* Custom Width : _____	176
<i>*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required</i>				<input type="checkbox"/>	EXTBK13	* Custom Height : _____	176
REAR WHEELS <i>Select one below</i>			9	<i>*Note: Please refer to the Seating Booklet for custom sizes</i>			
<input type="checkbox"/>	EXTRWT0	20" Solid Treaded Tire	N/C	SEATING - CUSHION <i>Select one below</i>			13
<input type="checkbox"/>	EXTRWT2	22" Solid Treaded Tire	N/C	W () x D ()			
<input type="checkbox"/>	EXTRWT4	24" Solid Treaded Tire	N/C	<input type="checkbox"/>	EXTC01	RELAX Soft Cushion	480
<i>(STF 13-19")</i> <input type="checkbox"/>	EXTRWT12	12" Solid Tire	N/C	<input type="checkbox"/>	EXTC02	RELAX Comfy Cushion	624
<input type="checkbox"/>	EXTRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	EXTC03	RELAX Firm Cushion	784
<input type="checkbox"/>	EXTRWA2	22" Air Tire	120	<input type="checkbox"/>	EXTCGL1	Half Overlay Gel (1 Piece)	116
<input type="checkbox"/>	EXTRWA4	24" Air Tire	120	<input type="checkbox"/>	EXTCGL2	Full Overlay Gel (2 Pieces)	232
<input type="checkbox"/>	EXTAQR	Axles – Quick Release	120	<input type="checkbox"/>	EXTCCL	Contoured Leg Channel (Bolster) Std.	180
<input type="checkbox"/>	EXTHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	EXTCCLA	Contoured Leg Channel (Bolster) with ABS	260
CASTER <i>Select one below</i>			10	<input type="checkbox"/>	EXTCCLC	* Contoured Leg Channel (Bolster) Custom	320
<input type="checkbox"/>	EXTFC1A	5" X 1.25" Caster	N/C	<i>*Note: Please refer to the Seating Booklet for custom sizes</i>			
<input type="checkbox"/>	EXTFC1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	EXTC00	Extra Cushion Cover	99
<input type="checkbox"/>	EXTFC1C	8" X 1.25" Caster	N/C	ACCESSORIES			14
<input type="checkbox"/>	EXTFC2A	5" X 1.75" Caster (Soft Roll)	80	<input type="checkbox"/>	EXTSBA	Seat Belt (Auto Style)	26
<input type="checkbox"/>	EXTFC2B	6" X 1.75" Caster	80	<input type="checkbox"/>	EXTSBAP	Seat Belt (Padded)	56
<input type="checkbox"/>	EXTFC2C	8" X 1.75" Caster	80	<input type="checkbox"/>	EXTSBA4P	Seat Belt (4 Point Padded)	130
<input type="checkbox"/>	EXTFC4B	6" X 2 Caster	80	<input type="checkbox"/>	EXTTRA	Full Padded Tray Locking	280
<input type="checkbox"/>	EXTFC4C	8" X 2 Caster	80	<input type="checkbox"/>	EXTTRB	Full Padded Tray Slide On	280
WHEEL LOCK <i>Select one below</i>			11	<input type="checkbox"/>	EXTTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/>	EXTALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	EXTTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/>	EXTALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	EXTTRD	Half Padded Left Tray	240
<input type="checkbox"/>	EXTEXTH	Steel Push To Lock	N/C	<input type="checkbox"/>	EXTTRE	Half Padded Right Tray	240
<input type="checkbox"/>	EXTEXTL	Steel Pull To Lock	N/C	<input type="checkbox"/>	EXTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	EXTBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	EXTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	EXTBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	EXTTDFB	Transport Tie Downs Front	85
<input type="checkbox"/>	EXTBRXB	9" Brake Extension	130	<input type="checkbox"/>	EXTTDRB	Transport Tie Downs Rear	85
<input type="checkbox"/>	EXTATTBRX	Attendant Brakes	240	<input type="checkbox"/>	EXTOXY	Oxygen Tank Holder	195
				<input type="checkbox"/>	EXTIVA	Adjustable I V Pole	130
				<input type="checkbox"/>	EXTCHPT-L	Left Arm Trough Tray	400
				<input type="checkbox"/>	EXTCHPT-R	Right Arm Trough Tray	400
				<input type="checkbox"/>	EXTLATB	Lateral Bite out	180
				<input type="checkbox"/>	EXTADD8	8" Adductor	164
				<input type="checkbox"/>	EXTADD10	10" Adductor	164