

EXT Extreme Wheelchair - Tilt Wheelchair



INFORMATION

Customer Name: _____
Address: _____ PO #: _____
City: _____ Phone #: _____ Date: _____
Account #: _____ Authorizer Name: _____

FRAME TYPE	Select one below	1
☐ EXTWANA	Dynamic Tilt Frame	\$1,440
FRAME COLOUR	Select one below	2
☐ EXTB	Frame Color - Black	\$3,239
☐ EXTS	Frame Color - Silver	\$3,239
☐ EXTB	Frame Color - Blue	\$3,239
☐ EXTB	Frame Color - Burgundy	\$3,239
☐ EXTC	Frame - Colour : _____	\$3,599
CLIENT WEIGHT	Select one below	3
☐ EXTWA	Up to 250Lbs	N/C
☐ EXTWAMF	Exceeds 250Lbs	\$876
☐ EXTWAMG	Exceeds 350Lbs	\$1,190
☐ EXTWAMH	Exceeds 400Lbs (Max 600Lbs)	\$1,332
SEAT WIDTH	Select one below	4
☐ EXTSW16	Seat Width 16"	N/C
☐ EXTSW17	Seat Width 17"	N/C
☐ EXTSW18	Seat Width 18"	N/C
☐ EXTSW19	Seat Width 19"	\$479
☐ EXTSW20	Seat Width 20"	\$479
☐ EXTSWC	Custom Seat Width: _____	\$479
SEAT DEPTH	Select one below	5
☐ EXTSD16	Seat Depth 16"	N/C
☐ EXTSD17	Seat Depth 17"	N/C
☐ EXTSD18	Seat Depth 18"	N/C
☐ EXTSD19	Seat Depth 19"	\$318
☐ EXTSD20	Seat Depth 20"	\$318
☐ EXTSDC	Custom Seat Depth: _____	\$318
FOOT REST AND FOOT PLATES	Select one below	6
☐ EXTF6A	60 Degree Hanger ☐ Pivot ☐ Pin	N/C
☐ EXTF7A	70 Degree Hanger ☐ Pivot ☐ Pin	N/C
☐ EXTF8B	80 Degree Hanger ☐ Pin	N/C
☐ EXTEA	Elevating Legrest ☐ Pivot ☐ Pin	\$485
☐ EXTAL	Articulating Leg Rests ☐ Pivot ☐ Pin	\$540
☐ EXTEB	Elevating Legrest System	\$540
☐ EXTPPL	Plastic Footplates	N/C
☐ EXTPXL	Extra Large Footplates	\$180
☐ EXTPAF	Angle Adjustable Footplates	\$288
☐ EXT1FB	One Piece Foot Board	\$96
☐ EXT1FBP	One Piece Foot Board Padding	\$53
☐ EXTNLS	Nylon Leg Strap Padded	\$72
☐ EXTNLSC	Nylon Leg Strap Custom Size ()	\$144
☐ EXTCSL	Contracture Sling ☐ Standard ☐ M Style	\$288
☐ EXTOFST	Offset Foot Rest Receivers	N/C

*Note: Pin Style elevating legrest is not available for 14" Seat to Floor or lower

BACK TYPE	Select one below	7
☐ EXTAB1	Angle Adjustable Back	N/C
☐ EXTENDO	Endomorph Canes W ()	N/C
☐ EXTRBK	Recline Kit	\$672

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE): 16" – 18"
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
DYNAMIC TILT ADJUSTMENTS 0 TO 25 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CABLES
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
☐	12.25"	20"	5"	XL SHORT
☐	13.00"	20"	5"	XL SHORT
☐	13.75"	20"	5"	XL SHORT
☐	14.50"	20"	6"	SHORT
☐	14.50"	22"	5"	SHORT
☐	14.50"	22"	6"	SHORT
☐	15.25"	20"	5"	SHT/LONG
☐	15.25"	20"	6"	SHT/LONG
☐	15.25"	22"	5"	SHT/LONG
☐	15.25"	22"	6"	SHT/LONG
☐	16.00"	22"	6"	LONG
(HIGH SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
☐	16.75"	22"	6"	LONG
☐	17.50"	22"	6"	LONG
☐	17.50"	22"	6"	LONG
☐	17.50"	22"	8"	LONG
☐	17.50"	24"	8"	LONG
☐	18.25"	24"	6" or 8"	LONG

For more Seat To Floor combinations, visit our website

HEADREST		Select one below	8	SEATING - BACKS		Custom: select back & fill custom	13	
☐	EXTMAHF	Multi-Axis Hardware Complete with Large Pad	\$614	☐	EXTBK01	W () x H () RELAX Gentle Fit Back	\$653	
☐	EXTMAHFXL	Multi-Axis Hardware Complete with Flex Large	\$614	☐	EXTBK02	RELAX Full Fit Back	\$730	
☐	EXTMAHKP	Multi-Axis Hardware Complete with Medium Pad	\$614	☐	EXTBK03	RELAX Snug Fit Back	\$840	
☐	EXTMAHFXS	Multi-Axis Hardware Complete with Flex Medium	\$614	☐	EXTBK04	RELAX Form Fit Back No Laterals	\$984	
☐	EXTMAHFT	Multi-Axis Hardware Complete with Small Flat	\$614	☐	EXTBK06	RELAX Form Fit Back Gentle Laterals	\$984	
☐	EXTUPHT219	Contour Headrest Complete with Swivel Hardware	\$520	☐	EXTBK07	RELAX Form Fit Back Full Laterals	\$984	
☐	EXTUPHT218	Contour Headrest with Std. Assembly (1/2" Tube)	\$350	☐	EXTBK05	RELAX Form Fit Back Snug Laterals	\$984	
☐	EXTUPHT217	Contour Headrest with Std. Assembly (3/4" Tube)	\$265	☐	EXTBK10	RELAX Spiral Fit Back Gentle Laterals	\$998	
☐	EXTUHRA	Universal Hardware Adaptor (Adaptor only)	\$68	☐	EXTBK09	RELAX Spiral Fit Back Snug Laterals	\$998	
ARM REST		Select one below	9	☐	EXTBK11	RELAX Spiral Fit Back Full Laterals	\$998	
(Ht. 10"-14")	☐	EXTTA	T-Arm Ht. Adj./ Removable	\$288	☐	EXTTAB	RELAX Tension Back Upto 18"	\$168
(Ht. 7"-10")	☐	EXTTAS	T-Arm Ht. Adj./ Removable Short	\$288	☐	EXTTAB+	RELAX Tension Back Over 18"	\$288
	☐	EXTGWAPF	* Gel Waterfall Arm Pads	\$216	☐	EXTBK12	* Custom Width :	\$211
* Note: Gel Waterfall Arm Pads are only compatible with locking laptrays and full laptrays if required				☐	EXTBK13	* Custom Height :	\$211	
REAR WHEEL		Select one below	10	*Note: Please refer to the Seating Booklet for custom sizes				
☐	EXTRWT0	20" Solid Treaded Tire	N/C	☐	EXTBGL1	Gel Full Overlay	\$696	
☐	EXTRWT2	22" Solid Treaded Tire	N/C	☐	EXTBGL2	Gel Center Vertebrae	\$418	
☐	EXTRWT4	24" Solid Treaded Tire	N/C	☐	EXTBGL3	Gel Left Lateral	\$139	
(STF 13-19")	☐	EXTRWT12	* 12" Solid Tire	N/C	☐	EXTBGL4	Gel Right Lateral	\$139
☐	EXTRWT16	* 16" Solid Tire	N/C	SEATING - CUSHIONS				
☐	EXTRWA2	22" Air Tire	\$144			Select one below	14	
☐	EXTRWA4	24" Air Tire	\$144			W () x D ()		
☐	EXTAQR	Axles – Quick Release	\$144	☐	EXTC01	RELAX Soft Cushion	\$576	
☐	EXTAQR	Standard Handrim Config. (3/4")	\$144	☐	EXTC02	RELAX Comfy Cushion	\$749	
☐	EXTHRN	Narrow Handrim Config. (1/2" less)	N/C	☐	EXTC03	RELAX Firm Cushion	\$941	
☐	EXTHRT	Tight handrim Config. (No space)	N/C	☐	EXTC04	RELAX AIR Flow Cushion	\$704	
*Note: Quick Release is not available for wheels smaller than 20"				☐	EXTCGL2	Full Overlay Gel	\$430	
CASTER		Select one below	11	☐	EXTCCL	Contoured Leg Channel (Bolster) Std.	\$216	
☐	EXTFC1A	5" X 1.25" Caster	N/C	☐	EXTCCLA	Contoured Leg Channel (Bolster) with ABS	\$312	
☐	EXTFC1B	6" X 1.25" Caster	N/C	☐	EXTCCCL	* Contoured Leg Channel (Bolster) Custom	\$384	
☐	EXTFC1C	8" X 1.25" Caster	N/C	*Note: Please refer to the Seating Booklet for custom sizes				
☐	EXTFC2A	5" X 1.75" Caster (Soft Roll)	\$96	☐	EXTC00	Extra Cushion Cover	\$119	
☐	EXTFC2B	6" X 1.75" Caster	\$96	ACCESSORIES				
☐	EXTFC2C	8" X 1.75" Caster	\$96	☐	EXTSBA	Seat Belt (Auto Style)	\$31	
☐	EXTFC4B	6" X 2 Caster (Requires Steel Fork)	\$96	☐	EXTSBAP	Seat Belt (Padded)	\$67	
☐	EXTFC4C	8" X 2 Caster (Requires Steel Fork)	\$96	☐	EXTSBA4P	Seat Belt (4 Point Padded)	\$156	
WHEEL LOCK		Select one below	12	☐	EXTTRA	Full Padded Tray Locking	\$336	
☐	EXTALPH	Aluminium Push To Lock	N/C	☐	EXTTRB	Full Padded Tray Slide On	\$336	
☐	EXTALPL	Aluminium Pull To Lock	N/C	☐	EXTTRCL	Full Plastic Tray Locking ☐ ABS Black ☐ Clear	\$336	
☐	EXTSTPH	Steel Push To Lock	N/C	☐	EXTTRC	Full Plastic Tray Slide On ☐ ABS Black ☐ Clear	\$336	
☐	EXTSTPL	Steel Pull To Lock	N/C	☐	EXTTRD	Half Padded Left Tray	\$288	
☐	EXTBRXC	4" Brake Extension	N/C	☐	EXTTRE	Half Padded Right Tray	\$288	
☐	EXTBRXA	6" Brake Extension	N/C	☐	EXTTRF	Half Left Tray ☐ ABS Black ☐ Clear	\$288	
☐	EXTBRXB	9" Brake Extension	\$156	☐	EXTTRG	Half Right Tray ☐ ABS Black ☐ Clear	\$288	
☐	EXTATTBRX	Attendant Brakes	\$288	☐	EXTTDFB	Transport Tie Downs Front	\$102	
☐	EXTULB	Unilateral Hand Brake ☐ Left ☐ Right	\$144	☐	EXTTDRB	Transport Tie Downs Rear	\$102	
				☐	EXTOXY	Oxygen Tank Holder	\$234	
				☐	EXTIVA	Adjustable I V Pole	\$156	
				☐	EXTCHPT-L	Left Arm Trough Tray	\$480	
				☐	EXTCHPT-R	Right Arm Trough Tray	\$480	
				☐	EXTLATB	Lateral Bite out	\$216	
				☐	EXTADD8	8" Adductor	\$197	
				☐	EXTADD10	10" Adductor	\$197	
				☐	EXTADD12	12" Adductor	\$197	
				☐	EXTFRCH	Padded Footplates	\$106	