

GLY Glyder Wheelchair (Manual Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

- GLYS Frame - Silver 2390
- GLYB Frame - Black 2390
- GLYC Frame - Colour : _____ 2790

CLIENT WEIGHT *Select one below* **2**

- GLYWA Upto 250Lbs N/C
- GLYWAMF Exceeds 250Lbs 730
- GLYWAMG Exceeds 350Lbs 992
- GLYWAMH Exceeds 400Lbs (Max 600Lbs) 1110

SEAT WIDTH *Select one below* **3**

- GLYSW14 Seat Width 14" N/C
- GLYSW15 Seat Width 15" N/C
- GLYSW16 Seat Width 16" N/C
- GLYSW17 Seat Width 17" N/C
- GLYSW18 Seat Width 18" N/C
- GLYSW19 Seat Width 19" 399
- GLYSW20 Seat Width 20" 399
- GLYSWS Custom Seat Width: _____ 399

SEAT DEPTH *Select one below* **4**

- GLYSD14 Seat Depth 14" 350
- GLYSD15 Seat Depth 15" 350
- GLYSD16 Seat Depth 16" N/C
- GLYSD17 Seat Depth 17" N/C
- GLYSD18 Seat Depth 18" N/C
- GLYSD19 Seat Depth 19" 265
- GLYSD20 Seat Depth 20" 265
- GLYSDS Custom Seat Depth: _____ 265

BACK HEIGHT *Select one below* **5**

- GLYBH15 Back Height 15" N/C
- GLYBH16 Back Height 16" N/C
- GLYBH17 Back Height 17" N/C
- GLYBH18 Back Height 18" N/C
- GLYBH19 Back Height 19" N/C
- GLYBH20 Back Height 20" N/C
- GLYBH21 Back Height 21" N/C
- GLYBHC Back Height Custom : _____ 120

Notes:

SPECIFICATIONS

- LIGHT WEIGHT ALUMINIUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- VERTICAL/HORIZONTAL REAR AXLE ADJUSTABILITY
- ANTI-TIP ROLLERS
- 6" BRAKE EXTENSIONS
- DOUBLE CROSS BRACE
- REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> SL 12.75"	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L 14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	

(STANDARD SEAT SETTING)

<input type="checkbox"/> MH 16.75	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	8"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	24"	8"	LONG
<input type="checkbox"/> A 18.25"	24"	8"	LONG

BACK TYPE *Select one below* **6**

<input type="checkbox"/> GLYAB1	Angle Adjustable Back	N/C
---------------------------------	-----------------------	-----

FOOT REST AND FOOT PLATES *Select one below* **7**

<input type="checkbox"/> GLYF6A	60 Degree Hanger	N/C
<input type="checkbox"/> GLYF7A	70 Degree Hanger	N/C
<input type="checkbox"/> GLYF8A	80 Degree Hanger	N/C
<input type="checkbox"/> GLYPPL	Plastic	N/C
<input type="checkbox"/> GLYPXL	Extra Large	150
<input type="checkbox"/> GLYPAFE	Angle Adjustable	240
<input type="checkbox"/> GLYEA	Elevating Legrest	404
<input type="checkbox"/> GLYAL	Articulating Leg Rests	450
<input type="checkbox"/> GLY1FB	One Piece Foot Board	80
<input type="checkbox"/> GLYNLS	Nylon Leg Strap Padded	60
<input type="checkbox"/> GLYCSL	Contracture Sling	240

ARM REST *Select one below* **8**

(Std 10"-14") <input type="checkbox"/> GLYFA	Flip Back Ht. Adj./ Removable	240
(Std 10"-14") <input type="checkbox"/> GLYTA	T-Arm Ht. Adj./ Removable	340
(Short 7"-11") <input type="checkbox"/> GLYFAS	Flip Back Ht. Adj./ Removable	240
(Short 7"-11") <input type="checkbox"/> GLYTAS	T-Arm Ht. Adj./ Removable	340
<input type="checkbox"/> GLYGWAPS	Gel Waterfall Arm Pads	180

REAR WHEELS *Select one below* **9**

<input type="checkbox"/> GLYRWT0	20" Solid Treaded Tire	N/C
<input type="checkbox"/> GLYRWT2	22" Solid Treaded Tire	N/C
<input type="checkbox"/> GLYRWT4	24" Solid Treaded Tire	N/C
<input type="checkbox"/> GLYRWS0	20" Spoked Treaded Tire	N/C
<input type="checkbox"/> GLYRWS2	22" Spoked Treaded Tire	N/C
<input type="checkbox"/> GLYRWS4	24" Spokedtreaded Tire	N/C
<input type="checkbox"/> GLYRWA2	22" Air Tire	120
<input type="checkbox"/> GLYRWA4	24" Air Tire	120
<input type="checkbox"/> GLYAQR	Axles – Quick Release	120

CASTER *Select one below* **10**

<input type="checkbox"/> GLYFC1A	5" X 1.25" Caster	N/C
<input type="checkbox"/> GLYFC1B	6" X 1.25" Caster	N/C
<input type="checkbox"/> GLYFC1C	8" X 1.25" Caster	N/C
<input type="checkbox"/> GLYFC2A	5" X 1.75" Caster (Soft Roll)	80
<input type="checkbox"/> GLYFC2B	6" X 1.75" Caster	80
<input type="checkbox"/> GLYFC2C	8" X 1.75" Caster	80
<input type="checkbox"/> GLYFC4B	6" X 2 Caster	80
<input type="checkbox"/> GLYFC4C	8" X 2 Caster	80

WHEEL LOCK *Select one below* **11**

<input type="checkbox"/> GLYALPH	Aluminium Push To Lock	N/C
<input type="checkbox"/> GLYALPL	Aluminium Pull To Lock	N/C
<input type="checkbox"/> GLYSTPH	Steel Push To Lock	N/C
<input type="checkbox"/> GLYSTPL	Steel Pull To Lock	N/C
<input type="checkbox"/> GLYBRXC	4" Brake Extension	N/C
<input type="checkbox"/> GLYBRXA	6" Brake Extension	N/C
<input type="checkbox"/> GLYBRXB	9" Brake Extension	130

SEATING - BACKS *For custom, choose one back then fill custom* **12**

W () x H ()	
<input type="checkbox"/> GLYBK01	RELAX Gentle Fit Back 544
<input type="checkbox"/> GLYBK02	RELAX Full Fit Back 608
<input type="checkbox"/> GLYBK03	RELAX Snug Fit Back 700
<input type="checkbox"/> GLYBK04	RELAX Form Fit Back No Laterals 820
<input type="checkbox"/> GLYBK06	RELAX Form Fit Back Gentle Laterals 820
<input type="checkbox"/> GLYBK07	RELAX Form Fit Back Full Laterals 820
<input type="checkbox"/> GLYBK05	RELAX Form Fit Back Snug Laterals 820
<input type="checkbox"/> GLYBK10	RELAX Spiral Fit Back Gentle Laterals 832
<input type="checkbox"/> GLYBK09	RELAX Spiral Fit Back Snug Laterals 832
<input type="checkbox"/> GLYBK11	RELAX Spiral Fit Back Full Laterals 832
<input type="checkbox"/> GLYTAB	RELAX Tension Back Upto 18" 140
<input type="checkbox"/> GLYTAB+	RELAX Tension Back Over 18" 240
<input type="checkbox"/> GLYBGL1	Half Overlay Gel (3 Piece) 348
<input type="checkbox"/> GLYBGL2	Full Overlay Gel (6 Pieces) 696
<input type="checkbox"/> GLYPGLY	Pin Extension (1" Extension) 40
<input type="checkbox"/> GLYBK12	*Custom Width : _____ +176
<input type="checkbox"/> GLYBK13	*Custom Height : _____ +176

*Note: Please refer to the Seating Booklet for custom sizes

SEATING - CUSHION *Select one below* **13**

W () x D ()	
<input type="checkbox"/> GLYC01	RELAX Soft Cushion 480
<input type="checkbox"/> GLYC02	RELAX Comfy Cushion 624
<input type="checkbox"/> GLYC03	RELAX Firm Cushion 784
<input type="checkbox"/> GLYCGL1	Half Overlay Gel (1 Piece) 116
<input type="checkbox"/> GLYCGL2	Full Overlay Gel (2 Pieces) 232
<input type="checkbox"/> GLYCCL	Contoured Leg Channel (Bolster) Std. 180
<input type="checkbox"/> GLYCCLC	Contoured Leg Channel (Bolster) Custom 320

*Note: Please refer to the Seating Booklet for custom sizes

<input type="checkbox"/> GLYC00	Extra Cushion Cover	99
---------------------------------	---------------------	----

ACCESSORIES **14**

<input type="checkbox"/> GLYSBA	Seat Belt (Auto Style)	26
<input type="checkbox"/> GLYSBAP	Seat Belt (Padded)	56
<input type="checkbox"/> GLYSBA4P	Seat Belt (4 Point Padded)	130
<input type="checkbox"/> GLYTRA	Padded Tray Locking	280
<input type="checkbox"/> GLYTRB	Padded Tray Slide On	280
<input type="checkbox"/> GLYTRC	Full Plastic Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/> GLYTRD	Half Padded Left Tray	240
<input type="checkbox"/> GLYTRE	Half Padded Right Tray	240
<input type="checkbox"/> GLYTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/> GLYTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/> GLYTDFB	Transport Tie Downs Front	85
<input type="checkbox"/> GLYTDRB	Transport Tie Downs Rear	85
<input type="checkbox"/> GLYOXY	Oxygen Tank Holder	195
<input type="checkbox"/> GLYIVA	Adjustable I V Pole	130
<input type="checkbox"/> GLYCHPT-L	Left Arm Trough Tray	400
<input type="checkbox"/> GLYCHPT-R	Right Arm Trough Tray	400
<input type="checkbox"/> GLYADD8	8" Adductor	164
<input type="checkbox"/> GLYADD10	10" Adductor	164