

MGP Magic Plus Wheelchair - Folding Wheelchair



INFORMATION

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR Select one below **1 ADP**

- MGPS Frame - Silver \$1,690 WM0001062
- MGPB Frame - Black \$1,690 WM0001062
- MGPC Frame - Colour : _____ \$2,090 WM0001062

CLIENT WEIGHT Select one below **2 ADP**

- MGPWA Up to 250Lbs N/C
- MGPWAMF Exceeds 250Lbs \$730 WAMF
- MGPWAMG Exceeds 350Lbs \$992 WAMG
- MGPWAMH Exceeds 400Lbs (Max 600Lbs) \$1,110 WAMH

SEAT WIDTH Select one below **3 ADP**

- MGPSW14 Seat Width 14" N/C
- MGPSW15 Seat Width 15" N/C
- MGPSW16 Seat Width 16" N/C
- MGPSW17 Seat Width 17" N/C
- MGPSW18 Seat Width 18" N/C
- MGPSW19 Seat Width 19" \$399 WAMA
- MGPSW20 Seat Width 20" \$399 WAMA
- MGPSWS Custom Seat Width: _____ \$399 WAMA

SEAT DEPTH Select one below **4 ADP**

- MGPSD14 Seat Depth 14" N/C
- MGPSD15 Seat Depth 15" N/C
- MGPSD16 Seat Depth 16" N/C
- MGPSD17 Seat Depth 17" N/C
- MGPSD18 Seat Depth 18" N/C
- MGPSD19 Seat Depth 19" \$265 WAMB
- MGPSD20 Seat Depth 20" \$265 WAMB
- MGPSDS Custom Seat Depth: _____ \$265 WAMB

BACK HEIGHT Select one below **5 ADP**

- MGPBH15 Back Height 15" N/C
- MGPBH16 Back Height 16" N/C
- MGPBH17 Back Height 17" N/C
- MGPBH18 Back Height 18" N/C
- MGPBH19 Back Height 19" N/C
- MGPBH20 Back Height 20" N/C
- MGPBH21 Back Height 21" N/C
- MGPBHC Back Height Custom : _____ \$120

Notes:

STANDARD NO CHARGE SPECIFICATIONS

- LIGHT WEIGHT ALUMINIUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINIUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- VERTICAL REAR AXLE ADJUSTIBILITY
- ANTI-TIP ROLLERS
- 6" BRAKE EXTENSIONS
- DOUBLE CROSS BRACE

SEAT TO FLOOR SPECIFICATIONS (± 0.5")

(LOW SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	12.75"	20"	5"	2H-N
<input type="checkbox"/>	13"	20"	5"	3H,2H,1H
<input type="checkbox"/>	13.75"	20"	5"	3H,2H
<input type="checkbox"/>	14.25"	20"	6"	5H,3H,2H
<input type="checkbox"/>	14.5"	22"	5"	5H,3H
<input type="checkbox"/>	14.25"	22"	6"	5H,3H,2H
<input type="checkbox"/>	15.25"	20"	5"	5H
<input type="checkbox"/>	15"	20"	6"	5H,3H
<input type="checkbox"/>	15"	22"	5"	5H
<input type="checkbox"/>	15"	22"	6"	5H,3H
<input type="checkbox"/>	16.25"	22"	6"	5H
<input type="checkbox"/>	16.5"	22"	8"	5H
<input type="checkbox"/>	17.25"	22"	8"	5H

(HEIGHT EXTENDED CASTER)				
	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	18.25"	24"	8"	5H,3H

For more Seat To Floor combinations, visit our website

FOOT REST AND FOOT PLATES		Select one below	5	ADP
<input type="checkbox"/>	MGPF6A	60 Degree Hanger Pivot	N/C	
<input type="checkbox"/>	MGPF7A	70 Degree Hanger Pivot	N/C	
<input type="checkbox"/>	MGPF8A	80 Degree Hanger Pivot	N/C	
<input type="checkbox"/>	MGPFA	Elevating Legrest Pivot	\$404	WAMK
<input type="checkbox"/>	MGPAL	Articulating Leg Rests Pivot	\$450	WAMK
<input type="checkbox"/>	MGPPPL	Plastic Footplates	N/C	
<input type="checkbox"/>	MGPPXL	Extra Large Footplates	\$150	SE0000322
<input type="checkbox"/>	MGPAPF	Angle Adjustable Footplates	\$240	WAMJ
<input type="checkbox"/>	MGP1FB	One Piece Foot Board	\$80	SEFND0015
<input type="checkbox"/>	MGP1FBP	One Piece Foot Board Padding	\$44	SEFND0035
<input type="checkbox"/>	MGPNLS	Nylon Leg Strap Padded	\$60	SEFND0010
<input type="checkbox"/>	MGPNLS	Nylon Leg Strap Custom Size ()	\$120	SEFCF107L
<input type="checkbox"/>	MGPNSL	Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style	\$240	SEFND0020

BACK TYPE		Select one below	6	ADP
<input type="checkbox"/>	MGPAB1	Angle Adjustable Back	N/C	

ARM REST		Select one below	8	ADP
(Ht. 10"-14")	<input type="checkbox"/>	MGPFA	Flip Back Ht. Adj./ Removable	\$340
(Ht. 10"-14")	<input type="checkbox"/>	MGPFA	T-Arm Ht. Adj./ Removable	\$240
(Ht. 7"-11")	<input type="checkbox"/>	MGPFA	Flip Back Ht. Adj./ Removable Short	N/C
(Ht. 7"-11")	<input type="checkbox"/>	MGPFA	T-Arm Ht. Adj./ Removable Short	\$340
<input type="checkbox"/>	MGPWAPS	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030

* Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required

REAR WHEEL		Select one below	9	ADP
<input type="checkbox"/>	MGRWTO	20" Solid Treaded Tire	N/C	
<input type="checkbox"/>	MGRWT2	22" Solid Treaded Tire	N/C	
<input type="checkbox"/>	MGRWT4	24" Solid Treaded Tire	N/C	
<input type="checkbox"/>	MGRW50	20" Spoked Treaded Tire	N/C	
<input type="checkbox"/>	MGRW52	22" Spoked Treaded Tire	N/C	
<input type="checkbox"/>	MGRW54	24" Spokedtreaded Tire	N/C	
<input type="checkbox"/>	MGRWA2	22" Air Tire	\$120	
<input type="checkbox"/>	MGRWA4	24" Air Tire	\$120	
<input type="checkbox"/>	MGPAQR	Axles – Quick Release	\$120	WAMU

CASTER		Select one below	10	ADP
<input type="checkbox"/>	MGPFC1A	5" X 1.25" Caster	N/C	
<input type="checkbox"/>	MGPFC1B	6" X 1.25" Caster	N/C	
<input type="checkbox"/>	MGPFC1C	8" X 1.25" Caster	N/C	
<input type="checkbox"/>	MGPFC2A	5" X 1.75" Caster (Soft Roll)	\$80	
<input type="checkbox"/>	MGPFC2B	6" X 1.75" Caster	\$80	
<input type="checkbox"/>	MGPFC2C	8" X 1.75" Caster	\$80	
<input type="checkbox"/>	MGPFC4B	6" X 2 Caster	\$80	
<input type="checkbox"/>	MGPFC4C	8" X 2 Caster	\$80	

WHEEL LOCK		Select one below	11	ADP
<input type="checkbox"/>	MGPALPH	Aluminium Push To Lock	N/C	
<input type="checkbox"/>	MGPALPL	Aluminium Pull To Lock	N/C	
<input type="checkbox"/>	MGPSTPH	Steel Push To Lock	N/C	
<input type="checkbox"/>	MGPSTPL	Steel Pull To Lock	N/C	
<input type="checkbox"/>	MGPBRXC	4" Brake Extension	N/C	
<input type="checkbox"/>	MGPBRXA	6" Brake Extension	N/C	
<input type="checkbox"/>	MGPBRXB	9" Brake Extension	\$130	

SEATING -BACKS		Custom: select back & fill custom	12	ADP
W () x H ()				
<input type="checkbox"/>	MGPBK01	RELAX Gentle Fit Back	\$544	SEBWM0020
<input type="checkbox"/>	MGPBK02	RELAX Full Fit Back	\$608	SEBWM0010
<input type="checkbox"/>	MGPBK03	RELAX Snug Fit Back	\$700	SEBWM0027
<input type="checkbox"/>	MGPBK04	RELAX Form Fit Back No Laterals	\$820	SEBWM0026
<input type="checkbox"/>	MGPBK06	RELAX Form Fit Back Gentle Laterals	\$820	SEBWM0026
<input type="checkbox"/>	MGPBK07	RELAX Form Fit Back Full Laterals	\$820	SEBWM0026
<input type="checkbox"/>	MGPBK05	RELAX Form Fit Back Snug Laterals	\$820	SEBWM0026
<input type="checkbox"/>	MGPBK10	RELAX Spiral Fit Back Gentle Laterals	\$832	SEBWM0030
<input type="checkbox"/>	MGPBK09	RELAX Spiral Fit Back Snug Laterals	\$832	SEBWM0030
<input type="checkbox"/>	MGPBK11	RELAX Spiral Fit Back Full Laterals	\$832	SEBWM0030
<input type="checkbox"/>	MGPBTA	RELAX Tension Back Upto 18"	\$140	WAMD
<input type="checkbox"/>	MGPBTA+	RELAX Tension Back Over 18"	\$240	WAME
<input type="checkbox"/>	MGPMPGP	Pin Extension (1" Extension)	\$40	
<input type="checkbox"/>	MGPBK12	* Custom Width : _____	\$176	SEBWM0035
<input type="checkbox"/>	MGPBK13	* Custom Height : _____	\$176	SEBWM0040

*Note: Please refer to the Seating Booklet for custom sizes

<input type="checkbox"/>	MGPBGL1	Gel Full Overlay	\$348	SEBND0020 X 3
<input type="checkbox"/>	MGPBGL2	Gel Center Vertebrae	\$116	SEBND0020
<input type="checkbox"/>	MGPBGL3	Gel Left Lateral	\$116	SEBND0020
<input type="checkbox"/>	MGPBGL4	Gel Right Lateral	\$116	SEBND0020

Maximum 1 piece additional pressure relief gel

SEATING - CUSHIONS		Select one below	13	ADP
W () x D ()				
<input type="checkbox"/>	MGPC01	RELAX Soft Cushion	\$480	SESWM0050
<input type="checkbox"/>	MGPC02	RELAX Comfy Cushion	\$624	SESWM0030
<input type="checkbox"/>	MGPC03	RELAX Firm Cushion	\$784	SESWM0040
<input type="checkbox"/>	MGPC04	RELAX AIR FLOW Cushion	\$587	SE0001669
<input type="checkbox"/>	MGPCGL2	Full Overlay Gel	\$258	SESND1045
<input type="checkbox"/>	MGPCCL	Contoured Leg Channel (Bolster) Std.	\$180	SEMCF203L
<input type="checkbox"/>	MGPCCLA	Contoured Leg Channel (Bolster) with ABS	\$260	SESND1035 SESND1015
<input type="checkbox"/>	MGPCCLC	* Contoured Leg Channel (Bolster) Custom	\$320	SEMCF209L

*Note: Please refer to the Seating Booklet for custom sizes

<input type="checkbox"/>	MGPC00	Extra Cushion Cover	\$99	SESND1055
--------------------------	--------	---------------------	------	-----------

ACCESSORIES			14	ADP
<input type="checkbox"/>	MGPSBA	Seat Belt (Auto Style)	\$26	
<input type="checkbox"/>	MGPSBAP	Seat Belt (Padded)	\$56	SERND0030
<input type="checkbox"/>	MGPSBA4P	Seat Belt (4 Point Padded)	\$130	SERND0020
<input type="checkbox"/>	MGPTRA	Full Padded Tray Locking	\$280	SETND0005 SETND0035
<input type="checkbox"/>	MGPTRB	Full Padded Tray Slide On	\$280	SETND0005 SETND0035
<input type="checkbox"/>	MGPTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	MGPTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	MGPTRD	Half Padded Left Tray	\$240	SETND0015 SETND0035
<input type="checkbox"/>	MGPTRR	Half Padded Right Tray	\$240	SETND0015 SETND0035
<input type="checkbox"/>	MGPTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>	MGPTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>	MGPTRFB	Transport Tie Downs Front	\$85	
<input type="checkbox"/>	MGPTRRB	Transport Tie Downs Rear	\$85	
<input type="checkbox"/>	MGPPOXY	Oxygen Tank Holder	\$195	WEPX
<input type="checkbox"/>	MGPPIVA	Adjustable I V Pole	\$130	
<input type="checkbox"/>	MGPCHPT-L	Left Arm Trough Tray	\$400	SETND0015 SETND0035
<input type="checkbox"/>	MGPCHPT-R	Right Arm Trough Tray	\$400	SETND0030 SEMCF305L
<input type="checkbox"/>	MGPPLATB	Lateral Bite out	\$180	SELND0015
<input type="checkbox"/>	MGPADD8	8" Adductor	\$164	SELND0005
<input type="checkbox"/>	MGPADD10	10" Adductor	\$164	SELND0005
<input type="checkbox"/>	MGPADD12	12" Adductor	\$164	SELND0005
<input type="checkbox"/>	MGPFRCH	Padded Footplates	\$88	SEMCF803L