

MGP Magic Plus Wheelchair (Manual Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

- MGPS Frame - Silver 1590
- MGPB Frame - Black 1590
- MGPC Frame - Colour : _____ 1990

CLIENT WEIGHT *Select one below* **2**

- MGPWA Upto 250Lbs N/C
- MGPWAMF Exceeds 250Lbs 730
- MGPWAMG Exceeds 350Lbs 992
- MGPWAMH Exceeds 400Lbs (Max 600Lbs) 1110

SEAT WIDTH *Select one below* **3**

- MGPSW14 Seat Width 14" N/C
- MGPSW15 Seat Width 15" N/C
- MGPSW16 Seat Width 16" N/C
- MGPSW17 Seat Width 17" N/C
- MGPSW18 Seat Width 18" N/C
- MGPSW19 Seat Width 19" 399
- MGPSW20 Seat Width 20" 399
- MGPSWS Custom Seat Width: _____ 399

SEAT DEPTH *Select one below* **4**

- MGPSD14 Seat Depth 14" 350
- MGPSD15 Seat Depth 15" 350
- MGPSD16 Seat Depth 16" N/C
- MGPSD17 Seat Depth 17" N/C
- MGPSD18 Seat Depth 18" N/C
- MGPSD19 Seat Depth 19" 265
- MGPSD20 Seat Depth 20" 265
- MGPSDS Custom Seat Depth: _____ 265

BACK HEIGHT *Select one below* **5**

- MGPBH15 Back Height 15" N/C
- MGPBH16 Back Height 16" N/C
- MGPBH17 Back Height 17" N/C
- MGPBH18 Back Height 18" N/C
- MGPBH19 Back Height 19" N/C
- MGPBH20 Back Height 20" N/C
- MGPBH21 Back Height 21" N/C
- MGPBHC Back Height Custom : _____ 120

Notes:

SPECIFICATIONS

- LIGHT WEIGHT ALUMINIUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- VERTICAL REAR AXLE ADJUSTIBILITY
- ANTI-TIP ROLLERS
- 6" BRAKE EXTENSIONS
- DOUBLE CROSS BRACE

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	SL 12.75"	20"	5"	XL SHORT
<input type="checkbox"/>	SL 13.00"	20"	5"	XL SHORT
<input type="checkbox"/>	LL 13.75"	20"	5"	XL SHORT
<input type="checkbox"/>	L 14.50"	20"	6"	SHORT
<input type="checkbox"/>	L 14.50"	22"	5"	SHORT
<input type="checkbox"/>	L14.50"	22"	6"	SHORT
<input type="checkbox"/>	LH 15.25"	20"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	20"	6"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	5"	SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	6"	SHT/LONG
<input type="checkbox"/>	SH 16.00"	22"	6"	LONG

(STANDARD SEAT SETTING)

<input type="checkbox"/>	MH 16.75"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	22"	8"	LONG
<input type="checkbox"/>	H 17.50"	22"	6"	LONG
<input type="checkbox"/>	H 17.50"	24"	8"	LONG
<input type="checkbox"/>	A 18.25"	24"	8"	LONG

BACK TYPE		6
<input type="checkbox"/>	MGPAB1 Angle Adjustable Back	N/C
FOOT REST AND FOOT PLATES <i>Select one below</i>		7
<input type="checkbox"/>	MGPF6A 60 Degree Hanger	N/C
<input type="checkbox"/>	MGPF7A 70 Degree Hanger	N/C
<input type="checkbox"/>	MGPF8A 80 Degree Hanger	N/C
<input type="checkbox"/>	MGPPPL Plastic	N/C
<input type="checkbox"/>	MGPPXL Extra Large	150
<input type="checkbox"/>	MGPPAFE Angle Adjustable	240
<input type="checkbox"/>	MGPEA Elevating Legrest	404
<input type="checkbox"/>	MGPAL Articulating Leg Rests	450
<input type="checkbox"/>	MGP1FB One Piece Foot Board	80
<input type="checkbox"/>	MGPNLS Nylon Leg Strap Padded	60
<input type="checkbox"/>	MGPCSL Contracture Sling	240
ARM REST <i>Select one below</i>		8
(Std 10"-14") <input type="checkbox"/>	MGPFA Flip Back Ht. Adj./ Removable	240
(Std 10"-14") <input type="checkbox"/>	MGPTA T-Arm Ht. Adj./ Removable	340
(Short 7"-11") <input type="checkbox"/>	MGPFAS Flip Back Ht. Adj./ Removable	240
(Short 7"-11") <input type="checkbox"/>	MGPTAS T-Arm Ht. Adj./ Removable	340
<input type="checkbox"/>	MGPGWAPF Gel Waterfall Arm Pads	180
REAR WHEELS <i>Select one below</i>		9
<input type="checkbox"/>	MGPRWT0 20" Solid Treaded Tire	N/C
<input type="checkbox"/>	MGPRWT2 22" Solid Treaded Tire	N/C
<input type="checkbox"/>	MGPRWT4 24" Solid Treaded Tire	N/C
<input type="checkbox"/>	MGPpRWS0 20" Spoked Treaded Tire	N/C
<input type="checkbox"/>	MGPpRWS2 22" Spoked Treaded Tire	N/C
<input type="checkbox"/>	MGPpRWS4 24" Spokedtreaded Tire	N/C
<input type="checkbox"/>	MGPRWA2 22" Air Tire	120
<input type="checkbox"/>	MGPRWA4 24" Air Tire	120
<input type="checkbox"/>	MGPAQR Axles – Quick Release	120
CASTER <i>Select one below</i>		10
<input type="checkbox"/>	MGPFC1A 5" X 1.25" Caster	N/C
<input type="checkbox"/>	MGPFC1B 6" X 1.25" Caster	N/C
<input type="checkbox"/>	MGPFC1C 8" X 1.25" Caster	N/C
<input type="checkbox"/>	MGPFC2A 5" X 1.75" Caster (Soft Roll)	80
<input type="checkbox"/>	MGPFC2B 6" X 1.75" Caster	80
<input type="checkbox"/>	MGPFC2C 8" X 1.75" Caster	80
<input type="checkbox"/>	MGPFC4B 6" X 2 Caster	80
<input type="checkbox"/>	MGPFC4C 8" X 2 Caster	80
WHEEL LOCK <i>Select one below</i>		11
<input type="checkbox"/>	MGPALPH Aluminium Push To Lock	N/C
<input type="checkbox"/>	MGPALPL Aluminium Pull To Lock	N/C
<input type="checkbox"/>	MGPSTPH Steel Push To Lock	N/C
<input type="checkbox"/>	MGPSTPL Steel Pull To Lock	N/C
<input type="checkbox"/>	MGPBRXC 4" Brake Extension	N/C
<input type="checkbox"/>	MGPBRXA 6" Brake Extension	N/C
<input type="checkbox"/>	MGPBRXB 9" Brake Extension	130

SEATING - BACKS <i>For custom, choose one back then fill custom</i>		12
W () x H ()		
<input type="checkbox"/>	MGPBK01 RELAX Gentle Fit Back	544
<input type="checkbox"/>	MGPBK02 RELAX Full Fit Back	608
<input type="checkbox"/>	MGPBK03 RELAX Snug Fit Back	700
<input type="checkbox"/>	MGPBK04 RELAX Form Fit Back No Laterals	820
<input type="checkbox"/>	MGPBK06 RELAX Form Fit Back Gentle Laterals	820
<input type="checkbox"/>	MGPBK07 RELAX Form Fit Back Full Laterals	820
<input type="checkbox"/>	MGPBK05 RELAX Form Fit Back Snug Laterals	820
<input type="checkbox"/>	MGPBK10 RELAX Spiral Fit Back Gentle Laterals	832
<input type="checkbox"/>	MGPBK09 RELAX Spiral Fit Back Snug Laterals	832
<input type="checkbox"/>	MGPBK11 RELAX Spiral Fit Back Full Laterals	832
<input type="checkbox"/>	MGPTAB RELAX Tension Back Upto 18"	140
<input type="checkbox"/>	MGPTAB+ RELAX Tension Back Over 18"	240
<input type="checkbox"/>	MGPBGL1 Half Overlay Gel (3 Piece)	348
<input type="checkbox"/>	MGPBGL2 Full Overlay Gel (6 Pieces)	696
<input type="checkbox"/>	MGPMMGP Pin Extension (1" Extension)	40
	MGPBK12 *Custom Width : _____	+176
	MGPBK13 *Custom Height : _____	+176

*Note: Please refer to the Seating Booklet for custom sizes

SEATING - CUSHION <i>Select one below</i>		13
W () x D ()		
<input type="checkbox"/>	MGPC01 RELAX Soft Cushion	480
<input type="checkbox"/>	MGPC02 RELAX Comfy Cushion	624
<input type="checkbox"/>	MGPC03 RELAX Firm Cushion	784
<input type="checkbox"/>	MGPCGL1 Half Overlay Gel (1 Piece)	116
<input type="checkbox"/>	MGPCGL2 Full Overlay Gel (2 Pieces)	232
<input type="checkbox"/>	MGPCCL Contoured Leg Channel (Bolster) Std.	180
<input type="checkbox"/>	MGPCCL Contoured Leg Channel (Bolster) with ABS	260
<input type="checkbox"/>	MGPCCLC Contoured Leg Channel (Bolster) Custom	320
*Note: Please refer to the Seating Booklet for custom sizes		
<input type="checkbox"/>	MGPC00 Extra Cushion Cover	99

ACCESSORIES		14
<input type="checkbox"/>	MGPSBA Seat Belt (Auto Style)	26
<input type="checkbox"/>	MGPSBAP Seat Belt (Padded)	56
<input type="checkbox"/>	MGPSBA4P Seat Belt (4 Point Padded)	130
<input type="checkbox"/>	MGPTRA Padded Tray Locking	280
<input type="checkbox"/>	MGPTRB Padded Tray Slide On	280
<input type="checkbox"/>	MGPTRC Full Plastic Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/>	MGPTRD Half Padded Left Tray	240
<input type="checkbox"/>	MGPTRE Half Padded Right Tray	240
<input type="checkbox"/>	MGPTRF Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	MGPTRG Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	MGPTDFB Transport Tie Downs Front	85
<input type="checkbox"/>	MGPTDRB Transport Tie Downs Rear	85
<input type="checkbox"/>	MGPOXY Oxygen Tank Holder	195
<input type="checkbox"/>	MGPIVA Adjustable I V Pole	130
<input type="checkbox"/>	MGPCHPT-L Left Arm Trough Tray	400
<input type="checkbox"/>	MGPCHPT-R Right Arm Trough Tray	400
<input type="checkbox"/>	MGPADD8 8" Adductor	164
<input type="checkbox"/>	MGPADD10 10" Adductor	164
<input type="checkbox"/>	MGPADD12 12" Adductor	164