

STP Super Tilt Plus Wheelchair - Tilt Wheelchair



INFORMATION

Customer Name: _____
Address: _____ PO #: _____
City: _____ Phone #: _____ Date: _____
Account #: _____ Authorizer Name: _____



FRAME COLOUR		Select one below	1
☐ MP5B	Tilt Frame - Black		\$3,239
☐ MP5S	Tilt Frame - Silver		\$3,239
☐ MP5C	Tilt Frame - Custom : _____		\$3,719

CLIENT WEIGHT		Select one below	2
☐ STPWA	Up to 250Lbs		N/C
☐ STPWAMF	Exceeds 250Lbs		\$876
☐ STPWAMG	Exceeds 350Lbs		\$1,190
☐ STPWAMH	Exceeds 400Lbs (Max 600Lbs)		\$1,332

SEAT WIDTH		Select one below	3
☐ STPSW16	Seat Width 16"		N/C
☐ STPSW17	Seat Width 17"		N/C
☐ STPSW18	Seat Width 18"		N/C
☐ STPSW19	Seat Width 19"		\$479
☐ STPSW20	Seat Width 20"		\$479
☐ STPSWC	Custom Seat Width: _____		\$479

SEAT DEPTH		Select one below	4
☐ STPSD16	Seat Depth 16"		N/C
☐ STPSD17	Seat Depth 17"		N/C
☐ STPSD18	Seat Depth 18"		N/C
☐ STPSD19	Seat Depth 19"		\$318
☐ STPSD20	Seat Depth 20"		\$318
☐ STPSDC	Custom Seat Depth: _____		\$318

FOOT REST AND FOOT PLATES		Select one below	5
☐ STPOFST	Offset Foot Rest		N/C
☐ STPF6A	60 Degree Hanger ☐ Pivot ☐ Pin		N/C
☐ STPF7A	70 Degree Hanger ☐ Pivot ☐ Pin		N/C
☐ STPF8A	80 Degree Hanger ☐ Pivot ☐ Pin		N/C
☐ STPEA	Elevating Legrest ☐ Pivot ☐ Pin		\$485
☐ STPAL	Articulating L/R ☐ Pivot ☐ Pin		\$540
☐ STPEB	Elevating Legrest System		\$540
☐ STPPPL	Plastic Footplates		N/C
☐ STPPXL	Extra Large Footplates		\$180
☐ STPPAF	Angle Adjustable Footplates		\$288
☐ STP1FB	One Piece Foot Board		\$96
☐ STP1FBP	One Piece Foot Board Padding		\$53
☐ STPNLS	Nylon Leg Strap Padded		\$72
☐ STPNLSC	Nylon Leg Strap Custom Size ()		\$144
☐ STPCSL	Contracture Sling ☐ Standard ☐ M Style		\$288

BACK TYPE		Select one below	7
☐ STPAB1	Angle Adjustable Back		N/C
☐ STPENDO	Endomorph Canes W ()		N/C
☐ STPRBK	Recline Kit		\$672

STANDARD NO CHARGE SPECIFICATIONS	
SEAT WIDTH: 14" – 18"	
SEAT DEPTH (ADJUSTABLE)	
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"	
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]	
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]	
WEIGHT CAPACITY: 250 LBS	
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)	
TILT ADJUSTMENTS 0 TO 50 (DEGREES)	
TWO-WAY LOCKING GAS STRUT	
SINGLE RELEASE, DUAL CYLINDERS	
STANDARD LOCKOUT COLLARS	
STANDARD T STYLE ARMS	
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS	

TILT OPTION		Select one below	6
☐ STPINTS	Independant Power Tilt System		\$2,400

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
☐	13.00"	20"	5"	XL SHORT
☐	13.75"	20"	5"	XL SHORT
☐	14.50"	20"	6"	SHORT
☐	14.50"	22"	5"	SHORT
☐	14.50"	22"	6"	SHORT
☐	15.25"	20"	5"	SHT/LONG
☐	15.25"	20"	6"	SHT/LONG
☐	15.25"	22"	5"	SHT/LONG
☐	15.25"	22"	6"	SHT/LONG
☐	16.00"	22"	6"	LONG
(HIGH SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
☐	16.75"	22"	6"	LONG
☐	17.50"	22"	6"	LONG
☐	17.50"	22"	6"	LONG
☐	17.50"	24"	6"	LONG
☐	18.25"	24"	8"	LONG

For more Seat To Floor combinations, visit our website

HEADREST		Select one below	8
☐	STPMAHF	Multi-Axis Hardware Complete with Large Pad	\$614
☐	STPMAHFXL	Multi-Axis Hardware Complete with Flex Large	\$614
☐	STPMAHKP	Multi-Axis Hardware Complete with Medium Pad	\$614
☐	STPMAHFXS	Multi-Axis Hardware Complete with Flex Medium	\$614
☐	STPMAHFT	Multi-Axis Hardware Complete with Small Flat	\$614
☐	STPUPHT219	Contour Headrest Complete with Swivel Hardware	\$520
☐	STPUPHT218	Contour Headrest with Std. Assembly (1/2" Tube)	\$350
☐	STPUPHT217	Contour Headrest with Std. Assembly (3/4" Tube)	\$265
☐	STPUHRA	Universal Hardware Adaptor (Adaptor only)	\$68

ARM REST		Select one below	9
(Ht. 10"-14")	☐ STPTA	T-Arm Ht. Adj./ Removable	\$288
(Ht. 10"-14")	☐ STPTAS	T-Arm Ht. Adj./ Removable Short	\$288
	☐ STPGWAPF	* Gel Waterfall Arm Pads	\$216

*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required

REAR WHEEL		Select one below	10
☐	STPRWTO	20" Solid Treaded Tire	N/C
☐	STPRWT2	22" Solid Treaded Tire	N/C
☐	STPRWT4	24" Solid Treaded Tire	N/C
(STF 13-19")	☐ STPRWT12	* 12" Solid Tire	N/C
lower frame	☐ STPRWT12I	* 12" Solid Tire Inside Mount	\$144
☐	STPRWT16	* 16" Solid Tire	N/C
☐	STPRWA2	22" Air Tire	\$144
☐	STPRWA4	24" Air Tire	\$144
☐	STPAQR	Axles – Quick Release	\$144
☐	STPHRT	Narrow Handrim Config. (1/2" less)	N/C

*Note: Quick Release is not available for wheels smaller than 20"

CASTER		Select one below	11
☐	STPFC1A	5" X 1.25" Caster	N/C
☐	STPFC1B	6" X 1.25" Caster	N/C
☐	STPFC1C	8" X 1.25" Caster	N/C
☐	STPFC2A	5" X 1.75" Caster	\$96
☐	STPFC2B	6" X 1.75" Caster	\$96
☐	STPFC2C	8" X 1.75" Caster	\$96
☐	STPFC4B	6" X 2 Caster	\$96
☐	STPFC4C	8" X 2 Caster	\$96

WHEEL LOCK		Select one below	12
☐	STPALPH	Aluminium Push To Lock	N/C
☐	STPALPL	Aluminium Pull To Lock	N/C
☐	STPSTPH	Steel Push To Lock	N/C
☐	STPSTPL	Steel Pull To Lock	N/C
☐	STPBRXC	4" Brake Extension	N/C
☐	STPBRXA	6" Brake Extension	N/C
☐	STPBRXB	9" Brake Extension	\$156
☐	STPATTBRX	Attendant Brakes	\$288
☐	STPULB	Unilateral Hand Brake ☐ Left ☐ Right	\$144

SEATING - BACKS		Custom: select back & fill custom	13
	W () x H ()		
☐	STPBK01	RELAX Gentle Fit Back	\$653
☐	STPBK02	RELAX Full Fit Back	\$730
☐	STPBK03	RELAX Snug Fit Back	\$840
☐	STPBK04	RELAX Form Fit Back No Laterals	\$984
☐	STPBK06	RELAX Form Fit Back Gentle Laterals	\$984
☐	STPBK07	RELAX Form Fit Back Full Laterals	\$984
☐	STPBK05	RELAX Form Fit Back Snug Laterals	\$984
☐	STPBK10	RELAX Spiral Fit Back Gentle Laterals	\$998
☐	STPBK09	RELAX Spiral Fit Back Snug Laterals	\$998
☐	STPBK11	RELAX Spiral Fit Back Full Laterals	\$998
☐	STPTAB	RELAX Tension Back Upto 18"	\$168
☐	STPTAB+	RELAX Tension Back Over 18"	\$288
☐	STPBK12	* Custom Width : _____	\$211
☐	STPBK13	* Custom Height : _____	\$211

*Note: Please refer to the Seating Booklet for custom sizes

☐	STPBGL1	Gel Full Overlay	\$418
☐	STPBGL2	Gel Center Vertebrae	\$139
☐	STPBGL3	Gel Left Lateral	\$139
☐	STPBGL4	Gel Right Lateral	\$139

Maximum 1 piece additional pressure relief gel

SEATING - CUSHIONS		Select one below	14
	W () x D ()		
☐	STPC01	RELAX Soft Cushion	\$576
☐	STPC02	RELAX Comfy Cushion	\$749
☐	STPC03	RELAX Firm Cushion	\$941
☐	STPC04	RELAX AIR FLOW Cushion	\$704
☐	STPCGL2	Full Overlay Gel	\$310
☐	STPCCL	Contoured Leg Channel (Bolster) Std.	\$216
☐	STPCCLA	Contoured Leg Channel (Bolster) with ABS	\$312
☐	STPCCLC	* Contoured Leg Channel (Bolster) Custom	\$384

*Note: Please refer to the Seating Booklet for custom sizes

☐	STPC00	Extra Cushion Cover	\$119
--------------------------	--------	---------------------	-------

ACCESSORIES		15	
☐	STPSBA	Seat Belt (Auto Style)	\$31
☐	STPSBAP	Seat Belt (Padded)	\$67
☐	STPSBA4P	Seat Belt (4 Point Padded)	\$156
☐	STPTRA	Full Padded Tray Locking	\$336
☐	STPTRB	Full Padded Tray Slide On	\$336
☐	STPTRCL	Full Plastic Tray Locking ☐ ABS Black ☐ Clear	\$336
☐	STPTRC	Full Plastic Tray Slide On ☐ ABS Black ☐ Clear	\$336
☐	STPTRD	Half Padded Left Tray	\$288
☐	STPTRE	Half Padded Right Tray	\$288
☐	STPTRF	Half Left Tray ☐ ABS Black ☐ Clear	\$288
☐	STPTRG	Half Right Tray ☐ ABS Black ☐ Clear	\$288
☐	STPTDFB	Transport Tie Downs Front	\$102
☐	STPTDRB	Transport Tie Downs Rear	\$102
☐	STPOXY	Oxygen Tank Holder	\$234
☐	STPIVA	Adjustable I V Pole	\$156
☐	STPCHPT-L	Left Arm Trough Tray	\$480
☐	STPCHPT-R	Right Arm Trough Tray	\$480
☐	STPLATB	Lateral Bite out	\$216
☐	STPADD8	8" Adductor	\$197
☐	STPADD10	10" Adductor	\$197
☐	STPADD12	12" Adductor	\$197
☐	STPFRCH	Padded Footplates	\$106