

NVT NV Dynamic Tilt Wheelchair



INFORMATION

Customer Name: _____
Address: _____ PO #: _____
City: _____ Phone #: _____ Date: _____
Account #: _____ Authorizer Name: _____



DYNAMIC FRAME		1
☒ NVTB	Dynamic Tilt Frame - Black	\$3,239
☒ NVTWANA	Standard Manual Wheelchair Frame with Dynamic T	\$1,134
☒ NVTWAMV	Clothing Guard	\$119

*Note: All three ADP codes apply when you order the NV Tilt Wheelchair

FRAME COLOUR		2
☐ NVTS	Tilt Frame - Silver	\$3,239
☐ NVTBL	Tilt Frame - Blue	\$3,239
☐ NVTBU	Tilt Frame - Burgundy	\$3,239
☐ NVTC	Tilt Frame - Custom: _____	\$480

CLIENT WEIGHT		Select one below	3
☐ NVTWA	Up to 250Lbs		N/C
☐ NVTWAMF	Exceeds 250Lbs		\$876
☐ NVTWAMG	Exceeds 350Lbs		\$1,190
☐ NVTWAMH	Exceeds 400Lbs (Max 600Lbs)		\$1,332

SEAT WIDTH		Select one below	4
☐ NVTSW18	Seat Width 18"		N/C
☐ NVTSW19	Seat Width 19"		\$479
☐ NVTSW20	Seat Width 20"		\$479
☐ NVTSWC	Custom Seat Width: _____		\$479

SEAT DEPTH		Select one below	5
☐ NVTSD16	Seat Depth 16"		N/C
☐ NVTSD17	Seat Depth 17"		N/C
☐ NVTSD18	Seat Depth 18"		N/C
☐ NVTSD19	Seat Depth 19"		\$318
☐ NVTSD20	Seat Depth 20"		\$318
☐ NVTSDC	Custom Seat Depth: _____		\$318

FOOT REST AND FOOT PLATES		Select one below	6
☐ NVTOFST	Offset Foot Rest		N/C
☐ NVTF6A	60 Degree Hanger ☐ Pivot		N/C
☐ NVTF7A	70 Degree Hanger ☐ Pivot		N/C
☐ NVTF8A	80 Degree Hanger ☐ Pivot		N/C
☐ NVTEA	Elevating Legrest ☐ Pivot		\$485
☐ NVTAL	Articulating L/R ☐ Pivot		\$540
☐ NVTEB	Elevating Legrest System		\$540
☐ NVTPPL	Plastic Footplates		N/C
☐ NVTPXL	Extra Large Footplates		\$180
☐ NVTPAF	Angle Adjustable Footplates		\$288
☐ NVT1FB	One Piece Foot Board		\$96
☐ NVT1FBP	One Piece Foot Board Padding		\$53
☐ NVTNLS	Nylon Leg Strap Padded		\$72
☐ NVTNLSCL	Nylon Leg Strap Custom Size ()		\$144
☐ NVTCSL	Contracture Sling ☐ Standard ☐ M Style		\$288

BACK TYPE		Select one below	7
☐ NVTAB1	Angle Adjustable Back		N/C
☐ NVTENDO	Endomorph Canes W () ☐ Standard		N/C
☐ NVTRBK	Recline Kit		\$672

STANDARD NO CHARGE SPECIFICATIONS	
SEAT WIDTH: 14" – 18"	
SEAT DEPTH (ADJUSTABLE)	
SEAT HEIGHT (ADJUSTABLE): 16.75", 17.75" and 18.75"	
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]	
WEIGHT CAPACITY: 250 LBS	
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)	
TILT ADJUSTMENTS 0 TO 30 (DEGREES)	
TWO-WAY LOCKING GAS STRUT	
SINGLE RELEASE, DUAL CYLINDERS	
STANDARD LOCKOUT COLLARS	
STANDARD T STYLE ARMS	
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS	

TILT OPTION		Select one below	8
☐ NVTINTS	Independant Power Tilt System		\$2,400

SEAT TO FLOOR SPECIFICATIONS (± 0.5")			
(LOW SEAT SETTING)			
STF	R.WHEEL	CASTER	FORK

☐	16.75"	20"	6"	LONG
☐	17.75"	22"	6-8"	LONG
☐	18.75"	24"	6-8"	LONG

HEADREST		Select one below	9
☐	NVTMAHF	Multi-Axis Hardware Complete with Large Pad	\$614
☐	NVTMAHFXL	Multi-Axis Hardware Complete with Flex Large	\$614
☐	NVTMAHKP	Multi-Axis Hardware Complete with Medium Pad	\$614
☐	NVTMAHFXS	Multi-Axis Hardware Complete with Flex Medium	\$614
☐	NVTMAHFT	Multi-Axis Hardware Complete with Small Flat	\$614
☐	NVTUPHT219	Contour Headrest Complete with Swivel Hardware	\$520
☐	NVTUPHT218	Contour Headrest with Std. Assembly (1/2" Tube)	\$350
☐	NVTUPHT217	Contour Headrest with Std. Assembly (3/4" Tube)	\$265
☐	NVTUHRA	Universal Hardware Adaptor (Adaptor only)	\$103

ARM REST		Select one below	10
(Ht. 10"-14")	☐ NVTTA	T-Arm Ht. Adj./ Removable	\$288
(Ht. 10"-14")	☐ NVTTAS	T-Arm Ht. Adj./ Removable Short	\$288
	☐ NVTGWAPF	* Gel Waterfall Arm Pads	\$216

*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays and full laptrays if required

REAR WHEEL		Select one below	11
☐	NVTRWT0	20" Solid Treaded Tire	N/C
☐	NVTRWT2	22" Solid Treaded Tire	N/C
☐	NVTRWT4	24" Solid Treaded Tire	N/C
☐	NVTRWA2	22" Air Tire	\$144
☐	NVTRWA4	24" Air Tire	\$144
☐	NVTAQR	Axles – Quick Release	\$144
☐	NVTHRS	Standard Handrim Config. (3/4")	N/C
☐	NVTHRN	Narrow Handrim Config. (1/2" less)	N/C
☐	NVTHRT	Tight handrim Config. (No space)	N/C

*Note: Quick Release is not available for wheels smaller than 20"

CASTER		Select one below	12
☐	NVTF1B	6" X 1.25" Caster	N/C
☐	NVTF1C	8" X 1.25" Caster	N/C
☐	NVTF2B	6" X 1.75" Caster	\$96
☐	NVTF2C	8" X 1.75" Caster	\$96
☐	NVTF4B	6" X 2 Caster (Requires Steel Fork)	\$96
☐	NVTF4C	8" X 2 Caster (Requires Steel Fork)	\$96

WHEEL LOCK		Select one below	13
☐	NVTALPH	Aluminium Push To Lock	N/C
☐	NVTALPL	Aluminium Pull To Lock	N/C
☐	NVTBRXC	4" Brake Extension	N/C
☐	NVTBRXA	6" Brake Extension	N/C
☐	NVTBRXB	9" Brake Extension	\$156
☐	NVTATTBRX	Attendant Brakes	\$288
☐	NVTULB	Unilateral Hand Brake ☐ Left ☐ Right	\$144

SEATING - BACKS		Custom: select back & fill custom	14
		W () x H ()	
☐	NVTBK01	RELAX Gentle Fit Back	\$653
☐	NVTBK02	RELAX Full Fit Back	\$730
☐	NVTBK03	RELAX Snug Fit Back	\$840
☐	NVTBK04	RELAX Form Fit Back No Laterals	\$984
☐	NVTBK06	RELAX Form Fit Back Gentle Laterals	\$984
☐	NVTBK07	RELAX Form Fit Back Full Laterals	\$984
☐	NVTBK05	RELAX Form Fit Back Snug Laterals	\$984
☐	NVTBK10	RELAX Spiral Fit Back Gentle Laterals	\$998
☐	NVTBK09	RELAX Spiral Fit Back Snug Laterals	\$998
☐	NVTBK11	RELAX Spiral Fit Back Full Laterals	\$998
☐	NVTTAB	RELAX Tension Back Upto 18"	\$168
☐	NVTTAB+	RELAX Tension Back Over 18"	\$288
☐	NVTBK12	* Custom Width : _____	\$211
☐	NVTBK13	* Custom Height : _____	\$211

*Note: Please refer to the Seating Booklet for custom sizes

☐	NVTBGL1	Gel Full Overlay	\$696
☐	NVTBGL2	Gel Center Vertebrae	\$418
☐	NVTBGL3	Gel Left Lateral	\$139
☐	NVTBGL4	Gel Right Lateral	\$139

SEATING - CUSHIONS		Select one below	15
		W () x D ()	
☐	NVTC01	RELAX Soft Cushion	\$576
☐	NVTC02	RELAX Comfy Cushion	\$749
☐	NVTC03	RELAX Firm Cushion	\$941
☐	NVTC04	RELAX AIR FLOW Cushion	\$704
☐	NVTCGL2	Full Overlay Gel	\$430
☐	NVTCCL	Contoured Leg Channel (Bolster) Std.	\$216
☐	NVTCCLA	Contoured Leg Channel (Bolster) with ABS	\$312
☐	NVTCCLC	* Contoured Leg Channel (Bolster) Custom	\$384

*Note: Please refer to the Seating Booklet for custom sizes

☐	NVTC00	Extra Cushion Cover	\$119
--------------------------	--------	---------------------	-------

ACCESSORIES		16	
☐	NVTSBA	Seat Belt (Auto Style)	\$31
☐	NVTSBAP	Seat Belt (Padded)	\$67
☐	NVTSBA4P	Seat Belt (4 Point Padded)	\$156
☐	NVTTRA	Full Padded Tray Locking	\$336
☐	NVTTRB	Full Padded Tray Slide On	\$336
☐	NVTTRCL	Full Plastic Tray Locking ☐ ABS Black ☐ Clear	\$336
☐	NVTTRC	Full Plastic Tray Slide On ☐ ABS Black ☐ Clear	\$336
☐	NVTTRD	Half Padded Left Tray	\$288
☐	NVTTRE	Half Padded Right Tray	\$288
☐	NVTTRF	Half Left Tray ☐ ABS Black ☐ Clear	\$288
☐	NVTTRG	Half Right Tray ☐ ABS Black ☐ Clear	\$288
☐	NVTDFB	Transport Tie Downs Front	\$102
☐	NVTDRB	Transport Tie Downs Rear	\$102
☐	NVTOXY	Oxygen Tank Holder	\$234
☐	NVTIVA	Adjustable I V Pole	\$156
☐	NVTCHPT-L	Left Arm Trough Tray	\$480
☐	NVTCHPT-R	Right Arm Trough Tray	\$480
☐	NVTLATB	Lateral Bite out	\$216
☐	NVTADD8	8" Adductor	\$197
☐	NVTADD10	10" Adductor	\$197
☐	NVTADD12	12" Adductor	\$197
☐	NVTFRCH	Padded Footplates	\$106