

EXT Extreme Wheelchair (Tilt Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



The Extreme tilt wheelchair has a light weight frame and is depth adjustable. The tilt system offers users a change in orientation to redistribute pressure from one area to another. The front seat to floor height also remains constant during a 0 to 25 degrees tilt

SEAT TO FLOOR SPECIFICATIONS			
STF	R.WHEEL	CASTER	FORK
SL 13.00"	20"	5"	XL SHORT

SPECIFICATIONS						1
Item	CODE	Description	Item	CODE	Description	
Colour	EXTB	Dynamic Frame - Black	Caster	EXTFC1A	5" X 1.25" Caster	
Client Weight	EXTWA	Up to 250 lbs.	Foot Hanger	EXTF7A	70 Degree Hanger - Pivot Style	
Seat Width	EXTSWXX	<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"	Foot Plates	EXTPAF	Angle Adjustable Footplates	
Seat Depth	EXTSD18	Seat Depth 18"	Wheel Lock	EXTALPH	Aluminium Push To Lock	
Back Type	EXTAB1	Angle Adjustable Back	Wheel Lock Extension	EXTBRXA	6" Brake Extension	
Arm Rest	EXTTA	T-Arm Height Adjustable / Removable	Rear Wheels	EXTRWTO	20" Solid Treaded Tire	

ADDITIONAL OPTIONS 2

- EXTBK01 **RELAX** Gentle Fit Back 544
- EXTMAHF Multi-Axis Hardware Complete with Large Pad 512
- EXTMAHFXL Multi-Axis Hardware Complete with Flex Large 512
- EXTMAHKP Multi-Axis Hardware Complete with Medium Pad 512
- EXTMAHFXS Multi-Axis Hardware Complete with Flex Medium 512
- EXTMAHFT Multi-Axis Hardware Complete with Small Flat 512

***** PLEASE NOTE: *****

ONLY "ADDITIONAL OPTIONS" MENTIONED ABOVE CAN BE ADDED TO THE QUICK SHIP ORDER PRODUCTS. IF YOU REQUIRE FURTHER MODIFICATION PLEASE FILL OUT OUR REGULAR EXTREME ORDER FORM.

MINIMUM OF 2 UNITS ARE REQUIRED FOR SAME DAY OR NEXT DAY PROCESSING AND SHIPPING.

NOTES:

Number of 16" Seat Width Wheelchairs: _____
Number of 18" Seat Width Wheelchairs: _____
Number of 20" Seat Width Wheelchairs: _____
Total Wheelchairs Required: _____