

STP STP Wheelchair (Tilt Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

- STPB Tilt Frame - Black 2699
- STPS Tilt Frame - Silver 2699
- STPBL Tilt Frame - Blue 2699
- STPBU Tilt Frame - Burgundy 2699
- STPC Tilt Frame - Custom : _____ 3099

CLIENT WEIGHT *Select one below* **2**

- STPWA Up to 250Lbs N/C
- STPWAMF Exceeds 250Lbs 730
- STPWAMG Exceeds 350Lbs 992
- STPWAMH Exceeds 400Lbs (Max 600Lbs) 1110

SEAT WIDTH *Select one below* **3**

- STPSW16 Seat Width 16" N/C
- STPSW17 Seat Width 17" N/C
- STPSW18 Seat Width 18" N/C
- STPSW19 Seat Width 19" 399
- STPSW20 Seat Width 20" 399
- STPSWC Custom Seat Width: _____ 399

SEAT DEPTH *Select one below* **4**

- STPSD16 Seat Depth 16" N/C
- STPSD17 Seat Depth 17" N/C
- STPSD18 Seat Depth 18" N/C
- STPSD19 Seat Depth 19" 265
- STPSD20 Seat Depth 20" 265
- STPSDC Custom Seat Depth: _____ 265

FOOT REST AND FOOT PLATES *Select one below* **5**

- STPOFST Offset Foot Rest N/C
- STPF6A 60 Degree Hanger Pivot Pin N/C
- STPF7A 70 Degree Hanger Pivot Pin N/C
- STPF8A 80 Degree Hanger Pivot Pin N/C
- STPEA Elevating Legrest Pivot Pin 404
- STPAL Articulating L/R Pivot Pin 450
- STPEB Elevating Legrest System 450
- STPPPL Plastic N/C
- STPPXL Extra Large 150
- STPPAF Angle Adjustable Footplates 240
- STPFRCH Padded Footplates 88
- STP1FB One Piece Foot Board 80
- STP1FBP One Piece Foot Board Padding 44
- STPNLS Nylon Leg Strap Padded 60
- STPNLSC Nylon Leg Strap Custom Size () 120
- STPCSCL Contracture Sling Standard M Style 240

STANDARD NO CHARGE SPECIFICATIONS

SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS
STANDARD T STYLE ARMS
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

Tilt Options **6**

- STPINTS Independant Power Tilt System 2000

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	
(STANDARD SEAT SETTING)				
<input type="checkbox"/> MH 16.75"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	8"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	24"	6"	LONG	
<input type="checkbox"/> A 18.25"	24"	8"	LONG	

Notes:

BACK TYPE		7	SEATING - BACKS		13	
<input type="checkbox"/>	STPAB1	Angle Adjustable Back	N/C	W () x H ()		
<i>Specify Width</i> <input type="checkbox"/>	STPENDO	Endomorph Canes	N/C	<input type="checkbox"/>	STPBK01 RELAX Gentle Fit Back 544	
<input type="checkbox"/>	STPRBK	Recline Kit	560	<input type="checkbox"/>	STPBK02 RELAX Full Fit Back 608	
HEADREST		8	<i>For custom, choose one back then fill custom</i>			
<input type="checkbox"/>	STPMAHF	Multi-Axis Hardware Complete with Large Pad	512	<input type="checkbox"/>	STPBK03 RELAX Snug Fit Back 700	
<input type="checkbox"/>	STPMAHFXL	Multi-Axis Hardware Complete with Flex Large	512	<input type="checkbox"/>	STPBK04 RELAX Form Fit Back No Laterals 820	
<input type="checkbox"/>	STPMAHKP	Multi-Axis Hardware Complete with Medium Pad	512	<input type="checkbox"/>	STPBK06 RELAX Form Fit Back Gentle Laterals 820	
<input type="checkbox"/>	STPMAHFXS	Multi-Axis Hardware Complete with Flex Medium	512	<input type="checkbox"/>	STPBK07 RELAX Form Fit Back Full Laterals 820	
<input type="checkbox"/>	STPMAHFT	Multi-Axis Hardware Complete with Small Flat	512	<input type="checkbox"/>	STPBK05 RELAX Form Fit Back Snug Laterals 820	
<input type="checkbox"/>	STPUPHT219	Contour Headrest Complete with Swivel Hardware	433	<input type="checkbox"/>	STPBK10 RELAX Spiral Fit Back Gentle Laterals 832	
<i>(1/2" Tubing)</i> <input type="checkbox"/>	STPUPHT218	Contour Headrest with Std. Assembly (1/2")	292	<input type="checkbox"/>	STPBK09 RELAX Spiral Fit Back Snug Laterals 832	
<i>(3/4" Tubing)</i> <input type="checkbox"/>	STPUPHT217	Contour Headrest with Std. Assembly (3/4")	292	<input type="checkbox"/>	STPBK11 RELAX Spiral Fit Back Full Laterals 832	
<input type="checkbox"/>	STPUHA	Universal Hardware Adaptor (Adaptor only)	57	<input type="checkbox"/>	STPTAB RELAX Tension Back Upto 18" 140	
<input type="checkbox"/>	STPMAHGA-P	Multi Axis with Headrest Plate (No Pad)	255	<input type="checkbox"/>	STPTAB+ RELAX Tension Back Over 18" 240	
ARM REST		9	<input type="checkbox"/>	STPBGL1	Half Overlay Gel (3 Piece) 348	
<i>(Std 10"-14")</i> <input type="checkbox"/>	STPTA	T-Arm Ht. Adj./ Removable	N/C	<input type="checkbox"/>	STPBGL2	Full Overlay Gel (6 Pieces) 696
<i>(Std 10"-14")</i> <input type="checkbox"/>	STPFA	Flip Back Ht. Adj./ Removable	340	<input type="checkbox"/>	STPPSTP	Pin Extension (1" Extension) 40
<i>(Short 7"-11")</i> <input type="checkbox"/>	STPTAS	T-Arm Ht. Adj./ Removable	N/C	<input type="checkbox"/>	STPBK12	* Custom Width : _____ +176
<i>(Short 7"-11")</i> <input type="checkbox"/>	STPFAS	Flip Back Ht. Adj./ Removable	340	<input type="checkbox"/>	STPBK13	* Custom Height : _____ +176
<input type="checkbox"/>	STPGWAPP	* Gel Waterfall Arm Pads	180	*Note: Please refer to the Seating Booklet for custom sizes		
*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required			SEATING - CUSHION		14	
REAR WHEELS		10	<i>Select one below</i>			
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C	<input type="checkbox"/>	STPC01 RELAX Soft Cushion 480	
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C	<input type="checkbox"/>	STPC02 RELAX Comfy Cushion 624	
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C	<input type="checkbox"/>	STPC03 RELAX Firm Cushion 784	
<i>(STF 13-19")</i> <input type="checkbox"/>	STPRWT12	12" Solid Tire	N/C	<input type="checkbox"/>	STPCGL1	Half Overlay Gel (1 Piece) 116
<i>(STF 17-19")</i> <input type="checkbox"/>	STPRWT12I	12" Solid Tire Inside Mount	120	<input type="checkbox"/>	STPCGL2	Full Overlay Gel (2 Pieces) 232
<input type="checkbox"/>	STPRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	STPCCL	Contoured Leg Channel (Bolster) Std. 180
<input type="checkbox"/>	STPRWA2	22" Air Tire	120	<input type="checkbox"/>	STPCCLA	Contoured Leg Channel (Bolster) with ABS 260
<input type="checkbox"/>	STPRWA4	24" Air Tire	120	<input type="checkbox"/>	STPCCLC	* Contoured Leg Channel (Bolster) Custom 320
<input type="checkbox"/>	STPAQR	Axles – Quick Release	120	*Note: Please refer to the Seating Booklet for custom sizes		
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	STPC00	Extra Cushion Cover 99
CASTER		11	ACCESSORIES		15	
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style) 26
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBAP	Seat Belt (Padded) 56
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBA4P	Seat Belt (4 Point Padded) 130
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster (Soft Roll)	80	<input type="checkbox"/>	STPTRA	Full Padded Tray Locking 280
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	80	<input type="checkbox"/>	STPTRB	Full Padded Tray Slide On 280
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	80	<input type="checkbox"/>	STPTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 280
<input type="checkbox"/>	STPFC4B	6" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	STPTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 280
<input type="checkbox"/>	STPFC4C	8" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	STPTRD	Half Padded Left Tray 240
WHEEL LOCK		12	<input type="checkbox"/>	STPTRE	Half Padded Right Tray 240	
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	STPTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	STPTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C	<input type="checkbox"/>	STPTDFB	Transport Tie Downs Front 85
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C	<input type="checkbox"/>	STPTDRB	Transport Tie Downs Rear 85
<input type="checkbox"/>	STPBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	STPOXY	Oxygen Tank Holder 195
<input type="checkbox"/>	STPBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	STPIVA	Adjustable I V Pole 130
<input type="checkbox"/>	STPBRXB	9" Brake Extension	130	<input type="checkbox"/>	STPCHPT-L	Left Arm Trough Tray 400
<input type="checkbox"/>	STPATTBRX	Attendant Brakes	240	<input type="checkbox"/>	STPCHPT-R	Right Arm Trough Tray 400
				<input type="checkbox"/>	STPLATB	Lateral Bite out 180
				<input type="checkbox"/>	STPADD8	8" Adductor 164
				<input type="checkbox"/>	STPADD10	10" Adductor 164
				<input type="checkbox"/>	STPADD12	12" Adductor 164