

TNP Transit Wheelchair (Manual Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR		Select one below	1
<input type="checkbox"/>	TNPS	Frame - Grey Steel Frame	1000
<input type="checkbox"/>	TNPB	Frame - Black	1000

CLIENT WEIGHT		Select one below	2
<input type="checkbox"/>	TNPWA	Upto 250 lbs	N/C
<input type="checkbox"/>	TNPWAMF	Exceeds 250 lbs	300
<input type="checkbox"/>	TNPWAMG	Exceeds 350 lbs (Max 600 lbs)	500

SEAT WIDTH		Select one below	3
<input type="checkbox"/>	TNPSW16	Seat Width 16"	N/C
<input type="checkbox"/>	TNPSW18	Seat Width 18"	N/C
<input type="checkbox"/>	TNPSW20	Seat Width 20"	300
<input type="checkbox"/>	TNPSW22	Seat Width 22"	300

SEAT DEPTH		Select one below	4
<input type="checkbox"/>	TNPSD18	Seat Depth 18"	N/C
<input type="checkbox"/>	TNPSD20	Seat Depth 20"	300

FOOT REST AND FOOT PLATES		Select one below	7
<input type="checkbox"/>	TNPF7A	70 Degree Hanger	N/C
<input type="checkbox"/>	TNPEA	Elevating Legrest	404
<input type="checkbox"/>	TNPAL	Articulating Leg Rests	450

FOOT REST TYPE		Select one below	8
<input type="checkbox"/>	TNPFL	Fixed Legrest	N/C
<input type="checkbox"/>	TNPRL	Removable Legrest	N/C

ARM REST		Select one below	9
<input type="checkbox"/>	TNPFAF	Fixed Height - Flip Back	N/C
<input type="checkbox"/>	TNPFAFN	Fixed Height - Non Flip Back	N/C

WHEEL LOCK		Select one below	10
<input type="checkbox"/>	TNPSTPH	Steel Push To Lock	N/C
<input type="checkbox"/>	TNPSTPL	Steel Pull To Lock	N/C
<input type="checkbox"/>	TNPBRXC	4" Brake Extension	N/C
<input type="checkbox"/>	TNPBRXA	6" Brake Extension	N/C
<input type="checkbox"/>	TNPBRXB	9" Brake Extension	130

SEATING - Upholstery		Select one below	11
<input type="checkbox"/>	TNPUP01	Blue Vinyl Upholstery	N/C
<input type="checkbox"/>	TNPUP02	Pink Vinyl Upholstery	N/C
<input type="checkbox"/>	TNPUP03	Green Vinyl Upholstery	N/C
<input type="checkbox"/>	TNPUP04	Black Vinyl Upholstery	N/C
<input type="checkbox"/>	TNPUP05	Purple Vinyl Upholstery	N/C
<input type="checkbox"/>	TNPUPC	Custom Vinyl Upholstery	N/C

Specify: _____

Logo			12
<input type="checkbox"/>	TNPLOGO	Custom Company Logo	100

SPECIFICATIONS			
STEEL FRAME AND CROSS BRACES			
STANDARD 18" BACK HEIGHT			
EXTRUDED ALUMINUM FORKS			
CHROME HAND RIMS			
VINYL UPHOLSTERY			
POWDER COATED FRAME			
FIXED HEIGHT - FLIP BACK OR FIXED HEIGHT - NON FLIP BACK			
SWING AWAY FOOTREST WITH HEEL LOOPS			
PRECISION BEARINGS			
MOULDED PLASTIC FOOTPLATES			
ANTI-TIP ROLLERS			
6" BRAKE EXTENSIONS			
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS			

SEAT TO FLOOR SPECIFICATIONS			
(STANDARD SEAT SETTING)			
<input type="checkbox"/>	H 15.50"	22"	6" LONG
<input type="checkbox"/>	H 17.50"	24"	8" LONG

ACCESSORIES		13	
<input type="checkbox"/>	TNPOXYIV	Oxygen Tank and Adjustable I V Pole Combo	325

Notes: