

VGP Voyager Wheelchair - Folding Wheelchair



INFORMATION

Customer Name: _____
Address: _____ PO #: _____
City: _____ Phone #: _____ Date: _____
Account #: _____ Authorizer Name: _____

FRAME COLOUR		Select one below	1
☐	VGPS	Frame - Silver	\$2,028
☐	VGPB	Frame - Black	\$2,028
☐	VGPBL	Frame - Blue	\$2,028
☐	VGPBU	Frame - Burgundy	\$2,028
☐	VGPC	Frame - Colour : _____	\$2,508

CLIENT WEIGHT		Select one below	2
☐	VGPWA	Up to 250Lbs	N/C
☐	VGPWAMF	Exceeds 250Lbs	\$876
☐	VGPWAMG	Exceeds 350Lbs	\$1,190
☐	VGPWAMH	Exceeds 400Lbs (Max 600Lbs)	\$1,332

SEAT WIDTH		Select one below	3
☐	VGPSW14	Seat Width 14"	N/C
☐	VGPSW15	Seat Width 15"	N/C
☐	VGPSW16	Seat Width 16"	N/C
☐	VGPSW17	Seat Width 17"	N/C
☐	VGPSW18	Seat Width 18"	N/C
☐	VGPSW19	Seat Width 19"	\$479
☐	VGPSW20	Seat Width 20"	\$479
☐	VGPSWS	Custom Seat Width: _____	\$479

SEAT DEPTH		Select one below	4
☐	VGPSD14	Seat Depth 14"	N/C
☐	VGPSD15	Seat Depth 15"	N/C
☐	VGPSD16	Seat Depth 16"	N/C
☐	VGPSD17	Seat Depth 17"	N/C
☐	VGPSD18	Seat Depth 18"	N/C
☐	VGPSD19	Seat Depth 19"	\$318
☐	VGPSD20	Seat Depth 20"	\$318
☐	VGPSDS	Custom Seat Depth: _____	\$318

BACK HEIGHT		Select one below	5
☐	VGPBH15	Back Height 15"	N/C
☐	VGPBH16	Back Height 16"	N/C
☐	VGPBH17	Back Height 17"	N/C
☐	VGPBH18	Back Height 18" (Standard)	N/C
☐	VGPBH19	Back Height 19"	N/C
☐	VGPBH20	Back Height 20"	N/C
☐	VGPBH21	Back Height 21"	N/C
☐	VGPBHC	Back Height Custom : _____	\$144

STANDARD NO CHARGE SPECIFICATIONS
HEAVY DUTY STEEL FRAME AND CROSS BRACES
STANDARD 18" BACK HEIGHT
20", 22", OR 24" NYLON MOULDED REAR WHEELS
EXTRUDED ALUMINUM FORKS
CHROME HAND RIMS
5", 6", OR 8" CASTERS
BLACK NYLON UPHOLSTERY
POWDER COATED FRAME
FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
PRECISION BEARINGS
MOULDED PLASTIC FOOTPLATES
SEAT WIDTH 14" TO 18"
VERTICAL/HORIZONTAL REAR AXLE PLATE ADJUSTABILITY
ANTI-TIP ROLLERS
6" BRAKE EXTENSIONS
DOUBLE CROSS BRACE
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
	STF	R.WHEEL	CASTER	FORK
☐	12.75"	20"	5"	2H-N
☐	13"	20"	5"	3H,2H,1H
☐	13.75"	20"	5"	3H,2H
☐	14.25"	20"	6"	5H,3H,2H
☐	14.5"	22"	5"	5H,3H
☐	14.25"	22"	6"	5H,3H,2H
☐	15.25"	20"	5"	5H
☐	15"	20"	6"	5H,3H
☐	15"	22"	5"	5H
☐	15"	22"	6"	5H,3H
☐	16.25"	22"	6"	5H
☐	16.5"	22"	8"	5H
☐	17.25"	22"	8"	5H
(HEIGHT EXTENDED CASTER)				
	STF	R.WHEEL	CASTER	FORK
☐	18.25"	24"	8"	5H,3H

For more Seat To Floor combinations, visit our website

Notes:

FOOT REST AND FOOT PLATES		Select one below	6	SEATING -BACKS		Custom: select back & fill custom	12		
☐	VGPF6A	60 Degree Hanger Pivot	N/C		W () x H ()				
☐	VGPF7A	70 Degree Hanger Pivot	N/C	☐	VGPBK01	RELAX Gentle Fit Back	\$653		
☐	VGPF8A	80 Degree Hanger Pivot	N/C	☐	VGPBK02	RELAX Full Fit Back	\$730		
☐	VGPEA	Elevating Legrest Pivot	\$485	☐	VGPBK03	RELAX Snug Fit Back	\$840		
☐	VGPAL	Articulating Leg Rests Pivot	\$540	☐	VGPBK04	RELAX Form Fit Back No Laterals	\$984		
☐	VGPPPL	Plastic	N/C	☐	VGPBK06	RELAX Form Fit Back Gentle Laterals	\$984		
☐	VGPPXL	Extra Large Foot Rest	\$180	☐	VGPBK07	RELAX Form Fit Back Full Laterals	\$984		
☐	VGPPAF	Angle Adjustable Footplates	\$288	☐	VGPBK05	RELAX Form Fit Back Snug Laterals	\$984		
☐	VGPIFB	One Piece Foot Board	\$96	☐	VGPBK10	RELAX Spiral Fit Back Gentle Laterals	\$998		
☐	VGPIFBP	One Piece Foot Board Padding	\$53	☐	VGPBK09	RELAX Spiral Fit Back Snug Laterals	\$998		
☐	VGPNLS	Nylon Leg Strap Padded	\$72	☐	VGPBK11	RELAX Spiral Fit Back Full Laterals	\$998		
☐	VGPNLSC	Nylon Leg Strap Custom Size ()	\$144	☐	VGPTAB	RELAX Tension Back Upto 18"	\$168		
☐	VGPCSL	Contracture Sling ☐ Standard ☐ M Style	\$288	☐	VGPTAB+	RELAX Tension Back Over 18"	\$288		
BACK TYPE		Select one below	7	☐	VGPBK12	* Custom Width : _____	\$211		
☐	VGAB1	Angle Adjustable Back	N/C	☐	VGPBK13	* Custom Height : _____	\$211		
ARM REST		Select one below	8	*Note: Please refer to the Seating Booklet for custom sizes					
(Ht. 10"-14")	☐	VGPF6A	Flip Back Ht. Adj./ Removable	\$408	☐	VGPBGL1	Gel Full Overlay	\$696	
(Ht. 10"-14")	☐	VGPTA	T-Arm Ht. Adj./ Removable	\$288	☐	VGPBGL2	Gel Center Vertebrae	\$418	
(Ht. 7"-11")	☐	VGPFAS	Flip Back Ht. Adj./ Removable Short	N/C	☐	VGPBGL3	Gel Left Lateral	\$139	
(Ht. 7"-11")	☐	VGPTAS	T-Arm Ht. Adj./ Removable Short	\$408	☐	VGPBGL4	Gel Right Lateral	\$139	
(Ht. 10"-14")	☐	VGPP3PA	3 Point Arm. Removable/No Flip	\$408	SEATING - CUSHIONS			Select one below	13
(Ht. 7"-11")	☐	VGPP3PAS	3 Point Arm. Removable/No Flip Short	\$408		W () x D ()			
☐	VGPGWAPF	* Gel Waterfall Arm Pads	\$216	☐	VGPC01	RELAX Soft Cushion	\$576		
* Note: Gel Waterfall Arm Pads are only compatible with locking laptray and full laptrays if required				☐	VGPC02	RELAX Comfy Cushion	\$749		
REAR WHEEL		Select one below	9	☐	VGPC03	RELAX Firm Cushion	\$941		
☐	VGPRWT0	20" Solid Treaded Tire	N/C	☐	VGPC04	RELAX AIR FLOW Cushion	\$704		
☐	VGPRWT2	22" Solid Treaded Tire	N/C	☐	VGPCGL2	Full Overlay Gel (2 Pieces)	\$430		
☐	VGPRWT4	24" Solid Treaded Tire	N/C	☐	VGPCCL	Contoured Leg Channel (Bolster) Std.	\$216		
☐	VGPRWA2	22" Air Tire	\$144	☐	VGPCCLA	Contoured Leg Channel (Bolster) with ABS	\$312		
☐	VGPRWA4	24" Air Tire	\$144	☐	VGPCCLC	* Contoured Leg Channel (Bolster) Custom	\$384		
☐	VGPAQR	Axles – Quick Release	\$144	*Note: Please refer to the Seating Booklet for custom sizes					
CASTER		Select one below	10	☐	VGPC00	Extra Cushion Cover	\$119		
☐	VGPF1A	5" X 1.25" Caster	N/C	ACCESSORIES			14		
☐	VGPF1B	6" X 1.25" Caster	N/C	☐	VGPSBA	Seat Belt (Auto Style)	\$31		
☐	VGPF1C	8" X 1.25" Caster	N/C	☐	VGPSBAP	Seat Belt (Padded)	\$67		
☐	VGPF2A	5" X 1.75" Caster (Soft Roll)	\$96	☐	VGPSBA4P	Seat Belt (4 Point Padded)	\$156		
☐	VGPF2B	6" X 1.75" Caster	\$96	☐	VGPTRA	Full Padded Tray Locking	\$336		
☐	VGPF2C	8" X 1.75" Caster	\$96	☐	VGPTRB	Full Padded Tray Slide On	\$336		
☐	VGPF4B	6" X 2 Caster (Requires Steel Fork)	\$96	☐	VGPTRCL	Full Plastic Tray Locking ☐ ABS Black ☐ Clear	\$336		
☐	VGPF4C	8" X 2 Caster (Requires Steel Fork)	\$96	☐	VGPTRC	Full Plastic Tray Slide On ☐ ABS Black ☐ Clear	\$336		
WHEEL LOCK		Select one below	11	☐	VGPTRD	Half Padded Left Tray	\$288		
☐	VGPA1PH	Aluminium Push To Lock	N/C	☐	VGPTRE	Half Padded Right Tray	\$288		
☐	VGPA1PL	Aluminium Pull To Lock	N/C	☐	VGPTRF	Half Left Tray ☐ ABS Black ☐ Clear	\$288		
☐	VGPS1PH	Steel Push To Lock	N/C	☐	VGPTRG	Half Right Tray ☐ ABS Black ☐ Clear	\$288		
☐	VGPS1PL	Steel Pull To Lock	N/C	☐	VGPTDFB	Transport Tie Downs Front	\$102		
☐	VGPB1RX	4" Brake Extension	N/C	☐	VGPTDRB	Transport Tie Downs Rear	\$102		
☐	VGPB1RXA	6" Brake Extension	N/C	☐	VGPOXY	Oxygen Tank Holder	\$234		
☐	VGPB1RXB	9" Brake Extension	\$156	☐	VGPIVA	Adjustable I V Pole	\$156		
				☐	VGPC1PT-L	Left Arm Trough Tray	\$480		
				☐	VGPC1PT-R	Right Arm Trough Tray	\$480		
				☐	VGPLATB	Lateral Bite out	\$216		
				☐	VGPA1DD8	8" Adductor	\$197		
				☐	VGPA1DD10	10" Adductor	\$197		
				☐	VGPA1DD12	12" Adductor	\$197		
				☐	VGPF1RCH	Padded Footplates	\$106		