

RELAX

SPIRAL FIT BACK



Rear Back Angle



5 Points of Tension Adjustability

The **Relax Spiral Fit** by Power Plus Mobility is designed to meet a variety of positioning needs. This is an aluminum light weight formed back with hardware that allows for multiple adjustments. The Relax Spiral Fit with its wide range of laterals is suited for users needing complex postural support and spinal adjustments. Accommodates for curvature of the spine allowing for comfortable posterior pelvic tilt.

Quick reference guide

RELAX SERIES - SPIRAL FIT BACK

REGULAR SIZES ADP: SEBWM0030

Width = 16" - 18"

Part Number: Refer below

Height = 18"

CUSTOM ADP: WIDTH SEBWM0035, HEIGHT SEBWM0040

Width = 15" or <, 19" or >

Part Number: Refer below

Height = 17" or <, 19" or >

"<" means less and ">" means more

PART NUMBERS

KTBF (W X H) **No Laterals**

KTBS (W X H) **Gentle Laterals**

KTBM (W X H) **Snug Laterals**

KTBD (W X H) **Full Laterals**

E.g. For a W = 15" and H = 17" Back; Part number would be KTB1517


RELAX

RELAX SPIRAL FIT COMES WITH THE FOLLOWING BACK SHELL OPTIONS




NO LATERAL DIMENSIONS

X = Standard shell height 18"




GENTLE SHELL LATERAL DIMENSIONS

X = Standard shell height 18"
Y = Lateral depth 3" finished
Z = Lateral height 10"



FULL SHELL LATERAL DIMENSIONS

X = Standard shell height 18"
Y = Lateral depth 5" finished
Z = Lateral height 10"




SNUG SHELL LATERAL DIMENSIONS

W = Standard shell height 18"
X = Lateral depth 5" finished
Y = Space below lateral 4"
Z = Lateral height 6" above cutout

KEY FEATURES

- Wide range of laterals with cutout that allows for curvature of the spine
- Comfortable High Resilience Foam
- 3lbs Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Adjustable/Removable Lumbar Support
- Easy to Install and Remove
- Height Adjustable

OVERLAYS



ADP: SEBND0015 SEBND0015 SEBND0020

* Up-charge items

SPECIFICATIONS	RELAX SERIES - SPIRAL FIT BACK
ADP	SEBWM0030 , CUSTOM W SEBWM0035 , H SEBWM0040
STANDARD WIDTH RANGE	16" - 18"
CUSTOM WIDTH RANGE	15" or less and 19" or more
STANDARD HEIGHT RANGE	18"
CUSTOM HEIGHT RANGE	17" or less and 19" or more
AVAILABLE HARDWARE STYLES	One style which is universal to all our backs
ANGLE ADJUSTABLE	25° anterior adjustability, 20° of posterior adjustability
WIDTH ADJUSTABLE	2" of adjustability through the pins on the back
FITS TUBING SIZE	7/8" tubing or 1"
HEADREST MOUNT	Backs pre-drilled for headrest supports and lateral supports