

# EXTREME

TILT WHEELCHAIR



ADP CODE: WM0001041

AADL CODE: W305

The **EXTREME** tilt Wheelchair from Power Plus Mobility has a light weight frame and is depth adjustable. The **EXTREME** tilt system provides a change of orientation to redistribute pressure from one area to another. The front seat to floor height also remains constant during 0 To 25 degrees of tilt.




Seat To Floor Height  
Adjustability



Zero Maintenance Soft Set Cylinders



Offset Front Hangers  
(Optional)



Dual Cylinders Allows For Consistent Tilt For All Weight Capacities



Standard Angle Adjustable Stroller Handles



Angle Adjustable Front Journals

## ABOUT EXTREME

Our wheelchair **EXTREME** lives up to its name. Exceeding the boundaries of what you would find in a wheelchair, the **EXTREME** is in a category all on its own. When Power Plus Mobility first began the process of creating **EXTREME**, it was obvious that this product would deliver extraordinary results. With its ability to be fully adjustable and easily customizable, the **EXTREME** is perfect for a variety of users. We built it with our customers in mind and every detail has been fully analyzed and researched. The **EXTREME** is a prime example of the hard work we put into our products.

## CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# EXTREME



## STANDARD FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts for Easy Tilt (Up To 25° )
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Relax Seating System
- New Gas Cylinder Release Mechanism
- Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals
- New Recessed Axles
- Maintains Front STF Height In Tilt
- Overall Wheelchair Width: 29" (18" Seat Width)

SPECIFICATIONS	EXTREME
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	(ADJUSTABLE): 16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 12.25" TO 18.25"
Tilt ANGLE:	0° TO 25°
BACK ANGLE:	(ADJUSTABLE): 0°, 5°, 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	DYNAMIC TILT ADJUSTMENTS 0° TO 25°
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS

**EXTREME**  
TILT WHEELCHAIR



Learn more about this product by scanning the QR Code