

HORIZON TILT

TILT COMMODORE

STAINLESS STEEL



AADL CODE: G014 / G067

The HORIZON TILT gives users the ability to comfortably carry out their daily personal care routine safely and effectively. The HORIZON TILT serves similar functions as the Horizon regular commode but has the added bonus of the tilt feature. This option gives users a risk-free change in seating positioning.

HORIZON TILT 1

All Features Of Horizon Are The Following :

- Weight capacity 250 lbs.
- 18" x 18" Easy to clean urethane seat
- Seat width by frame 18"
- Overall width 23.5"
- Rear axle receiver ready
- Silver powder coated finish
- Black urethane contour headrest
- 30° of tilt
- Black urethane easy to clean back rest
- Black molded sanitary cover
- Flip back arms
- Black urethane easy to clean arm pads
- 6" x 10" x 14" Large bucket
- 20" Back height
- 4" x 5" Rust resistant nylon casters (3 total locking casters and 1 directional lock caster)
- 70° Footrest
- Weight of chair: ≈ 46 lbs.

UP CHARGE OPTIONS

- 14" to 20" Seat depth
- 21" to 26" Seat depth
- 14" to 17" Seat width
- 23" to 26" Seat width
- Custom back height
- Capacity 750 to 1000 lbs.
- Tension back
- Lock down flip back arms
- Custom soft seat 18" to 22" (closed oval)
- Custom soft seat 18" to 22" (key hole)
- 60° Footrest
- 80° Footrest
- 90° Footrest
- Elevating legrest - not water proof
- Angle adjustable footplates - not water proof
- One piece footboard
- Seat belt (plastic)
- Height adjustable casters
- Custom seat to floor height
- Custom soft seats
- Self propelling wheels

Seat Heights

- 24" Wheels X 5" Casters = 21" Seat to Floor Height

Wheel Option: Self Propelling Wheels

- 24" Wheels/ 5" Casters

Warranty: 5 Years On Frame

- 1 Year on components



Learn more about this product by scanning the QR Code

HORIZON TILT 2

All Features Of Horizon 1 Plus The Following:

- Weight capacity 350 lbs.
- 18" X 18" Easy to clean urethane seat
- Seat width by frame 20"
- Overall width 25.5"
- Heavy duty construction
- Weight of chair: ≈ 49 lbs.

HORIZON TILT 3

All Features Of Horizon 1 Plus The Following:

- Weight capacity 450 lbs.
- 21" X 18" Easy to clean urethane seat
- Seat width by frame 22"
- Overall width 27.5"
- Extra heavy duty construction
- Weight of chair: ≈ 51 lbs.

