

RELAX SERIES

RELAX

SNUG FIT BACK



MSRP: \$846



The **Relax Snug Fit** is an aluminum light weight formed back with hardware that allows for multiple adjustments. The Relax Snug Fit with its open hip laterals is suited for users needing complex postural support and spinal alignment. Provides 7" of mid Thoracic support.

Quick Reference Guide

RELAX SERIES - SNUG FIT BACK		
STANDARD ADP: SEBWM0027	CUSTOM WIDTH ADP: SEBWM0027 & SEBWM0035	CUSTOM HEIGHT ADP: SEBWM0027 & SEBWM0040
Height: 18" Width: 16" & 18"	Height: 18" Width: All except 16" & 18"	Height: All except 18" Width: 16" & 18"
Available Heights: 16", 18" & 20"		

CUSTOM
ADP: **SEBWM0027** W= **SEBWM0035** & H= **SEBWM0040**


Part Number: Refer below

PART NUMBERS

MDU(W X H) **Snug Laterals**

Fill in width (W) and height (H) to get part number.

E.g. For a W = 15" and H = 17" Back; Part number would be MDU1517




W = Standard shell height 18"
 X = Lateral depth 5" finished
 Y = Space below lateral 4"
 Z = Lateral height 6" above cutout




Rear Back Angle



Inner Foam and Incontinent Cover



Reposition pin hardware back by 1" on the shell.
 Optional 1" longer pins to mount a smaller back on a wider chair



Accommodates Pelvic rotation and will allow for full back support with 25° of posterior and 20° of anterior adjustability.

KEY FEATURES

- Provides complex posterial support with open hip laterals
- Comfortable High Resilience Foam
- 3 lbs. Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Adjustable/Removable Lumbar Support
- Easy to Install and Remove
- Height Adjustable
- TEK4 Waerproof material at the bottom of back for easy clean up

OVERLAYS

* Foam Pad	* Visco	* Gel
ADP: SEBND0015	ADP: SEBND0015	ADP: SEBND0020

* Up-charge items

SPECIFICATIONS	RELAX SERIES - SNUG FIT BACK
ADP	SEBWM0027 W= SEBWM0035 & H= SEBWM0040
STANDARD WIDTH RANGE	16" - 18"
CUSTOM WIDTH RANGE	17" or less and 19" or more
STANDARD HEIGHT RANGE	18"
CUSTOM HEIGHT RANGE	17" or less and 19" or more
AVAILABLE HARDWARE STYLES	One style which is universal to all our backs
ANGLE ADJUSTABLE	25° anterior adjustability, 20° of posterior adjustability
WIDTH ADJUSTABLE	2" of adjustability through the pins on the back
FITS TUBING SIZE	7/8" tubing or 1"
HEADREST MOUNT	Backs pre-drilled for headrest supports and lateral supports
FINISHED HEIGHT	Back Shell Size +2"