

SNUG FORM FIT BACK



RELAX
SERIES



The Relax Form Fit Back with *Snug* laterals is an aluminum light weight formed back with hardware that allows for multiple adjustments. The Relax Form Fit with its wide range of laterals is suited for complex seating and positioning needs. Cubes are removable to custom fit for spinal anomalies.

Quick Reference Guide

RELAX SERIES - SNUG FORM BACK		
STANDARD	CUSTOM WIDTH	CUSTOM HEIGHT
ADP: SEBWM0026	ADP: SEBWM0026 & SEBWM0035	ADP: SEBWM0026 & SEBWM0040
Height: 18" Width: 16" & 18"	Height: 18" Width: All except 16" & 18"	Height: All except 18" Width: 16" & 18"
Available Heights: 16", 18" & 20"		

CUSTOM
ADP: **SEBWM0026** W= **SEBWM0035** & H= **SEBWM0040**

Part Number: Refer below

PART NUMBERS

SHUM (W X H) **Snug Laterals**

Fill in width (W) and height (H) to get part number.


E.g. For a W = 15" and H = 17" Back; Part number would be SHUM1517

Poly Foam Cubes Part Numbers:

QUB6x6	Poly Foam Cubes 6x6
QUB12x12	Poly Foam Cubes 12x12


OVERLAYS

* Foam Pad * Visco * Gel



ADP: SEBND0015 SEBND0015 SEBND0020

* Up-charge items



W = Standard shell height 18"
 X = Lateral depth 5" finished
 Y = Space below lateral 4"
 Z = Lateral height 6" above cutout

CUSTOMIZE ANYTHING


Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

KEY FEATURES

- Snug fit back laterals with removable cubes for custom fit
- Comfortable High Resilience Foam
- 3 lbs. Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Easy to Install and Remove
- Height Adjustable
- TEK4 Waterproof material at the bottom of back for easy clean up



Removable Cube Foam



Accommodates Pelvic rotation and will allow for full back support with 25° of posterior and 20° of anterior adjustability.

SPECIFICATIONS	RELAX SERIES - FORM FIT BACK
ADP	SEBWM0026 W= SEBWM0035 & H= SEBWM0040
STANDARD WIDTH RANGE	16" - 18"
CUSTOM WIDTH RANGE	15" or less and 19" or more
STANDARD HEIGHT RANGE	18"
CUSTOM HEIGHT RANGE	17" or less and 19" or more
AVAILABLE HARDWARE STYLES	One style which is universal to all our backs
ANGLE ADJUSTABLE	25° anterior adjustability, 20° of posterior adjustability
WIDTH ADJUSTABLE	2" of adjustability through the pins on the back
FITS TUBING SIZE	7/8" tubing or 1"
HEADREST MOUNT	Backs pre-drilled for headrest supports and lateral supports
FINISHED HEIGHT	Back Shell Size +2"