

# SOFT CUSHION



# RELAX SERIES



Waterproof zipper seals in the cushion and prevents fluids from entering inside of the cover.



Soft Cushion Layered

The **Relax Soft cushion** is formed to create excellent comfort through a one piece open cell foam. The base is a firm ultra cell polyurethane foam that provides a solid foundation for your cushion. This cushion provides good support and pressure right through the pelvic and thigh areas.

### Quick Reference Guide

#### RELAX SERIES - SOFT CUSHION

REGULAR SIZES ADP: **SESWM0050**

Width = 14" - 20"

Part Number: Refer below

Length = 14" - 20"

#### PART NUMBERS

Part Number: RS(W X L)



Fill in width (W) and length (L) to get part number.

E.g. For a W = 15 and L = 17 Cushion; Part number would be RS1517

#### COVER OPTIONS

Part Number	Description	ADP Codes
RS-COV	Standard Cover	<b>SESND1050</b>

Pressure Reading Scale in mm-Hg



Pressure Map	Relax Soft
Duration	1 hr
Weight	125 lbs./57 Kg
Average Pressure	30 mm-Hg

Pressure Map	Relax Soft
Duration	1 hr
Weight	200 lbs./57 Kg
Average Pressure	30 mm-Hg

**CUSTOMIZE ANYTHING**

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

**KEY FEATURES**

- Superior pressure reduction for users
- Light weight cushion 2.3 lbs.
- Excellent contouring for positioning
- Water fall front for thigh support and comfort
- Adductor for comfort and proper thigh positioning
- Incontinent inner cover
- Easy clean, incontinent, four way stretch outer cover
- Waterproof zipper on standard cushions

**OVERLAYS**



ADP: SEBND0015      SEBND0020

\* Up-charge items

SPECIFICATIONS	RELAX SERIES - SOFT CUSHION
ADP	<b>SESWM0050</b>
PRODUCT WEIGHT	2.3 lbs. (16" - 20")
WIDTH	14" - 20"
DEPTH	14" - 20"
HEIGHT	2.5"
WEIGHT CAPACITY	250 lbs.