

TENSION ADJUSTABLE BACK



RELAX
SERIES

Our **Tension Back** Upholstery provides trunk support and positioning needs. Additional overlays can be added such as Poly Foam, Visco Foam or Gel.

Quick Reference Guide			
Accessories - Tension Back Upholstery			
Part Number	Description	Dia	Product ADP Code
TAB	Tension Back Upholstery Nylon 18" or Less	1	WAMD
TABC	Tension Back Upholstery Custom Nylon 19" or Greater	1	WAME
TABD	Tension Back Upholstery Darlex 18" or Less	6	WAMD
TABCD	Tension Back Upholstery Custom Darlex 19" or Greater	6	WAME
TABE	Tension Back Upholstery Egg Crate 18" or Less	5	WAMD & SEBND0015
TABCE	Tension Back Upholstery Custom Egg Crate 19" or Greater	5	WAME & SEBND0015
TABG	Tension Back Upholstery Gel 18" or Less	4	WAMD & SEBND0020
TABCG	Tension Back Upholstery Custom Gel 19" or Greater	4	WAME & SEBND0020
PRG-P	Additional Pressure Relief Gel (6)	4	SEBND0020
PRF-P	Additional Pressure Relief Foam (6)	2,3	SEBND0015

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.



Rear Back Angle



Rear straps can be adjusted to the curve of the spine

OVERLAYS




Diagram 2, 3, 4, 5 and 6 * Up-charge items

KEY FEATURES

- Comfortable High Resilience Foam
- Easy to Install and Remove
- Height Adjustable
- Adjustable Straps at the Rear

SPECIFICATIONS	RELAX SERIES - TENSION BACK
ADP	WAMD & WAME
STANDARD WIDTH RANGE	16" - 18"
CUSTOM WIDTH RANGE	15" or less and 19" or more
STANDARD HEIGHT RANGE	18"
CUSTOM HEIGHT RANGE	17" or less and 19" or more