

**VGP VOYAGER Wheelchair (Manual Wheelchair)**



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_

| FRAME COLOUR             |         | Select one below   | 1    |
|--------------------------|---------|--------------------|------|
| <input type="checkbox"/> | VOYAGER | PLUS FRAME - Black | 1450 |

| SEAT WIDTH               |           | Select one below | 2   |
|--------------------------|-----------|------------------|-----|
| <input type="checkbox"/> | SW14"-20" | Specify _____    | N/C |

SW14"-20" Light Weight Aluminum Frame

|                          |           |               |     |
|--------------------------|-----------|---------------|-----|
| <input type="checkbox"/> | SW21"-28" | Specify _____ | N/C |
|--------------------------|-----------|---------------|-----|

SW21"-28" Must Include Plus Package 1, 2 Or 3

| SEAT DEPTH |  | Select one below | 3 |
|------------|--|------------------|---|
|------------|--|------------------|---|

|                          |           |                               |     |
|--------------------------|-----------|-------------------------------|-----|
| <input type="checkbox"/> | SDF18"    | (By Frame)                    | N/C |
| <input type="checkbox"/> | SD16"-18" | Specify _____ (By Upholstery) | N/C |
| <input type="checkbox"/> | SDF20"    | (By Frame)                    | 350 |
| <input type="checkbox"/> | SD18"-20" | Specify _____ (By Upholstery) | N/C |
| <input type="checkbox"/> | SDF22"    | (By Frame)                    | 350 |
| <input type="checkbox"/> | SD20"-22" | Specify _____ (By Upholstery) | N/C |

| BACK |  | Select one below | 4 |
|------|--|------------------|---|
|------|--|------------------|---|

|                          |           |                       |     |
|--------------------------|-----------|-----------------------|-----|
| <input type="checkbox"/> | BH15"-20" | Specify _____         | N/C |
| <input type="checkbox"/> | BH21"+    | Specify _____         | 120 |
| <input type="checkbox"/> | AB        | Angle Adj. Back       | N/C |
| <input type="checkbox"/> | ENDO*     | Endomorph Canes W ( ) | N/C |

Specify Width

Endomorph Canes require TA or TAS Arm Rest

| STROLLER HANDLE |  | Select one below | 5 |
|-----------------|--|------------------|---|
|-----------------|--|------------------|---|

|                          |     |                                  |     |
|--------------------------|-----|----------------------------------|-----|
| <input type="checkbox"/> | SHS | Removable Stroller Handle        | N/C |
| <input type="checkbox"/> | SHR | Removable Raised Stroller Handle | N/C |

| FOOT REST |  | Select one below | 6 |
|-----------|--|------------------|---|
|-----------|--|------------------|---|

|                          |    |                           |     |
|--------------------------|----|---------------------------|-----|
| <input type="checkbox"/> | F6 | 60 Degree Hanger          | N/C |
| <input type="checkbox"/> | F7 | 70 Degree Hanger          | N/C |
| <input type="checkbox"/> | F8 | 80 Degree Hanger          | N/C |
| <input type="checkbox"/> | EA | Elevating Legrest (A2)    | 404 |
| <input type="checkbox"/> | AL | Articulating Legrest (A2) | 450 |

| FOOT REST AND FOOT PLATES |  | Select one below | 7 |
|---------------------------|--|------------------|---|
|---------------------------|--|------------------|---|

|                          |     |                                  |     |
|--------------------------|-----|----------------------------------|-----|
| <input type="checkbox"/> | PPL | Plastic                          | N/C |
| <input type="checkbox"/> | PAF | Angle Adjustable - Extra Large   | 150 |
| <input type="checkbox"/> | PXL | Extra Large - Alum (F6A,F7A,F8A) | 150 |
| <input type="checkbox"/> | 1FB | One Piece Foot Board             | 80  |

| ARM REST |  | Select one below | 8 |
|----------|--|------------------|---|
|----------|--|------------------|---|

|                          |       |                                    |     |
|--------------------------|-------|------------------------------------|-----|
| <input type="checkbox"/> | FA    | Flip Back Ht. Adj./ Removable      | 240 |
| <input type="checkbox"/> | TA*   | T-Arm Ht. Adj./ Removable          | 240 |
| <input type="checkbox"/> | FAS   | Flip Back (Short) Ht. Adj./ Remove | 240 |
| <input type="checkbox"/> | TAS*  | T-Arm (Short) Ht. Adj./ Removable  | 240 |
| <input type="checkbox"/> | GWAPS | Gel Waterfall Arm Pads (A1)        | 150 |

14" - 20" Wide Light Weight Aluminum Frame Only

Please Note : A1 = Clients Responsibility

Please Note : A2 = Clinical Justification

**SPECIFICATIONS**

- LIGHT WEIGHT ALUMINUM FRAME AND CROSS BRACES
- STANDARD 18" BACK HEIGHT
- 20", 22", OR 24" NYLON MOULDED REAR WHEELS
- EXTRUDED ALUMINUM FORKS
- CHROME HAND RIMS
- 5", 6", OR 8" CASTERS
- BLACK NYLON UPHOLSTERY
- POWDER COATED FRAME
- FLIP BACK HEIGHT ADJUSTABLE REMOVABLE ARMS
- DETACHABLE SWING AWAY FOOTREST WITH HEEL LOOPS
- PRECISION BEARINGS
- MOULDED PLASTIC FOOTPLATES
- SEAT WIDTH 14" TO 18"
- BLACK PADDED ARM PADS
- DOUBLE CROSS BRACE

**SEAT TO FLOOR SPECIFICATIONS**

|                          | STF        | R.WHEEL | CASTER | FORK     |
|--------------------------|------------|---------|--------|----------|
| <input type="checkbox"/> | SL 12.25   | 20"     | 5"     | XL SHORT |
| <input type="checkbox"/> | SL 13.00"  | 20"     | 5"     | XL SHORT |
| <input type="checkbox"/> | LL 13.75"  | 20"     | 5"     | XL SHORT |
| <input type="checkbox"/> | L 14.50"   | 20"     | 6"     | SHORT    |
| <input type="checkbox"/> | L 14.50"   | 22"     | 5"     | SHORT    |
| <input type="checkbox"/> | L14.50"    | 22"     | 6"     | SHORT    |
| <input type="checkbox"/> | LH 15.25"  | 20"     | 5"     | SHT/LONG |
| <input type="checkbox"/> | LH 15.25"  | 20"     | 6"     | SHT/LONG |
| <input type="checkbox"/> | LH 15.25"  | 22"     | 5"     | SHT/LONG |
| <input type="checkbox"/> | LH 15.25"  | 22"     | 6"     | SHT/LONG |
| <input type="checkbox"/> | SH 16.00"  | 22"     | 6"     | LONG     |
| <input type="checkbox"/> | MH 16.75   | 22"     | 6"     | LONG     |
| <input type="checkbox"/> | H 17.50"   | 22"     | 6"     | LONG     |
| <input type="checkbox"/> | H 17.50"   | 22"     | 8"     | LONG     |
| <input type="checkbox"/> | H 17.50"   | 22"     | 6"     | LONG     |
| <input type="checkbox"/> | H 17.50"   | 24"     | 8"     | LONG     |
| <input type="checkbox"/> | TH18.25"   | 24"     | 8"     | LONG     |
| <input type="checkbox"/> | TT 19.00"* | 24"     | 8"     | LONG     |

\*TT 19.00' requires taller caster journal

| REAR WHEELS              |      |                                 | Select one below | 9 |
|--------------------------|------|---------------------------------|------------------|---|
| <input type="checkbox"/> | RWT0 | 20" MAG Solid Treaded Tire      | N/C              |   |
| <input type="checkbox"/> | RWT2 | 22" MAG Solid Treaded Tire      | N/C              |   |
| <input type="checkbox"/> | RWT4 | 24" MAG Solid Treaded Tire      | N/C              |   |
| <input type="checkbox"/> | RWP2 | 22" Pyramid Tire                | 80               |   |
| <input type="checkbox"/> | RWP4 | 24" Pyramid Tire                | 80               |   |
| <input type="checkbox"/> | RWA2 | 22" Air Tire                    | 120              |   |
| <input type="checkbox"/> | RWA4 | 24" Air Tire                    | 120              |   |
| <input type="checkbox"/> | AQR  | Axles – Quick Release (349 Lbs) | N/C              |   |
| <input type="checkbox"/> | PCHR | Plastic Coated Hand Rims        | N/C              |   |

14" - 20" Wide Light Weight Aluminium Frame Only

| CASTER                   |      |                       | Select one below | 10 |
|--------------------------|------|-----------------------|------------------|----|
| <input type="checkbox"/> | FC1A | 5" X 1.25" Caster     | N/C              |    |
| <input type="checkbox"/> | FC1B | 6" X 1.25" Caster     | N/C              |    |
| <input type="checkbox"/> | FC1C | 8" X 1.25" Caster     | N/C              |    |
| <input type="checkbox"/> | FC2A | 5" X 1.75" Caster     | 80               |    |
| <input type="checkbox"/> | FC2B | 6" X 1.75" Caster     | 80               |    |
| <input type="checkbox"/> | FC2C | 8" X 1.75" Caster     | 80               |    |
| <input type="checkbox"/> | FC3A | 5" X 1.25" Air Caster | 80               |    |
| <input type="checkbox"/> | FC3B | 6" X 1.25" Air Caster | 80               |    |
| <input type="checkbox"/> | FC3C | 8" X 1.25" Air Caster | 80               |    |

| WHEEL LOCK               |      |                        | Select one below | 11 |
|--------------------------|------|------------------------|------------------|----|
| <input type="checkbox"/> | ALPL | Aluminium Pull To Lock | N/C              |    |
| <input type="checkbox"/> | ALPH | Aluminium Push To Lock | N/C              |    |
| <input type="checkbox"/> | STPL | Steel Pull To Lock     | N/C              |    |
| <input type="checkbox"/> | STPH | Steel Push To Lock     | N/C              |    |
| <input type="checkbox"/> | BRXA | 6" Brake Extension     | 40               |    |
| <input type="checkbox"/> | BRXB | 9" Brake Extension     | 60               |    |

| SEATING - BACKS          |     |                         | Select one below | 12 |
|--------------------------|-----|-------------------------|------------------|----|
|                          |     |                         | W ( ) x H ( )    |    |
| <input type="checkbox"/> | TAB | Tension Adjustable Back | 180              |    |

| ACCESSORIES              |         |                             | 14    |
|--------------------------|---------|-----------------------------|-------|
| <input type="checkbox"/> | SBA     | Seat Belt (Auto Style)      | 23    |
| <input type="checkbox"/> | SBB     | Seat Belt (Velcro Style)    | 15.75 |
| <input type="checkbox"/> | CNC     | Crutch And Cane Holder (A1) | 60    |
| <input type="checkbox"/> | IVA     | Adjustable I V Pole (A1)    | 130   |
| <input type="checkbox"/> | NLSP    | Nylon Leg Strap Padded (A1) | 50    |
| <input type="checkbox"/> | VGPAMPL | Amputee Support Left        | 45    |
| <input type="checkbox"/> | VGPAMPR | Amputee Support Right       | 45    |
| <input type="checkbox"/> | VGPOXY  | Oxygen Tank Holder          | 130   |
| <input type="checkbox"/> | VGPATR  | Anti Tippers                | 88    |

Please Note : A1 = Clients Responsibility

| PLUS PACKAGES  |  | Select one below  | 15 |
|--|--|---|----|
| <input type="checkbox"/>   | <b>PLUS PACKAGE 1: AADL CODE: W101</b> | <b>CLIENT 250 TO 349 LBS - \$880.00</b>   |    |
|  | INCLUDES THE FOLLOWING:                | STEEL FRAME, DOUBLE CROSS BRACES WITH HEAVY DUTY BACK AND SEAT UPHOLSTERY       |    |
|  |  | SEAT WIDTH INCREASE FROM 16" TO 24"   |    |
|  |  | SEAT DEPTH INCREASE FROM 16" TO 22"   |    |
|  |  | 5,6 OR 8 X 1 1/4" HD CASTERS  |    |
|  |  | *** WELDED CASTER JOURNAL ***   |    |
| <input type="checkbox"/>   | <b>PLUS PACKAGE 2: AADL CODE: W102</b> | <b>CLIENT 350 TO 449 LBS - \$1140.00</b>  |    |
|  | INCLUDES THE FOLLOWING:                | STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY |    |
|  |  | SEAT WIDTH INCREASE FROM 18" TO 26"   |    |
|  |  | SEAT DEPTH INCREASE FROM 16" TO 22"   |    |
|  |  | 5,6 OR 8 X 1 1/4" HD CASTERS  |    |
|  |  | HD STEEL FORKS  |    |
|  |  | *** WELDED CASTER JOURNAL ***   |    |
| <input type="checkbox"/>   | <b>PLUS PACKAGE 3: AADL CODE: W103</b> | <b>CLIENT 450 TO 600 LBS - \$1300.00</b>  |    |
|  | INCLUDES THE FOLLOWING:                | STEEL FRAME, DOUBLE CROSS BRACES WITH EXTRA HEAVY DUTY BACK AND SEAT UPHOLSTERY |    |
|  |  | SEAT WIDTH INCREASE FROM 20" TO 28"   |    |
|  |  | SEAT DEPTH INCREASE FROM 16" TO 22"   |    |
|  |  | 5,6 OR 8 X 1 1/4" HD CASTERS  |    |
|  |  | HD STEEL FORKS  |    |
|  |  | *** WELDED CASTER JOURNAL ***   |    |
| <input type="checkbox"/>   | FIXED REINFORCED ONE PIECE BACK POST   | \$260.00  |    |
|  | (SPECIFY ANGLE _____)                  |   |    |
| <input type="checkbox"/>   | HEAVY DUTY TENSION ADJ. BACK UPH.      | \$220.00  |    |
| <b>Please Note : Reference Page 12 of the AADL Wheelchair APL for options and catalogue numbers for A wheelchairs.</b> |  |   |    |
| <b>Notes:</b>  |  |   |    |