

**EXT Extreme Wheelchair (Tilt Wheelchair)**



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_

FRAME COLOUR	Select one below	1
<input type="checkbox"/> EXTB	Dynamic Frame - Black	3899
<input type="checkbox"/> EXTC	Dynamic Frame - Colour : _____	4299

CLIENT WEIGHT	Select one below	2
<input type="checkbox"/> EXTWA	Up to 250Lbs	N/C
<input type="checkbox"/> EXTWAMF	Exceeds 250Lbs	730
<input type="checkbox"/> EXTWAMG	Exceeds 350Lbs	992
<input type="checkbox"/> EXTWAMH	Exceeds 400Lbs (Max 600Lbs)	1110




SEAT WIDTH	Select one below	3
<input type="checkbox"/> EXTSW16	Seat Width 16"	N/C
<input type="checkbox"/> EXTSW17	Seat Width 17"	N/C
<input type="checkbox"/> EXTSW18	Seat Width 18"	N/C
<input type="checkbox"/> EXTSW19	Seat Width 19"	399
<input type="checkbox"/> EXTSW20	Seat Width 20"	399
<input type="checkbox"/> EXTSWC	Custom Seat Width: _____	399

SEAT DEPTH	Select one below	4
<input type="checkbox"/> EXTSD16	Seat Depth 16"	N/C
<input type="checkbox"/> EXTSD17	Seat Depth 17"	N/C
<input type="checkbox"/> EXTSD18	Seat Depth 18"	N/C
<input type="checkbox"/> EXTSD19	Seat Depth 19"	265
<input type="checkbox"/> EXTSD20	Seat Depth 20"	265
<input type="checkbox"/> EXTSDC	Custom Seat Depth: _____	265

FOOT REST AND FOOT PLATES	Select one below	5
<input type="checkbox"/> EXTF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	404
<input type="checkbox"/> EXTAL	Articulating Leg Rests <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	450
<input type="checkbox"/> EXTEB	Elevating Legrest System	450
<input type="checkbox"/> EXTPPL	Plastic	N/C
<input type="checkbox"/> EXTPXL	Extra Large	150
<input type="checkbox"/> EXTPAF	Angle Adjustable Footplates	240
<input type="checkbox"/> EXTFRCH	Padded Footplates	88
<input type="checkbox"/> EXT1FB	One Piece Foot Board	80
<input type="checkbox"/> EXTNLS	Nylon Leg Strap Padded	60
<input type="checkbox"/> EXTNLS C	Nylon Leg Strap Custom Size ( )	120
<input type="checkbox"/> EXTC SL	Contracture Sling	240
<input type="checkbox"/> EXTOFST	Offset Foot Rest Receivers	N/C

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE): 16" – 18"
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
DYNAMIC TILT ADJUSTMENTS 0 TO 25 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS

SEAT TO FLOOR SPECIFICATIONS				
(LOW SEAT SETTING)				
STF	R. WHEEL	CASTER	FORK	
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	
(STANDARD SEAT SETTING)				
<input type="checkbox"/> MH 16.75"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	8"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	24"	8"	LONG	
<input type="checkbox"/> A 18.25"	24"	8"	LONG	

BACK TYPE		6	SEATING - BACKS		12
<input type="checkbox"/>	EXTAB1	Angle Adjustable Back	N/C	<b>W ( ) x H ( )</b>	
<i>Specify Width</i> <input type="checkbox"/>	EXTENDO	Endomorph Canes	N/C	<input type="checkbox"/>	EXTBK01 <b>RELAX</b> Gentle Fit Back 544
<input type="checkbox"/>	EXTRBK	Recline Kit	560	<input type="checkbox"/>	EXTBK02 <b>RELAX</b> Full Fit Back 608
HEADREST		7	SEATING - BACKS		12
<i>Select one below</i>			<i>Select one below</i>		
<input type="checkbox"/>	EXTMAHFLAT	Multi Axis Complete Flat Pad	460	<input type="checkbox"/>	EXTBK03 <b>RELAX</b> Snug Fit Back 700
<input type="checkbox"/>	EXTMAHF	Multi Axis Complete With Large Pad	460	<input type="checkbox"/>	EXTBK04 <b>RELAX</b> Form Fit Back No Laterals 820
<input type="checkbox"/>	EXTMAHKP	Multi Axis Complete With Small Pad	460	<input type="checkbox"/>	EXTBK06 <b>RELAX</b> Form Fit Back Gentle Laterals 820
 <input type="checkbox"/>	EXTMAHFX	Multi Axis FLEX Complex Complete Large	525	<input type="checkbox"/>	EXTBK07 <b>RELAX</b> Form Fit Back Full Laterals 820
<i>(1/2" Tubing)</i> <input type="checkbox"/>	EXTUPHT218	Std. vertical horizontal assembly (1/2")	256	<input type="checkbox"/>	EXTBK05 <b>RELAX</b> Form Fit Back Snug Laterals 820
<i>(3/4" Tubing)</i> <input type="checkbox"/>	EXTUPHT217	Std. vertical horizontal assembly (3/4")	256	<input type="checkbox"/>	EXTBK10 <b>RELAX</b> Spiral Fit Back Gentle Laterals 832
<input type="checkbox"/>	EXTUHA	Universal Hardware Adaptor	100	<input type="checkbox"/>	EXTBK09 <b>RELAX</b> Spiral Fit Back Snug Laterals 832
 <input type="checkbox"/>	EXTMAHGA-P	Multi Axis with Headrest Plate	224	<input type="checkbox"/>	EXTBK11 <b>RELAX</b> Spiral Fit Back Full Laterals 832
ARM REST		8	<input type="checkbox"/>	EXTTAB	<b>RELAX</b> Tension Back Upto 18" 140
<i>Select one below</i>			<input type="checkbox"/>	EXTTAB+	<b>RELAX</b> Tension Back Over 18" 240
<i>(Std 10"-14")</i> <input type="checkbox"/>	EXTTA	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	EXTBGL1 Half Overlay Gel (3 Piece) 348
<i>(Short 7"-11")</i> <input type="checkbox"/>	EXTTAS	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	EXTBGL2 Full Overlay Gel (6 Pieces) 696
<input type="checkbox"/>	EXTGWAPS	Gel Waterfall Arm Pads	180	<input type="checkbox"/>	EXTPEXT Pin Extension (1" Extension) 40
REAR WHEELS		9	<input type="checkbox"/>	EXTBK12	*Custom Width : _____ 176
<i>Select one below</i>			<input type="checkbox"/>	EXTBK13	*Custom Height : _____ 176
<input type="checkbox"/>	EXTRWT0	20" Solid Treaded Tire	N/C	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>	
<input type="checkbox"/>	EXTRWT2	22" Solid Treaded Tire	N/C	SEATING - CUSHION	
<input type="checkbox"/>	EXTRWT4	24" Solid Treaded Tire	N/C	<i>Select one below</i>	
<i>(STF 13-19")</i> <input type="checkbox"/>	EXTRWT12	12" Solid Tire	N/C	<b>W ( ) x D ( )</b>	
<input type="checkbox"/>	EXTRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	EXTC01 <b>RELAX</b> Soft Cushion 480
<input type="checkbox"/>	EXTRWA2	22" Air Tire	120	<input type="checkbox"/>	EXTC02 <b>RELAX</b> Comfy Cushion 624
<input type="checkbox"/>	EXTRWA4	24" Air Tire	120	<input type="checkbox"/>	EXTC03 <b>RELAX</b> Firm Cushion 784
<input type="checkbox"/>	EXTAQR	Axles – Quick Release	120	<input type="checkbox"/>	EXTCGL1 Half Overlay Gel (1 Piece) 116
<input type="checkbox"/>	EXTHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	EXTCGL2 Full Overlay Gel (2 Pieces) 232
CASTER		10	<input type="checkbox"/>	EXTCCL	Contoured Leg Channel (Bolster) Std. 180
<i>Select one below</i>			<input type="checkbox"/>	EXTCCLC	Contoured Leg Channel (Bolster) Custom 320
<input type="checkbox"/>	EXTFC1A	5" X 1.25" Caster	N/C	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>	
<input type="checkbox"/>	EXTFC1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	EXTC00 Extra Cushion Cover 99
<input type="checkbox"/>	EXTFC1C	8" X 1.25" Caster	N/C	ACCESSORIES	
<input type="checkbox"/>	EXTFC2A	5" X 1.75" Caster (Soft Roll)	80	<i>Select one below</i>	
<input type="checkbox"/>	EXTFC2B	6" X 1.75" Caster	80	<input type="checkbox"/>	EXTSBA Seat Belt (Auto Style) 26
<input type="checkbox"/>	EXTFC2C	8" X 1.75" Caster	80	<input type="checkbox"/>	EXTSBAP Seat Belt (Padded) 56
<input type="checkbox"/>	EXTFC4B	6" X 2 Caster	80	<input type="checkbox"/>	EXTSBA4P Seat Belt (4 Point Padded) 130
<input type="checkbox"/>	EXTFC4C	8" X 2 Caster	80	<input type="checkbox"/>	EXTTTRA Full Padded Tray Locking 280
WHEEL LOCK		11	<input type="checkbox"/>	EXTTRB	Full Padded Tray Slide On 280
<i>Select one below</i>			<input type="checkbox"/>	EXTTRC	Full Plastic Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 280
<input type="checkbox"/>	EXTALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	EXTTRD Half Padded Left Tray 240
<input type="checkbox"/>	EXTALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	EXTTRE Half Padded Right Tray 240
<input type="checkbox"/>	EXTEXTH	Steel Push To Lock	N/C	<input type="checkbox"/>	EXTTRF Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	EXTEXTL	Steel Pull To Lock	N/C	<input type="checkbox"/>	EXTTRG Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear 240
<input type="checkbox"/>	EXTBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	EXTTDFB Transport Tie Downs Front 85
<input type="checkbox"/>	EXTBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	EXTTDRB Transport Tie Downs Rear 85
<input type="checkbox"/>	EXTBRXB	9" Brake Extension	130	<input type="checkbox"/>	EXTOXY Oxygen Tank Holder 195
 <input type="checkbox"/>	STPATTBRX	Attendant Brakes	240	<input type="checkbox"/>	EXTIVA Adjustable I V Pole 130
<input type="checkbox"/>				<input type="checkbox"/>	EXTCHPT-L Left Arm Trough Tray 400
<input type="checkbox"/>				<input type="checkbox"/>	EXTCHPT-R Right Arm Trough Tray 400
<input type="checkbox"/>				<input type="checkbox"/>	EXTLATB Lateral Bite out 180
<input type="checkbox"/>				<input type="checkbox"/>	EXTADD8 8" Adductor 164
<input type="checkbox"/>				<input type="checkbox"/>	EXTADD10 10" Adductor 164
<input type="checkbox"/>				<input type="checkbox"/>	EXTADD12 12" Adductor 164