

EXT Extreme Wheelchair - Tilt Wheelchair



INFORMATION

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME TYPE	Select one below	1	ADP
<input type="checkbox"/> EXTWANA Dynamic Tilt Frame		\$1,200	WANA

FRAME COLOUR	Select one below	1	ADP
<input type="checkbox"/> EXTB Frame Color - Black		\$2,699	WM0001041
<input type="checkbox"/> EXTC Frame - Colour : _____		\$2,999	WM0001041

CLIENT WEIGHT	Select one below	2	ADP
<input type="checkbox"/> EXTWA Up to 250Lbs		N/C	
<input type="checkbox"/> EXTWAMF Exceeds 250Lbs		\$730	WAMF
<input type="checkbox"/> EXTWAMG Exceeds 350Lbs		\$992	WAMG
<input type="checkbox"/> EXTWAMH Exceeds 400Lbs (Max 600Lbs)		\$1,110	WAMH

SEAT WIDTH	Select one below	3	ADP
<input type="checkbox"/> EXTSW16 Seat Width 16"		N/C	
<input type="checkbox"/> EXTSW17 Seat Width 17"		N/C	
<input type="checkbox"/> EXTSW18 Seat Width 18"		N/C	
<input type="checkbox"/> EXTSW19 Seat Width 19"		\$399	WAMA
<input type="checkbox"/> EXTSW20 Seat Width 20"		\$399	WAMA
<input type="checkbox"/> EXTSWC Custom Seat Width: _____		\$399	WAMA

SEAT DEPTH	Select one below	4	ADP
<input type="checkbox"/> EXTS16 Seat Depth 16"		N/C	
<input type="checkbox"/> EXTS17 Seat Depth 17"		N/C	
<input type="checkbox"/> EXTS18 Seat Depth 18"		N/C	
<input type="checkbox"/> EXTS19 Seat Depth 19"		\$265	WAMB
<input type="checkbox"/> EXTS20 Seat Depth 20"		\$265	WAMB
<input type="checkbox"/> EXTSDC Custom Seat Depth: _____		\$265	WAMB

FOOT REST AND FOOT PLATES	Select one below	5	ADP
<input type="checkbox"/> EXTF6A 60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/> EXTF7A 70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/> EXTF8A 80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		N/C	
<input type="checkbox"/> EXTEA Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		\$404	WAMK
<input type="checkbox"/> EXTAL Articulating Leg Rests <input type="checkbox"/> Pivot <input type="checkbox"/> Pin		\$450	WAMK
<input type="checkbox"/> EXTEB Elevating Legrest System		\$450	WAMK
<input type="checkbox"/> EXTPPL Plastic Footplates		N/C	
<input type="checkbox"/> EXTPXL Extra Large Footplates		\$150	SE0000322
<input type="checkbox"/> EXTPAF Angle Adjustable Footplates		\$240	WAMJ
<input type="checkbox"/> EXT1FB One Piece Foot Board		\$80	SEFND0015
<input type="checkbox"/> EXT1FBP One Piece Foot Board Padding		\$44	SEFND0035
<input type="checkbox"/> EXTNLS Nylon Leg Strap Padded		\$60	SEFND0010
<input type="checkbox"/> EXTNLSLSC Nylon Leg Strap Custom Size ()		\$120	SEFCF107L
<input type="checkbox"/> EXTCSL Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style		\$240	SEFND0020
<input type="checkbox"/> EXTOFST Offset Foot Rest Receivers		N/C	

BACK TYPE	Select one below	6	ADP
<input type="checkbox"/> EXTAB1 Angle Adjustable Back		N/C	
<input type="checkbox"/> EXTENDO Endomorph Canes W ()		N/C	
<input type="checkbox"/> EXTRBK Recline Kit		\$560	WAMI

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE): 16" – 18"
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
DYNAMIC TILT ADJUSTMENTS 0 TO 25 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CABLES
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/> 12.25"	20"	5"	XL SHORT	
<input type="checkbox"/> 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> 14.50"	20"	6"	SHORT	
<input type="checkbox"/> 14.50"	22"	5"	SHORT	
<input type="checkbox"/> 14.50"	22"	6"	SHORT	
<input type="checkbox"/> 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> 16.00"	22"	6"	LONG	
(HIGH SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/> 16.75"	22"	5"	LONG	
<input type="checkbox"/> 17.50"	22"	5"	LONG	
<input type="checkbox"/> 17.50"	22"	5"	LONG	
<input type="checkbox"/> 17.50"	22"	8"	LONG	
<input type="checkbox"/> 17.50"	24"	8"	LONG	
<input type="checkbox"/> 18.25"	24"	6" or 8"	LONG	

For more Seat To Floor combinations, visit our website

HEADREST		Select one below	7	ADP	SEATING - BACKS		Custom: select back & fill custom	12	ADP
<input type="checkbox"/>	EXTMAHF	Multi-Axis Hardware Complete with Large Pad	\$512	SE0001542	<input type="checkbox"/>	EXTBK01	RELAX Gentle Fit Back	\$544	SEBWM0020
<input type="checkbox"/>	EXTMAHFXL	Multi-Axis Hardware Complete with Flex Large	\$512	SE0001542	<input type="checkbox"/>	EXTBK02	RELAX Full Fit Back	\$608	SEBWM0010
<input type="checkbox"/>	EXTMAHKP	Multi-Axis Hardware Complete with Medium Pad	\$512	SE0001542	<input type="checkbox"/>	EXTBK03	RELAX Snug Fit Back	\$700	SEBWM0027
<input type="checkbox"/>	EXTMAHFXS	Multi-Axis Hardware Complete with Flex Medium	\$512	SE0001542	<input type="checkbox"/>	EXTBK04	RELAX Form Fit Back No Laterals	\$820	SEBWM0026
<input type="checkbox"/>	EXTMAHFT	Multi-Axis Hardware Complete with Small Flat	\$512	SE0001542	<input type="checkbox"/>	EXTBK06	RELAX Form Fit Back Gentle Laterals	\$820	SEBWM0026
<input type="checkbox"/>	EXTUPHT219	Contour Headrest Complete with Swivel Hardware	\$433	SEHND0005 SEHCF050L	<input type="checkbox"/>	EXTBK07	RELAX Form Fit Back Full Laterals	\$820	SEBWM0026
<input type="checkbox"/>	EXTUPHT218	Contour Headrest with Std. Assembly (1/2" Tube)	\$292	SEHND0020	<input type="checkbox"/>	EXTBK05	RELAX Form Fit Back Snug Laterals	\$820	SEBWM0026
<input type="checkbox"/>	EXTUPHT217	Contour Headrest with Std. Assembly (3/4" Tube)	\$221	SEHWM0005	<input type="checkbox"/>	EXTBK10	RELAX Spiral Fit Back Gentle Laterals	\$832	SEBWM0030
<input type="checkbox"/>	EXTUHRA	Universal Hardware Adaptor (Adaptor only)	\$57	SEHND0020	<input type="checkbox"/>	EXTBK09	RELAX Spiral Fit Back Snug Laterals	\$832	SEBWM0030
ARM REST		Select one below	8	ADP	<input type="checkbox"/>	EXTBK11	RELAX Spiral Fit Back Full Laterals	\$832	SEBWM0030
(Ht. 10"-14")	<input type="checkbox"/>	EXTTA	T-Arm Ht. Adj./ Removable	\$240	<input type="checkbox"/>	EXTTAB	RELAX Tension Back Upto 18"	\$140	WAMD
(Ht. 10"-14")	<input type="checkbox"/>	EXTTAS	T-Arm Ht. Adj./ Removable Short	\$240	<input type="checkbox"/>	EXTTAB+	RELAX Tension Back Over 18"	\$240	WAME
<input type="checkbox"/>	EXTGWAFP	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030	<input type="checkbox"/>	EXTBK12	* Custom Width : _____	\$176	SEBWM0035
* Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required									
REAR WHEEL		Select one below	9	ADP	<input type="checkbox"/>	EXTBK13	* Custom Height : _____	\$176	SEBWM0040
<input type="checkbox"/>	EXTRWT0	20" Solid Treaded Tire	N/C		*Note: Please refer to the Seating Booklet for custom sizes				
<input type="checkbox"/>	EXTRWT2	22" Solid Treaded Tire	N/C		<input type="checkbox"/>	EXTBGL1	Gel Full Overlay	\$348	SEBND0020 X 3
<input type="checkbox"/>	EXTRWT4	24" Solid Treaded Tire	N/C		<input type="checkbox"/>	EXTBGL2	Gel Center Vertebrae	\$116	SEBND0020
(STF 13-19")	<input type="checkbox"/>	EXTRWT12	* 12" Solid Tire	N/C	<input type="checkbox"/>	EXTBGL3	Gel Left Lateral	\$116	SEBND0020
<input type="checkbox"/>	EXTRWT16	* 16" Solid Tire	N/C		<input type="checkbox"/>	EXTBGL4	Gel Right Lateral	\$116	SEBND0020
<input type="checkbox"/>	EXTRWA2	22" Air Tire	\$120		SEATING - CUSHIONS Select one below 13 ADP				
<input type="checkbox"/>	EXTRWA4	24" Air Tire	\$120		<input type="checkbox"/>	EXTC01	RELAX Soft Cushion	\$480	SESWM0050
<input type="checkbox"/>	EXTAQR	Axles – Quick Release	\$120	WAMU	<input type="checkbox"/>	EXTC02	RELAX Comfy Cushion	\$624	SESWM0030
<input type="checkbox"/>	EXTHRT	Narrow Handrim Config. (1/2" less)	N/C		<input type="checkbox"/>	EXTC03	RELAX Firm Cushion	\$784	SESWM0040
*Note: Quick Release is not available for wheels smaller than 20"									
CASTER		Select one below	10	ADP	<input type="checkbox"/>	EXTC04	RELAX AIR FLOW Cushion	\$587	SE0001669
<input type="checkbox"/>	EXTFC1A	5" X 1.25" Caster	N/C		<input type="checkbox"/>	EXTCGL2	Full Overlay Gel	\$258	SESND1045
<input type="checkbox"/>	EXTFC1B	6" X 1.25" Caster	N/C		<input type="checkbox"/>	EXTCCL	Contoured Leg Channel (Bolster) Std.	\$180	SEMCF203L
<input type="checkbox"/>	EXTFC1C	8" X 1.25" Caster	N/C		<input type="checkbox"/>	EXTCCLA	Contoured Leg Channel (Bolster) with ABS	\$260	SESDN1035 SESDN1015
<input type="checkbox"/>	EXTFC2A	5" X 1.75" Caster (Soft Roll)	\$80		<input type="checkbox"/>	EXTCCLC	* Contoured Leg Channel (Bolster) Custom	\$320	SEMCF209L
<input type="checkbox"/>	EXTFC2B	6" X 1.75" Caster	\$80		*Note: Please refer to the Seating Booklet for custom sizes				
<input type="checkbox"/>	EXTFC2C	8" X 1.75" Caster	\$80		<input type="checkbox"/>	EXTC00	Extra Cushion Cover	\$99	SESND1055
<input type="checkbox"/>	EXTFC4B	6" X 2 Caster	\$80		ACCESSORIES		14	ADP	
<input type="checkbox"/>	EXTFC4C	8" X 2 Caster	\$80		<input type="checkbox"/>	EXTSBA	Seat Belt (Auto Style)	\$26	
WHEEL LOCK		Select one below	11	ADP	<input type="checkbox"/>	EXTSBAP	Seat Belt (Padded)	\$56	SERND0030
<input type="checkbox"/>	EXTALPH	Aluminium Push To Lock	N/C		<input type="checkbox"/>	EXTSBA4P	Seat Belt (4 Point Padded)	\$130	SERND0020
<input type="checkbox"/>	EXTALPL	Aluminium Pull To Lock	N/C		<input type="checkbox"/>	EXTTRA	Full Padded Tray Locking	\$280	SETND0005 SETND0035 SETND0005 SETND0035
<input type="checkbox"/>	EXTSTPH	Steel Push To Lock	N/C		<input type="checkbox"/>	EXTTRB	Full Padded Tray Slide On	\$280	SETND0010
<input type="checkbox"/>	EXTSTPL	Steel Pull To Lock	N/C		<input type="checkbox"/>	EXTTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	EXTBRXC	4" Brake Extension	N/C		<input type="checkbox"/>	EXTTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010
<input type="checkbox"/>	EXTBRXA	6" Brake Extension	N/C		<input type="checkbox"/>	EXTTRD	Half Padded Left Tray	\$240	SETND0015 SETND0035 SETND0015 SETND0035
<input type="checkbox"/>	EXTBRXB	9" Brake Extension	\$130		<input type="checkbox"/>	EXTTRE	Half Padded Right Tray	\$240	SETND0015
<input type="checkbox"/>	EXTATTBRX	Attendant Brakes	\$240		<input type="checkbox"/>	EXTTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>	EXTULB	Unilateral Hand Brake <input type="checkbox"/> Left <input type="checkbox"/> Right	\$120		<input type="checkbox"/>	EXTTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015
<input type="checkbox"/>					<input type="checkbox"/>	EXTTDFB	Transport Tie Downs Front	\$85	
<input type="checkbox"/>					<input type="checkbox"/>	EXTTDRB	Transport Tie Downs Rear	\$85	
<input type="checkbox"/>					<input type="checkbox"/>	EXTOXY	Oxygen Tank Holder	\$195	WEPX
<input type="checkbox"/>					<input type="checkbox"/>	EXTIVA	Adjustable I V Pole	\$130	
<input type="checkbox"/>					<input type="checkbox"/>	EXTCHPT-L	Left Arm Trough Tray	\$400	SETND0015 SETND0035 SETND0030
<input type="checkbox"/>					<input type="checkbox"/>	EXTCHPT-R	Right Arm Trough Tray	\$400	SELDND0015
<input type="checkbox"/>					<input type="checkbox"/>	EXTLATB	Lateral Bite out	\$180	SELDND0015
<input type="checkbox"/>					<input type="checkbox"/>	EXTADD8	8" Adductor	\$164	SELDND0005
<input type="checkbox"/>					<input type="checkbox"/>	EXTADD10	10" Adductor	\$164	SELDND0005
<input type="checkbox"/>					<input type="checkbox"/>	EXTADD12	12" Adductor	\$164	SELDND0005
<input type="checkbox"/>					<input type="checkbox"/>	EXTFRCH	Padded Footplates	\$88	SEMCF803L