

NVT NVT Wheelchair (Tilt Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

- NVTB Tilt Frame - Black 3299
- NVTS Tilt Frame - Silver 3299
- NVTC Tilt Frame - Custom : _____ 3699

CLIENT WEIGHT *Select one below* **2**

- NVTWA Up to 250lbs N/C
- NVTWAMF Exceeds 250lbs 730
- NVTWAMG Exceeds 350lbs 992
- NVTWAMH Exceeds 400lbs (Max 600lbs) 1110

SEAT WIDTH *Select one below* **3**

- NVTSW16 Seat Width 16" N/C
- NVTSW17 Seat Width 17" N/C
- NVTSW18 Seat Width 18" N/C
- NVTSW19 Seat Width 19" 399
- NVTSW20 Seat Width 20" 399
- NVTSWC Custom Seat Width: _____ 399

SEAT DEPTH *Select one below* **4**

- NVTSD16 Seat Depth 16" N/C
- NVTSD17 Seat Depth 17" N/C
- NVTSD18 Seat Depth 18" N/C
- NVTSD19 Seat Depth 19" 265
- NVTSD20 Seat Depth 20" 265
- NVTSDC Custom Seat Depth: _____ 265

FOOT REST AND FOOT PLATES *Select one below* **5**

- NVTOFST Offset Foot Rest N/C
- NVTF6A 60 Degree Hanger Pivot Pin N/C
- NVTF7A 70 Degree Hanger Pivot Pin N/C
- NVTF8A 80 Degree Hanger Pivot Pin N/C
- NVTEA Elevating Legrest Pivot Pin 404
- NVTAL Articulating L/R Pivot Pin 450
- NVTEB Elevating Legrest System N/C
- NVTPLP Plastic N/C
- NVTPXL Extra Large 150
- NVTPAF Angle Adjustable Footplates 240
- NVTFRCH Padded Footplates 88
- NVT1FB One Piece Foot Board 80
- NVTNLS Nylon Leg Strap Padded 60
- NVTNLSL Nylon Leg Strap Custom **Size ()** 120
- NVTCSL Contracture Sling 240

STANDARD NO CHARGE SPECIFICATIONS

SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

Tilt Options **6**

- NVTINTS Independant Power Tilt System 2000

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)

	STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	

(STANDARD SEAT SETTING)

<input type="checkbox"/> MH 16.75"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	22"	8"	LONG
<input type="checkbox"/> H 17.50"	22"	6"	LONG
<input type="checkbox"/> H 17.50"	24"	8"	LONG
<input type="checkbox"/> A 18.25"	24"	8"	LONG

Notes:

BACK TYPE			7	SEATING - BACKS			13
<input type="checkbox"/>	NVTAB1	Angle Adjustable Back	N/C	W () x H ()			
<i>Specify Width</i> <input type="checkbox"/>	NVTENDO	Endomorph Canes	N/C	<input type="checkbox"/>	NVTBK01	RELAX Gentle Fit Back	544
<input type="checkbox"/>	NVTRBK	Recline Kit	N/C	<input type="checkbox"/>	NVTBK02	RELAX Full Fit Back	608
HEADREST			8	<input type="checkbox"/>	NVTBK03	RELAX Snug Fit Back	700
<i>Select one below</i>				<input type="checkbox"/>	NVTBK04	RELAX Form Fit Back No Laterals	820
<input type="checkbox"/>	NVMAHF	Multi-Axis Hardware Complete with Large Pad	512	<input type="checkbox"/>	NVTBK06	RELAX Form Fit Back Gentle Laterals	820
<input type="checkbox"/>	NVMAHFXL	Multi-Axis Hardware Complete with Flex Large	512	<input type="checkbox"/>	NVTBK07	RELAX Form Fit Back Full Laterals	820
<input type="checkbox"/>	NVMAHKP	Multi-Axis Hardware Complete with Medium Pad	512	<input type="checkbox"/>	NVTBK05	RELAX Form Fit Back Snug Laterals	820
<input type="checkbox"/>	NVMAHFXS	Multi-Axis Hardware Complete with Flex Medium	512	<input type="checkbox"/>	NVTBK10	RELAX Spiral Fit Back Gentle Laterals	832
<input type="checkbox"/>	NVMAHFT	Multi-Axis Hardware Complete with Small Flat	512	<input type="checkbox"/>	NVTBK09	RELAX Spiral Fit Back Snug Laterals	832
<i>(1/2" Tubing)</i> <input type="checkbox"/>	NVUPHT219	Contour Headrest Complete with Swivel Hardware	433	<input type="checkbox"/>	NVTBK11	RELAX Spiral Fit Back Full Laterals	832
<i>(3/4" Tubing)</i> <input type="checkbox"/>	NVUPHT218	Contour Headrest with Std. Assembly (1/2")	292	<input type="checkbox"/>	NVTTAB	RELAX Tension Back Upto 18"	140
<input type="checkbox"/>	NVUPHT217	Contour Headrest with Std. Assembly (3/4")	292	<input type="checkbox"/>	NVTTAB+	RELAX Tension Back Over 18"	240
<input type="checkbox"/>	NVUHA	Universal Hardware Adaptor (Adaptor only)	57	<input type="checkbox"/>	NVTBGL1	Half Overlay Gel (3 Piece)	348
<input type="checkbox"/>	NVMAHGA-P	Multi Axis with Headrest Plate (No Pad)	255	<input type="checkbox"/>	NVTBGL2	Full Overlay Gel (6 Pieces)	696
ARM REST			9	<input type="checkbox"/>	NVTNPVT	Pin Extension (1" Extension)	40
<i>Select one below</i>				<input type="checkbox"/>	NVTBK12	*Custom Width : _____	+176
<i>(Std 10"-14")</i> <input type="checkbox"/>	NVTTA	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	NVTBK13	*Custom Height : _____	+176
<i>(Std 10"-14")</i> <input type="checkbox"/>	NVTFA	Flip Back Ht. Adj./ Removable	340	*Note: Please refer to the Seating Booklet for custom sizes			
<i>(Short 7"-11")</i> <input type="checkbox"/>	NVTTAS	T-Arm Ht. Adj./ Removable	240	SEATING - CUSHION			14
<i>(Short 7"-11")</i> <input type="checkbox"/>	NVTFAS	Flip Back Ht. Adj./ Removable	340	<i>Select one below</i>			
<input type="checkbox"/>	NVTGWAFP	* Gel Waterfall Arm Pads	180	W () x D ()			
*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required				<input type="checkbox"/>	NVTC01	RELAX Soft Cushion	480
REAR WHEELS			10	<input type="checkbox"/>	NVTC02	RELAX Comfy Cushion	624
<i>Select one below</i>				<input type="checkbox"/>	NVTC03	RELAX Firm Cushion	784
<input type="checkbox"/>	NVTRWTO	20" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTCGL1	Half Overlay Gel (1 Piece)	116
<input type="checkbox"/>	NVTRWT2	22" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTCGL2	Full Overlay Gel (2 Pieces)	232
<input type="checkbox"/>	NVTRWT4	24" Solid Treaded Tire	N/C	<input type="checkbox"/>	NVTCCL	Contoured Leg Channel (Bolster) Std.	180
<i>(STF 13-19")</i> <input type="checkbox"/>	NVTRWT12	12" Solid Tire	N/C	<input type="checkbox"/>	NVCCLA	Contoured Leg Channel (Bolster) with ABS	260
<i>(STF 17-19")</i> <input type="checkbox"/>	NVTRWT12I	12" Solid Tire Inside Mount	120	<input type="checkbox"/>	NVTCCLC	* Contoured Leg Channel (Bolster) Custom	320
<input type="checkbox"/>	NVTRWT16	16" Solid Tire	N/C	*Note: Please refer to the Seating Booklet for custom sizes			
<input type="checkbox"/>	NVTRWA2	22" Air Tire	120	<input type="checkbox"/>	NVTC00	Extra Cushion Cover	99
<input type="checkbox"/>	NVTRWA4	24" Air Tire	120	ACCESSORIES			15
<input type="checkbox"/>	NVTAQR	Axles – Quick Release	120	<input type="checkbox"/>	NVTSBA	Seat Belt (Auto Style)	26
<input type="checkbox"/>	NVTHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	NVTSBAP	Seat Belt (Padded)	56
CASTER			11	<input type="checkbox"/>	NVTSBA4P	Seat Belt (4 Point Padded)	130
<i>Select one below</i>				<input type="checkbox"/>	NVTTA	Full Padded Tray Locking	280
<input type="checkbox"/>	NVTF1A	5" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTTB	Full Padded Tray Slide On	280
<input type="checkbox"/>	NVTF1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTTCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/>	NVTF1C	8" X 1.25" Caster	N/C	<input type="checkbox"/>	NVTTCLC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280
<input type="checkbox"/>	NVTF2A	5" X 1.75" Caster (Soft Roll)	80	<input type="checkbox"/>	NVTTD	Half Padded Left Tray	240
<input type="checkbox"/>	NVTF2B	6" X 1.75" Caster	80	<input type="checkbox"/>	NVTTRE	Half Padded Right Tray	240
<input type="checkbox"/>	NVTF2C	8" X 1.75" Caster	80	<input type="checkbox"/>	NVTTTF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	NVTF4B	6" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	NVTTTG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240
<input type="checkbox"/>	NVTF4C	8" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	NVTTDFB	Transport Tie Downs Front	85
WHEEL LOCK			12	<input type="checkbox"/>	NVTTDRB	Transport Tie Downs Rear	85
<i>Select one below</i>				<input type="checkbox"/>	NVTOXY	Oxygen Tank Holder	195
<input type="checkbox"/>	NVTALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	NVTIVA	Adjustable I V Pole	130
<input type="checkbox"/>	NVTALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	NVTCHPT-L	Left Arm Trough Tray	400
<input type="checkbox"/>	NVTNVTH	Steel Push To Lock	N/C	<input type="checkbox"/>	NVTCHPT-R	Right Arm Trough Tray	400
<input type="checkbox"/>	NVTNVTL	Steel Pull To Lock	N/C	<input type="checkbox"/>	NVTLATB	Lateral Bite out	180
<input type="checkbox"/>	NVTBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	NVTADD8	8" Adductor	164
<input type="checkbox"/>	NVTBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	NVTADD10	10" Adductor	164
<input type="checkbox"/>	NVTBRXB	9" Brake Extension	130	<input type="checkbox"/>	NVTADD12	12" Adductor	164
<input type="checkbox"/>	NVTATTBRX	Attendant Brakes	240				