

**STP Super Tilt Plus Wheelchair - Tilt Wheelchair**



**INFORMATION**

Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



FRAME COLOUR		Select one below	1	ADP
<input type="checkbox"/>	STPB	Tilt Frame - Black	2699	WM0001040
<input type="checkbox"/>	STPS	Tilt Frame - Silver	2699	WM0001040
<input type="checkbox"/>	STPBL	Tilt Frame - Blue	2699	WM0001040
<input type="checkbox"/>	STPBUR	Tilt Frame - Burgundy	2699	WM0001040
<input type="checkbox"/>	STPC	Tilt Frame - Custom : _____	3099	WM0001040

CLIENT WEIGHT		Select one below	2	ADP
<input type="checkbox"/>	STPWA	Up to 250Lbs	N/C	
<input type="checkbox"/>	STPWAMF	Exceeds 250Lbs	\$730	WAMF
<input type="checkbox"/>	STPWAMG	Exceeds 350Lbs	\$992	WAMG
<input type="checkbox"/>	STPWAMH	Exceeds 400Lbs (Max 600Lbs)	\$1,110	WAMH

SEAT WIDTH		Select one below	3	ADP
<input type="checkbox"/>	STPSW16	Seat Width 16"	N/C	
<input type="checkbox"/>	STPSW17	Seat Width 17"	N/C	
<input type="checkbox"/>	STPSW18	Seat Width 18"	N/C	
<input type="checkbox"/>	STPSW19	Seat Width 19"	\$399	WAMA
<input type="checkbox"/>	STPSW20	Seat Width 20"	\$399	WAMA
<input type="checkbox"/>	STPSWC	Custom Seat Width: _____	\$399	WAMA

SEAT DEPTH		Select one below	4	ADP
<input type="checkbox"/>	STPSD16	Seat Depth 16"	N/C	
<input type="checkbox"/>	STPSD17	Seat Depth 17"	N/C	
<input type="checkbox"/>	STPSD18	Seat Depth 18"	N/C	
<input type="checkbox"/>	STPSD19	Seat Depth 19"	\$265	WAMB
<input type="checkbox"/>	STPSD20	Seat Depth 20"	\$265	WAMB
<input type="checkbox"/>	STPSDC	Custom Seat Depth: _____	\$265	WAMB

FOOT REST AND FOOT PLATES		Select one below	5	ADP
<input type="checkbox"/>	STPOFST	Offset Foot Rest	N/C	
<input type="checkbox"/>	STPF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	
<input type="checkbox"/>	STPF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	
<input type="checkbox"/>	STPF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C	
<input type="checkbox"/>	STPEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	\$404	WAMK
<input type="checkbox"/>	STPAL	Articulating L/R <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	\$450	WAMK
<input type="checkbox"/>	STPEB	Elevating Legrest System	\$450	WAMK
<input type="checkbox"/>	STPPPL	Plastic Footplates	N/C	
<input type="checkbox"/>	STPPXL	Extra Large Footplates	\$150	SE0000322
<input type="checkbox"/>	STPPAF	Angle Adjustable Footplates	\$240	WAMJ
<input type="checkbox"/>	STP1FB	One Piece Foot Board	\$80	SEFND0015
<input type="checkbox"/>	STP1FBP	One Piece Foot Board Padding	\$44	SEFND0035
<input type="checkbox"/>	STPNLS	Nylon Leg Strap Padded	\$60	SEFND0010
<input type="checkbox"/>	STPNLSC	Nylon Leg Strap Custom Size ( )	\$120	SEFCF107L
<input type="checkbox"/>	STPCSCL	Contracture Sling <input type="checkbox"/> Standard <input type="checkbox"/> M Style	\$240	SEFND0020

STANDARD NO CHARGE SPECIFICATIONS	
SEAT WIDTH: 14" – 18"	
SEAT DEPTH (ADJUSTABLE)	
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"	
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]	
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]	
WEIGHT CAPACITY: 250 LBS	
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)	
TILT ADJUSTMENTS 0 TO 50 (DEGREES)	
TWO-WAY LOCKING GAS STRUT	
SINGLE RELEASE, DUAL CYLINDERS	
STANDARD LOCKOUT COLLARS	
STANDARD T STYLE ARMS	
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS	

TILT OPTION		Select one below	6	ADP
<input type="checkbox"/>	STPINTS	Independent Power Tilt System	\$2,000	

SEAT TO FLOOR SPECIFICATIONS (± 0.5")				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/>	13.00"	20"	5"	XL SHORT
<input type="checkbox"/>	13.75"	20"	5"	XL SHORT
<input type="checkbox"/>	14.50"	20"	6"	SHORT
<input type="checkbox"/>	14.50"	22"	5"	SHORT
<input type="checkbox"/>	14.50"	22"	6"	SHORT
<input type="checkbox"/>	15.25"	20"	5"	SHT/LONG
<input type="checkbox"/>	15.25"	20"	6"	SHT/LONG
<input type="checkbox"/>	15.25"	22"	5"	SHT/LONG
<input type="checkbox"/>	15.25"	22"	6"	SHT/LONG
<input type="checkbox"/>	16.00"	22"	6"	LONG
(HIGH SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/>	16.75"	22"	6"	LONG
<input type="checkbox"/>	17.50"	22"	6"	LONG
<input type="checkbox"/>	17.50"	22"	6"	LONG
<input type="checkbox"/>	17.50"	24"	6"	LONG
<input type="checkbox"/>	18.25"	24"	8"	LONG

For more Seat To Floor combinations, visit our website

BACK TYPE		Select one below	7	ADP	SEATING -BACKS		Custom: select back & fill custom	13	ADP	
<input type="checkbox"/>	STPAB1	Angle Adjustable Back	N/C		<input type="checkbox"/>	STPBK01	W ( ) x H ( ) <b>RELAX</b> Gentle Fit Back	\$544	SEBWM0020	
Specify Width	<input type="checkbox"/>	STPENDO	Endomorph Canes	W ( )	N/C	<input type="checkbox"/>	STPBK02	<b>RELAX</b> Full Fit Back	\$608	SEBWM0010
<input type="checkbox"/>	STPRBK	Recline Kit	\$560	WAMI	<input type="checkbox"/>	STPBK03	<b>RELAX</b> Snug Fit Back	\$700	SEBWM0027	
HEADREST		Select one below	8	ADP	<input type="checkbox"/>	STPBK04	<b>RELAX</b> Form Fit Back No Laterals	\$820	SEBWM0026	
<input type="checkbox"/>	STPMAHF	Multi-Axis Hardware Complete with Large Pad	\$512	SE0001542	<input type="checkbox"/>	STPBK06	<b>RELAX</b> Form Fit Back Gentle Laterals	\$820	SEBWM0026	
<input type="checkbox"/>	STPMAHFXL	Multi-Axis Hardware Complete with Flex Large	\$512	SE0001542	<input type="checkbox"/>	STPBK07	<b>RELAX</b> Form Fit Back Full Laterals	\$820	SEBWM0026	
<input type="checkbox"/>	STPMAHKP	Multi-Axis Hardware Complete with Medium Pad	\$512	SE0001542	<input type="checkbox"/>	STPBK05	<b>RELAX</b> Form Fit Back Snug Laterals	\$820	SEBWM0026	
<input type="checkbox"/>	STPMAHFXS	Multi-Axis Hardware Complete with Flex Medium	\$512	SE0001542	<input type="checkbox"/>	STPBK10	<b>RELAX</b> Spiral Fit Back Gentle Laterals	\$832	SEBWM0030	
<input type="checkbox"/>	STPMAHFT	Multi-Axis Hardware Complete with Small Flat	\$512	SE0001542	<input type="checkbox"/>	STPBK09	<b>RELAX</b> Spiral Fit Back Snug Laterals	\$832	SEBWM0030	
<input type="checkbox"/>	STPUPHT219	Contour Headrest Complete with Swivel Hardware	\$433	SEHND0005 SEHCF050L	<input type="checkbox"/>	STPBK11	<b>RELAX</b> Spiral Fit Back Full Laterals	\$832	SEBWM0030	
<input type="checkbox"/>	STPUPHT218	Contour Headrest with Std. Assembly (1/2" Tube)	\$292	SEHND0020	<input type="checkbox"/>	STPTAB	<b>RELAX</b> Tension Back Upto 18"	\$140	WAMD	
<input type="checkbox"/>	STPUPHT217	Contour Headrest with Std. Assembly (3/4" Tube)	\$221	SEHWM0005	<input type="checkbox"/>	STPTAB+	<b>RELAX</b> Tension Back Over 18"	\$240	WAME	
<input type="checkbox"/>	STPUHRA	Universal Hardware Adaptor (Adaptor only)	\$57	SEHND0020	<input type="checkbox"/>	STPBK12	* Custom Width : _____	\$176	SEBWM0035	
					<input type="checkbox"/>	STPBK13	* Custom Height : _____	\$176	SEBWM0040	
							<b>*Note: Please refer to the Seating Booklet for custom sizes</b>			
ARM REST		Select one below	9	ADP	<input type="checkbox"/>	STPBGL1	Gel Full Overlay	\$348	SEBND0020 X 3	
(Ht. 10"-14")	<input type="checkbox"/>	STPTA	T-Arm Ht. Adj./ Removable	\$240	<input type="checkbox"/>	STPBGL2	Gel Center Vertebrae	\$116	SEBND0020	
(Ht. 10"-14")	<input type="checkbox"/>	STPTAS	T-Arm Ht. Adj./ Removable Short	\$240	<input type="checkbox"/>	STPBGL3	Gel Left Lateral	\$116	SEBND0020	
<input type="checkbox"/>	STPGWAPF	* Gel Waterfall Arm Pads	\$180	SEAND0010 SEAND0030	<input type="checkbox"/>	STPBGL4	Gel Right Lateral	\$116	SEBND0020	
							<b>Maximum 1 piece additional pressure relief gel</b>			
							<b>*Note: Gel Waterfall Arm Pads are only compatible with locking laptrays if required</b>			
REAR WHEEL		Select one below	10	ADP	SEATING - CUSHIONS		Select one below	14	ADP	
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C		<input type="checkbox"/>	STPC01	W ( ) x D ( ) <b>RELAX</b> Soft Cushion	\$480	SESWM0050	
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C		<input type="checkbox"/>	STPC02	<b>RELAX</b> Comfy Cushion	\$624	SESWM0030	
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C		<input type="checkbox"/>	STPC03	<b>RELAX</b> Firm Cushion	\$784	SESWM0040	
(STF 13-19")	<input type="checkbox"/>	STPRWT12	* 12" Solid Tire	N/C	<input type="checkbox"/>	STPC04	<b>RELAX AIR FLOW</b> Cushion	\$587	SE0001669	
lower frame	<input type="checkbox"/>	STPRWT12I	* 12" Solid Tire Inside Mount	\$120	<input type="checkbox"/>	STPCGL2	Full Overlay Gel	\$258	SESND1045	
<input type="checkbox"/>	STPRWT16	* 16" Solid Tire	N/C		<input type="checkbox"/>	STPCCL	Contoured Leg Channel (Bolster) Std.	\$180	SEMCF203L	
<input type="checkbox"/>	STPRWA2	22" Air Tire	\$120		<input type="checkbox"/>	STPCCLA	Contoured Leg Channel (Bolster) with ABS	\$260	SESN10335 SESN10115	
<input type="checkbox"/>	STPRWA4	24" Air Tire	\$120		<input type="checkbox"/>	STPCCLC	* Contoured Leg Channel (Bolster) Custom	\$320	SEMCF209L	
<input type="checkbox"/>	STPAQR	Axles – Quick Release	\$120	WAMI			<b>*Note: Please refer to the Seating Booklet for custom sizes</b>			
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C		<input type="checkbox"/>	STPC00	Extra Cushion Cover	\$99	SESN1055	
							<b>*Note: Quick Release is not available for wheels smaller than 20"</b>			
CASTER		Select one below	11	ADP	ACCESSORIES			15	ADP	
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C		<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style)	\$26		
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C		<input type="checkbox"/>	STPSBAP	Seat Belt (Padded)	\$56	SERN00030	
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C		<input type="checkbox"/>	STPSBA4P	Seat Belt (4 Point Padded)	\$130	SERN00020	
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster	\$80		<input type="checkbox"/>	STPTRA	Full Padded Tray Locking	\$280	SETND0005 SETND0035	
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	\$80		<input type="checkbox"/>	STPTRB	Full Padded Tray Slide On	\$280	SETND0005 SETND0035	
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	\$80		<input type="checkbox"/>	STPTRCL	Full Plastic Tray Locking <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010	
<input type="checkbox"/>	STPFC4B	6" X 2 Caster	\$80		<input type="checkbox"/>	STPTRC	Full Plastic Tray Slide On <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$280	SETND0010	
<input type="checkbox"/>	STPFC4C	8" X 2 Caster	\$80		<input type="checkbox"/>	STPTRD	Half Padded Left Tray	\$240	SETND0015 SETND0035	
WHEEL LOCK		Select one below	12	ADP	<input type="checkbox"/>	STPTRE	Half Padded Right Tray	\$240	SETND0015 SETND0035	
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C		<input type="checkbox"/>	STPTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015	
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C		<input type="checkbox"/>	STPTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	\$240	SETND0015	
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C		<input type="checkbox"/>	STPTDFB	Transport Tie Downs Front	\$85		
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C		<input type="checkbox"/>	STPTDRB	Transport Tie Downs Rear	\$85		
<input type="checkbox"/>	STPBRXC	4" Brake Extension	N/C		<input type="checkbox"/>	STPOXY	Oxygen Tank Holder	\$195	WEPX	
<input type="checkbox"/>	STPBRXA	6" Brake Extension	N/C		<input type="checkbox"/>	STPIVA	Adjustable I V Pole	\$130		
<input type="checkbox"/>	STPBRXB	9" Brake Extension	\$130		<input type="checkbox"/>	STPCHPT-L	Left Arm Trough Tray	\$400	SETND0015 SETND0035	
<input type="checkbox"/>	STPATTBRX	Attendant Brakes	\$240		<input type="checkbox"/>	STPCHPT-R	Right Arm Trough Tray	\$400	SETND0030	
<input type="checkbox"/>	STPULB	Unilateral Hand Brake <input type="checkbox"/> Left <input type="checkbox"/> Right	\$120		<input type="checkbox"/>	STPLATB	Lateral Bite out	\$180	SELND0015	
					<input type="checkbox"/>	STPADD8	8" Adductor	\$164	SELND0005	
					<input type="checkbox"/>	STPADD10	10" Adductor	\$164	SELND0005	
					<input type="checkbox"/>	STPADD12	12" Adductor	\$164	SELND0005	
					<input type="checkbox"/>	STPFRCH	Padded Footplates	\$88	SEMCF803L	