

**TNP Transit Wheelchair (Manual Wheelchair)**



**INFORMATION**

Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



**FRAME COLOUR** Select one below **1 ADP**

- TNPF Frame Type - Folding N/C
- TNPR Frame Type - Rigid N/C

**FRAME COLOUR** Select one below **1 ADP**

- TNPS Frame - Grey Steel Frame 1000
- TNPB Frame - Black 1000

**CLIENT WEIGHT** Select one below **2 ADP**

- TNPWA Up to 250Lbs N/C
- TNPWAMF Exceeds 250Lbs \$300
- TNPWAMG Exceeds 350Lbs (Max 600Lbs) \$500

**SEAT WIDTH** Select one below **3 ADP**

- TNPSW16 Seat Width 16" N/C
- TNPSW18 Seat Width 18" N/C
- TNPSW20 Seat Width 20" \$300
- TNPSW22 Seat Width 22" \$300

**SEAT DEPTH** Select one below **4 ADP**

- TNPSD18 Seat Depth 18" N/C
- TNPSD20 Seat Depth 20" \$300

**FOOT REST AND FOOT PLATES** Select one below **5 ADP**

- TNPF7A 70 Degree Hanger N/C
- TNPEA Elevating Legrest \$404
- TNPAL Articulating Leg Rests \$450

**FOOT REST TYPE** Select one below **5 ADP**

- TNPFL Fixed Legrest N/C
- TNPRL Removable Legrest N/C

**ARM REST** Select one below **9 ADP**

- TNPF7A Fixed Height - Flip Back N/C
- TNPF7B Fixed Height - Non Flip Back N/C

**WHEEL LOCK** Select one below **10 ADP**

- TNPSTPH Steel Push To Lock N/C
- TNPSTPL Steel Pull To Lock N/C
- TNPBRXC 4" Brake Extension N/C
- TNPBRXA 6" Brake Extension N/C
- TNPBRXB 9" Brake Extension \$130

**SEATING - Upholstery** Select one below **11 ADP**

- TNPUP01 Blue Vinyl Upholstery N/C
  - TNPUP02 Pink Vinyl Upholstery N/C
  - TNPUP03 Green Vinyl Upholstery N/C
  - TNPUP04 Black Vinyl Upholstery N/C
  - TNPUP05 Purple Vinyl Upholstery N/C
  - TNPUPC Custom Vinyl Upholstery N/C
- Specify: \_\_\_\_\_

**LOGO** Select one below **12 ADP**

- TNPLOGO Custom Company Logo \$100

**ACCESSORIES** Select one below **12 ADP**

- TNPOXYIV Oxygen Tank and Adjustable I V Pole Combo \$325

**STANDARD NO CHARGE SPECIFICATIONS**

STEEL FRAME AND CROSS BRACES
STANDARD 18" BACK HEIGHT
EXTRUDED ALUMINUM FORKS
CHROME HAND RIMS
VINYL UPHOLSTERY
POWDER COATED FRAME
FIXED HEIGHT - FLIP BACK OR FIXED HEIGHT - NON FLIP BACK
SWING AWAY FOOTREST WITH HEEL LOOPS
PRECISION BEARINGS
ALUMINIUM MOLDED FOOTPLATES
ANTI-TIP ROLLERS
6" BRAKE EXTENSIONS
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS

**SEAT TO FLOOR SPECIFICATIONS (± 0.5")**

(STANDARD SEAT SETTING)			
STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 15.50"	22"	6"	LONG
<input type="checkbox"/> 17.50"	24"	8"	LONG

**NOTES**