


# GLYDER

FOLDING WHEELCHAIR



**ADP CODE:** WM0001061

The GLYDER is one of our newest wheelchairs, designed and constructed with accessibility and maximum mobility in mind. It is structured with lightweight frame and cross braces, horizontal and vertical versatile axles for easy transport. The chair also offers a variety of customizing options to suit different needs.



Horizontal Rear Axle Adjustability  
Ranging From 12.5" To 18.25"




Angle Adjustable Back Canes 0°, 5°,  
10° Or 15°



Dual Cross Braces



Two Point Flip Back Arm



Height Adjustable Arm



Swing In And Swing Out Front Rigging

## CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# GLYDER



## STANDARD FEATURES

- Light Weight Aluminum Frame & Cross Braces
- Multi Axle Position - Vertical & Horizontal
- 20", 22" Or 24" Urethane Molded Rear Wheels
- New Light Weight 20", 22" Or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" Or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- T-Arm Height Adjustable/ Removable
- New Flip Back Arms Which Allow Angle Back Adjustments
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 20" with Aluminum Frame & 21" > with Robust Steel Frame
- Black Padded Arm Pads
- Anti-Tip Rollers
- 6" Brake Extensions
- Overall Wheelchair Width: 27" (18" Seat Width Unfolded)  
18" (18" Seat Width folded)

SPECIFICATIONS	GLYDER
SEAT WIDTH:	14" TO 20" ALUMINUM FRAME, 21" > ROBUST STEEL FRAME
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.5" TO 18.25"
REAR WHEEL:	20", 22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	4", 5", 6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWING AWAY, DUAL SWING AWAY OR ELEVATING LEG RESTS
BACK ANGLE (ADJUSTABLE):	0°, 5°, 10° OR 15°
BACK HEIGHT RANGE:	15" TO 24"
WEIGHT CAPACITY:	Up to 600 lbs. With client weight package options
WEIGHT OF CHAIR:	31 lbs.

**GLYDER**  
FOLDING WHEELCHAIR



Learn more about this product by scanning the QR Code

