

ADP CODE: WM0001040

AADL CODE: W304

The STP is one of our most popular chairs and is designed for adjustments, maneuverability and comfort. The tilt system on this chair provides a change in seating position to redistribute pressure from one area to another.





Power Tilt Option



Center Pivot Swing In And Swing Out Front Rigging



Dual Cylinders Allows For Consistent Tilt For All Weight Capacities



Narrow Configuration (Inside Wheel Mount)



Rear Axle Adjustability



Rear Depth Adjustable Frame

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to users' specifications from either a drawing, pictures or a sample. The custom parts we make can be simple or complex and can be produced in small or large quantities based on what is requested. Contact your local sales representative with your specifications.





STANDARD FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts For Easy Tilt (Up To 55°)
- Tilt Collar lock, locked out at 33°
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Soft Set Gas Cylinder Release Mechanism
- Silver Grey or Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals
- Recessed Axle Plates
- Overall Wheelchair Width: 28" (18" Seat Width)

SPECIFICATIONS	STP
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 12.25" TO 18.25"
TILT ANGLE:	0° TO 55°
BACK ANGLE:	(ADJUSTABLE): 0°, 5°, 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS (33°)

ILT WHEELCHAIR



Learn more about this product by scanning the QR Code

