

PERSONAL FIT BACK




RELAX
SERIES

The Relax Personal Fit is an aluminum light weight formed back with hardware that allows for multiple adjustments.

Quick Reference Guide		
RELAX SERIES - PERSONAL FIT BACK		
STANDARD ADP: SE0001590	CUSTOM WIDTH ADP: SEBWM0035 & SE0001590	CUSTOM HEIGHT ADP: SEBWM0040 & SE0001590
Height: 18" Width: 16" & 18"	Height: 18" Width: All except 16" & 18"	Height: All except 18" Width: 16" & 18"
CUSTOM HEIGHT & WIDTH ADP: SE0001590 , W= SEBWM0035 & H= SEBWM0040		
Height: All except 18"	Width: All except 16" & 18"	Part Number: Refer below
PART NUMBERS		
PFBF (W X H) No Laterals		
Fill in width (W) and height (H) to get part number.		
E.g. For a W = 15" and H = 17" Back; Part number would be PFBF1517		

OVERLAYS

* Foam Pad * Visco * Gel



ADP: SEBND0015 SEBND0015 SEBND0020


* Up-charge items

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

KEY FEATURES

- Comfortable High Resilience Foam
- 3 lbs. Aluminum, Light Weight Shell
- Pre Punched holes to mount headrest and laterals
- Easy to Install and Remove
- Height Adjustable
- TEK4 Waterproof material at the bottom of back for easy clean up



X = Standard shell height 18"



Rear Back Angle



Accommodates Pelvic rotation and will allow for full back support with 25° of posterior and 20° of anterior adjustability.

PERSONAL FIT BACK

RELAX

SPECIFICATIONS	RELAX SERIES - PERSONAL FIT BACK
ADP	SE0001590, CUSTOM W= SEBWM0035 & CUSTOM H= SEBWM0040
STANDARD WIDTH RANGE	16" - 18"
CUSTOM WIDTH RANGE	15" or less and 19" or more
STANDARD HEIGHT RANGE	18"
CUSTOM HEIGHT RANGE	17" or less and 19" or more
AVAILABLE HARDWARE STYLES	One style which is universal to all our backs
ANGLE ADJUSTABLE	25° anterior adjustability, 20° of posterior adjustability
WIDTH ADJUSTABLE	2" of adjustability through the pins on the back
FITS TUBING SIZE	7/8" tubing or 1"
HEADREST MOUNT	Backs pre-drilled for headrest supports and lateral supports
FINISHED HEIGHT	Back Shell Size +2"